

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3000164								
Shigella species	64	MOTL/ 0	3000260			3000400		
Enterobacter agglomerans	35	MOTL/ 85	Shigella species	91		Shigella species	97	
			Enterobacter agglomerans	8		Shigella dysenteriae	2	
3000165 -7								
Enterobacter agglomerans	99		3000261 -3			3000404		
			Enterobacter agglomerans	99		Shigella species	96	
3000170						Shigella dysenteriae	3	
Enterobacter agglomerans	76	MOTL/ 85	3000264			3000410		
Shigella species	23	MOTL/ 0	Shigella species	64	MOTL/ 0	Shigella species	90	
			Enterobacter agglomerans	35	MOTL/ 85	Shigella dysenteriae	9	
3000171 -7								
Enterobacter agglomerans	99		3000265 -7			3000414		
			Enterobacter agglomerans	99		Shigella species	77	
3000200						Shigella dysenteriae	13	
Ewingella americana	98		3000270			Escherichia coli inaktiv	8	
Shigella species	1		Enterobacter agglomerans	76	MOTL/ 85			
			Shigella species	23	MOTL/ 0	3000420 -4		
3000204 -6						Shigella species	99	
Enterobacter agglomerans	99		3000271 -2					
			Enterobacter agglomerans	99		3000430		
3000210						Shigella species	99	
Ewingella americana	98		3000273			3000434		
			Enterobacter agglomerans	80	MOTL/ 85	Shigella species	97	
3000211 -7			Klebsiella ozaenae	19	MOTL/ 0	Enterobacter agglomerans	2	
Enterobacter agglomerans	99					3000440		
			3000274 -6			Shigella species	98	
3000220			Enterobacter agglomerans	99		3000444		
Shigella species	95					Shigella species	97	
Enterobacter agglomerans	4		3000277			Escherichia coli inaktiv	1	
			Enterobacter agglomerans	88				
3000221 -7			Klebsiella ozaenae	11		3000445 -6		
Enterobacter agglomerans	99					Enterobacter agglomerans	99	
			3000300 -7					
3000230 -7			Enterobacter agglomerans	99		3000450		
Enterobacter agglomerans	99					Shigella species	94	
			3000310 -7			Shigella dysenteriae	2	
3000240			Enterobacter agglomerans	99		Escherichia coli inaktiv	2	
Shigella species	93							
Enterobacter agglomerans	6		3000320 -7			3000451 -3		
			Enterobacter agglomerans	99		Enterobacter agglomerans	99	
3000241 -3								
Enterobacter agglomerans	99		3000330 -7			3000454		
			Enterobacter agglomerans	99		Shigella species	77	
3000244						Escherichia coli inaktiv	16	
Shigella species	71	MOTL/ 0	3000340 -7			Shigella dysenteriae	3	
Enterobacter agglomerans	28	MOTL/ 85	Enterobacter agglomerans	99				
3000245 -7			3000350 -7			3000455 -7		
Enterobacter agglomerans	99		Enterobacter agglomerans	99		Enterobacter agglomerans	99	
3000250			3000360 -7			3000460 -4		
Enterobacter agglomerans	69	MOTL/ 85	Enterobacter agglomerans	99		Shigella species	99	
Shigella species	30	MOTL/ 0						
			3000370 -7					
3000251 -7			Enterobacter agglomerans	99				
Enterobacter agglomerans	99							

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3001074 Enterobacter agglomerans	99		3001510 Proteus vulgaris	99		3002034 -6 Enterobacter agglomerans	99	
3001100 Proteus vulgaris	99		3001540 -1 Providencia alcalifaciens	99		3002037 Klebsiella rhinoscleromatis	95	
3001150 -6 Enterobacter agglomerans	99		3001750 -4 Enterobacter agglomerans	99		Enterobacter agglomerans	4	
3001170 -7 Enterobacter agglomerans	99		3001770 -6 Enterobacter agglomerans	99		3002040 Shigella species	62	MOTL/ 0
3001250 -6 Enterobacter agglomerans	99		3002000 Enterobacter agglomerans	99		Enterobacter agglomerans	37	MOTL/ 85
3001270 -7 Enterobacter agglomerans	99		3002003 Klebsiella rhinoscleromatis	99		3002041 -7 Enterobacter agglomerans	99	
3001334 -4 Enterobacter agglomerans	99		3002004 -6 Enterobacter agglomerans	99		3002050 -7 Enterobacter agglomerans	99	
3001350 -7 Enterobacter agglomerans	99		3002007 Klebsiella rhinoscleromatis	99		3002060 Shigella species	54	MOTL/ 0
3001360 -6 Enterobacter agglomerans	99		3002010 -2 Enterobacter agglomerans	99		Enterobacter agglomerans	45	MOTL/ 85
3001370 -7 Enterobacter agglomerans	99		3002013 Klebsiella rhinoscleromatis	90		3002061 -7 Enterobacter agglomerans	99	
3001400 Providencia alcalifaciens	65		Enterobacter agglomerans	9		3002070 -7 Enterobacter agglomerans	99	
Providencia rustigianii	31		3002014 -6 Enterobacter agglomerans	99		3002100 -2 Enterobacter agglomerans	99	
Proteus vulgaris	3		3002017 Klebsiella rhinoscleromatis	93		3002103 Klebsiella rhinoscleromatis	92	
3001401 Providencia alcalifaciens	99		Enterobacter agglomerans	6		Enterobacter agglomerans	7	
3001402 Providencia stuartii	99		3002020 -2 Enterobacter agglomerans	99		3002104 -6 Enterobacter agglomerans	99	
3001410 -1 Providencia alcalifaciens	99		3002023 Klebsiella rhinoscleromatis	99		3002107 Klebsiella rhinoscleromatis	94	
3001440 -1 Providencia alcalifaciens	99		3002024 -6 Enterobacter agglomerans	99		Enterobacter agglomerans	5	
3001500 Proteus vulgaris	52	GELT/ 91	3002027 Klebsiella rhinoscleromatis	98		3002110 -6 Enterobacter agglomerans	99	
Providencia rustigianii	28	GELT/ 0	3002030 -2 Enterobacter agglomerans	99		3002117 Klebsiella rhinoscleromatis	76	MOTL/ 0
Providencia alcalifaciens	18		3002033 Klebsiella rhinoscleromatis	93		Enterobacter agglomerans	23	MOTL/ 85
3001501 Providencia alcalifaciens	99		Enterobacter agglomerans	6		3002120 -2 Enterobacter agglomerans	99	
3001502 Proteus vulgaris	99		3002033 Klebsiella rhinoscleromatis	93		3002123 Klebsiella rhinoscleromatis	94	
3001074			Enterobacter agglomerans	6		Enterobacter agglomerans	5	
						3002124 -6 Enterobacter agglomerans	99	
						3002124		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3004610 -7 Enterobacter agglomerans	99		3005014 -4 Enterobacter agglomerans	99		3005230 -7 Enterobacter agglomerans	99	
3004620 -7 Enterobacter agglomerans	99		3005040 -1 Providencia alcalifaciens	99		3005240 -7 Enterobacter agglomerans	99	
3004630 -7 Enterobacter agglomerans	99		3005044 -6 Enterobacter agglomerans	99		3005250 -7 Enterobacter agglomerans	99	
3004640 -7 Enterobacter agglomerans	99		3005050 -7 Enterobacter agglomerans	99		3005260 -7 Enterobacter agglomerans	99	
3004650 -7 Enterobacter agglomerans	99		3005060 -6 Enterobacter agglomerans	99		3005270 -7 Enterobacter agglomerans	99	
3004660 -7 Enterobacter agglomerans	99		3005070 -7 Enterobacter agglomerans	99		3005300 -6 Enterobacter agglomerans	99	
3004670 -7 Enterobacter agglomerans	99		3005100 Providencia alcalifaciens	88		3005310 -7 Enterobacter agglomerans	99	
3004700 -7 Enterobacter agglomerans	99		Proteus vulgaris	11		3005320 -6 Enterobacter agglomerans	99	
3004710 -7 Enterobacter agglomerans	99		3005101 Providencia alcalifaciens	99		3005330 -7 Enterobacter agglomerans	99	
3004720 -7 Enterobacter agglomerans	99		3005110 -4 Enterobacter agglomerans	99		3005340 -7 Enterobacter agglomerans	99	
3004730 -7 Enterobacter agglomerans	99		3005130 -7 Enterobacter agglomerans	99		3005350 -7 Enterobacter agglomerans	99	
3004740 -7 Enterobacter agglomerans	99		3005140 Providencia alcalifaciens	92		3005360 -7 Enterobacter agglomerans	99	
3004750 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	7		3005370 -7 Enterobacter agglomerans	99	
3004760 -7 Enterobacter agglomerans	99		3005141 Providencia alcalifaciens	95		3005400 -1 Providencia alcalifaciens	98	
3004770 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	4		3005402 -3 Providencia stuartii	99	
3005000 Providencia alcalifaciens	88		3005142 -7 Enterobacter agglomerans	99		3005404 -5 Providencia alcalifaciens	99	
Proteus mirabilis	11		3005150 -7 Enterobacter agglomerans	99		3005406 Providencia stuartii	99	
3005001 Providencia alcalifaciens	99		3005160 -7 Enterobacter agglomerans	99		3005410 -1 Providencia alcalifaciens	99	
3005002 Providencia stuartii	99		3005170 -7 Enterobacter agglomerans	99		3005422 Providencia stuartii	99	
3005010 -1 Providencia alcalifaciens	99		3005200 Ewingella americana	99				
			3005210 -4 Enterobacter agglomerans	99				

BIOTYPE NUMBER	RP	ADD.TEST	BIOTYPE NUMBER	RP	ADD.TEST	BIOTYPE NUMBER	RP	ADD.TEST
ORGANISM NAME	%	POS. %	ORGANISM NAME	%	POS. %	ORGANISM NAME	%	POS. %
3005440 -1			3005600 -1			3006070 -7		
Providencia alcalifaciens	99		Providencia alcalifaciens	99		Enterobacter agglomerans	99	
3005442			3005602			3006100 -5		
Providencia stuartii	99		Providencia stuartii	99		Enterobacter agglomerans	99	
3005450 -1			3005634 -4			3006106		
Providencia alcalifaciens	98		Enterobacter agglomerans	99		Enterobacter agglomerans	60	GELT/ 2
3005454 -6			3005650 -7			Serratia marcescens	39	GELT/ 90
Enterobacter agglomerans	99		3005660 -6			3006107		
3005470 -6			Enterobacter agglomerans	99		Serratia marcescens	61	GELT/ 90
Enterobacter agglomerans	99		3005670 -7			Enterobacter agglomerans	38	GELT/ 2
3005500			Enterobacter agglomerans	99		3006110 -7		
Providencia alcalifaciens	88		3005710 -6			Enterobacter agglomerans	99	
Proteus vulgaris	11		Enterobacter agglomerans	99		3006120 -7		
3005501			3005730 -6			Enterobacter agglomerans	99	
Providencia alcalifaciens	99		Enterobacter agglomerans	99		3006130 -7		
3005502			3005740 -7			Enterobacter agglomerans	99	
Providencia stuartii	99		Enterobacter agglomerans	99		3006140 -7		
3005510			3005750 -7			Enterobacter agglomerans	99	
Providencia alcalifaciens	89		Enterobacter agglomerans	99		3006150 -7		
Proteus vulgaris	10		3005760 -7			Enterobacter agglomerans	99	
3005511			Enterobacter agglomerans	99		3006160 -7		
Providencia alcalifaciens	99		3005770 -7			Enterobacter agglomerans	99	
3005534			Enterobacter agglomerans	99		3006170 -3		
Enterobacter agglomerans	99		3006000 -7			Enterobacter agglomerans	99	
3005540 -1			Enterobacter agglomerans	99		3006174		
Providencia alcalifaciens	99		3006010 -7			Enterobacter agglomerans	97	
3005544			Enterobacter agglomerans	99		Enterobacter cloacae	2	
Enterobacter agglomerans	99		3006020 -7			3006175		
3005550			Enterobacter agglomerans	99		Enterobacter agglomerans	97	
Providencia alcalifaciens	90		3006030 -7			Enterobacter cloacae	2	
Enterobacter agglomerans	9		Enterobacter agglomerans	99		3006176 -7		
3005551			3006040 -7			Enterobacter agglomerans	99	
Providencia alcalifaciens	94		Enterobacter agglomerans	99		3006200		
Enterobacter agglomerans	5		3006050 -7			Ewingella americana	72	XYLO/ 13
3005552 -7			Enterobacter agglomerans	99		Enterobacter agglomerans	27	XYLO/ 93
Enterobacter agglomerans	99		3006060 -7			3006201 -7		
3005560 -6			Enterobacter agglomerans	99		Enterobacter agglomerans	99	
Enterobacter agglomerans	99		3006210 -7			Enterobacter agglomerans	99	
3005570 -7			Enterobacter agglomerans	99		3006210 -7		
Enterobacter agglomerans	99		3006210 -7			Enterobacter agglomerans	99	
3005440			3006210 -7			Enterobacter agglomerans	99	
			Enterobacter agglomerans	99		3006210		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3006750 -7 Enterobacter agglomerans	99		3007170 -7 Enterobacter agglomerans	99		3007434 -4 Enterobacter agglomerans	99	
3006760 -7 Enterobacter agglomerans	99		3007200 -4 Enterobacter agglomerans	99		3007450 -7 Enterobacter agglomerans	99	
3006770 -6 Enterobacter agglomerans	99		3007210 -7 Enterobacter agglomerans	99		3007460 -6 Enterobacter agglomerans	99	
3006777 Enterobacter agglomerans Klebsiella oxytoca	66 33	MOTL/ 85 MOTL/ 0	3007220 -6 Enterobacter agglomerans	99		3007470 -7 Enterobacter agglomerans	99	
3007000 Proteus mirabilis	99		3007230 -7 Enterobacter agglomerans	99		3007500 Proteus vulgaris	99	
3007002 Providencia stuartii	99		3007240 -7 Enterobacter agglomerans	99		3007502 Providencia stuartii	99	
3007010 -4 Enterobacter agglomerans	99		3007250 -7 Enterobacter agglomerans	99		3007510 -6 Enterobacter agglomerans	99	
3007030 -7 Enterobacter agglomerans	99		3007260 -7 Enterobacter agglomerans	99		3007530 -6 Enterobacter agglomerans	99	
3007040 -7 Enterobacter agglomerans	99		3007270 -7 Enterobacter agglomerans	99		3007540 -7 Enterobacter agglomerans	99	
3007050 -7 Enterobacter agglomerans	99		3007300 -7 Enterobacter agglomerans	99		3007550 -7 Enterobacter agglomerans	99	
3007060 -7 Enterobacter agglomerans	99		3007310 -7 Enterobacter agglomerans	99		3007560 -7 Enterobacter agglomerans	99	
3007070 -7 Enterobacter agglomerans	99		3007320 -7 Enterobacter agglomerans	99		3007570 -7 Enterobacter agglomerans	99	
3007100 -4 Enterobacter agglomerans	99		3007330 -7 Enterobacter agglomerans	99		3007610 -6 Enterobacter agglomerans	99	
3007110 -7 Enterobacter agglomerans	99		3007340 -7 Enterobacter agglomerans	99		3007630 -6 Enterobacter agglomerans	99	
3007120 -6 Enterobacter agglomerans	99		3007350 -7 Enterobacter agglomerans	99		3007640 -7 Enterobacter agglomerans	99	
3007130 -7 Enterobacter agglomerans	99		3007360 -7 Enterobacter agglomerans	99		3007650 -7 Enterobacter agglomerans	99	
3007140 -7 Enterobacter agglomerans	99		3007370 -7 Enterobacter agglomerans	99		3007660 -7 Enterobacter agglomerans	99	
3007150 -7 Enterobacter agglomerans	99		3007400 -1 Providencia alcalifaciens	99		3007670 -7 Enterobacter agglomerans	99	
3007160 -7 Enterobacter agglomerans	99		3007402 -2 Providencia stuartii	99		3007704 Enterobacter agglomerans	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3007710 -7 Enterobacter agglomerans	99		3010061 -7 Enterobacter agglomerans	99		3010250 -7 Enterobacter agglomerans	99	
3007720 -6 Enterobacter agglomerans	99		3010070 Yersinia pseudotubercul. Enterobacter agglomerans	66	HOTL/ 0	3010260 -7 Enterobacter agglomerans	99	
3007730 -7 Enterobacter agglomerans	99		3010071 -7 Enterobacter agglomerans	99		3010270 -7 Enterobacter agglomerans	99	
3007740 -7 Enterobacter agglomerans	99		3010100 Proteus vulgaris	99		3010300 -6 Enterobacter agglomerans	99	
3007750 -7 Enterobacter agglomerans	99		3010104 Enterobacter agglomerans	99		3010310 -7 Enterobacter agglomerans	99	
3007760 -7 Enterobacter agglomerans	99		3010110 Proteus vulgaris	84		3010320 -7 Enterobacter agglomerans	99	
3007770 -7 Enterobacter agglomerans	99		3010111 -7 Enterobacter agglomerans	99		3010330 -7 Enterobacter agglomerans	99	
3010000 Proteus mirabilis Yersinia pseudotubercul. Proteus vulgaris	82 10 7		3010124 Enterobacter agglomerans	99		3010340 -7 Enterobacter agglomerans	99	
3010010 Yersinia pseudotubercul. Enterobacter agglomerans	98 1		3010130 -7 Enterobacter agglomerans	99		3010350 -7 Enterobacter agglomerans	99	
3010014 -6 Enterobacter agglomerans	99		3010140 -7 Enterobacter agglomerans	99		3010360 -7 Enterobacter agglomerans	99	
3010030 Yersinia pseudotubercul. Enterobacter agglomerans	91 8		3010150 -7 Enterobacter agglomerans	99		3010370 -7 Enterobacter agglomerans	99	
3010034 -6 Enterobacter agglomerans	99		3010160 -7 Enterobacter agglomerans	99		3010400 Proteus vulgaris Morganella morganii	67 32	GELT/ 91 GELT/ 0
3010040 Yersinia pseudotubercul. Shigella species Enterobacter agglomerans	59 36 4	ESCL/ 95 ESCL/ 0	3010170 -7 Enterobacter agglomerans	99		3010450 -7 Enterobacter agglomerans	99	
3010041 -6 Enterobacter agglomerans	99		3010204 Enterobacter agglomerans	99		3010464 Enterobacter agglomerans	99	
3010050 Yersinia pseudotubercul. Enterobacter agglomerans	91 8		3010210 -7 Enterobacter agglomerans	99		3010470 -7 Enterobacter agglomerans	99	
3010051 -7 Enterobacter agglomerans	99		3010224 Enterobacter agglomerans	99		3010500 -4 Proteus vulgaris	99	
3010060 Shigella species Enterobacter agglomerans	85 14		3010230 -7 Enterobacter agglomerans	99		3010510 Proteus vulgaris	99	
3007710			3010240 -7 Enterobacter agglomerans	99		3010514 Enterobacter agglomerans	99	
						3010530 -6 Enterobacter agglomerans	99	
						3010530		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3010540 Proteus vulgaris	91		3011002 Proteus mirabilis	65	MALT/ 0	3011510 -4 Proteus vulgaris	99	
Enterobacter agglomerans	8		Proteus vulgaris	34	MALT/ 97	3011520 Proteus vulgaris	99	
3010544 -6 Enterobacter agglomerans	99		3011004 Proteus mirabilis	99		3011540 -2 Proteus vulgaris	99	
3010550 -7 Enterobacter agglomerans	99		3011010 Proteus vulgaris	59	MALT/ 97	3011550 Proteus vulgaris	99	
3010560 -6 Enterobacter agglomerans	99		Proteus mirabilis	40	MALT/ 0	3011700 Proteus vulgaris	99	
3010570 -7 Enterobacter agglomerans	99		3011020 Proteus mirabilis	99		3011774 Enterobacter agglomerans	99	
3010614 Enterobacter agglomerans	99		3011100 Proteus vulgaris	88		3012000 Proteus mirabilis	99	
3010630 -6 Enterobacter agglomerans	99		Proteus penneri	10		3012004 Enterobacter agglomerans	99	
3010640 -6 Enterobacter agglomerans	99		3011101 -4 Proteus vulgaris	99		3012010 -7 Enterobacter agglomerans	99	
3010650 -7 Enterobacter agglomerans	99		3011110 -2 Proteus vulgaris	99		3012024 Enterobacter agglomerans	99	
3010660 -6 Enterobacter agglomerans	99		3011140 Proteus vulgaris	99		3012030 -7 Enterobacter agglomerans	99	
3010670 -7 Enterobacter agglomerans	99		3011154 -4 Enterobacter agglomerans	99		3012040 -7 Enterobacter agglomerans	99	
3010710 -7 Enterobacter agglomerans	99		3011200 Proteus mirabilis	99		3012050 -7 Enterobacter agglomerans	99	
3010724 Enterobacter agglomerans	99		3011254 -4 Enterobacter agglomerans	99		3012060 -7 Enterobacter agglomerans	99	
3010730 -7 Enterobacter agglomerans	99		3011350 -4 Enterobacter agglomerans	99		3012070 -7 Enterobacter agglomerans	99	
3010740 -7 Enterobacter agglomerans	99		3011370 -7 Enterobacter agglomerans	99		3012100 Proteus vulgaris	94	
3010750 -7 Enterobacter agglomerans	99		3011400 Proteus vulgaris	85		Enterobacter agglomerans	5	
3010760 -7 Enterobacter agglomerans	99		Morganella morganii	13		3012104 -6 Enterobacter agglomerans	99	
3010770 -7 Enterobacter agglomerans	99		3011401 -4 Proteus vulgaris	99		3012110 -7 Enterobacter agglomerans	99	
3011000 Proteus mirabilis	85		3011410 Proteus vulgaris	99		3012120 -7 Enterobacter agglomerans	99	
Proteus vulgaris	14		3011440 Proteus vulgaris	99				
3010540			3011500 -6 Proteus vulgaris	99				

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3012130 -7 Enterobacter agglomerans	99		3012360 -7 Enterobacter agglomerans	99		3012640 -7 Enterobacter agglomerans	99	
3012140 -7 Enterobacter agglomerans	99		3012370 -7 Enterobacter agglomerans	99		3012650 -7 Enterobacter agglomerans	99	
3012150 -7 Enterobacter agglomerans	99		3012400 Proteus vulgaris	99		3012660 -7 Enterobacter agglomerans	99	
3012160 -7 Enterobacter agglomerans	99		3012414 Enterobacter agglomerans	99		3012670 -7 Enterobacter agglomerans	99	
3012170 -7 Enterobacter agglomerans	99		3012430 -4 Enterobacter agglomerans	99		3012700 -6 Enterobacter agglomerans	99	
3012200 -6 Enterobacter agglomerans	99		3012440 -6 Enterobacter agglomerans	99		3012710 -7 Enterobacter agglomerans	99	
3012210 -7 Enterobacter agglomerans	99		3012450 -7 Enterobacter agglomerans	99		3012720 -6 Enterobacter agglomerans	99	
3012220 -7 Enterobacter agglomerans	99		3012460 -6 Enterobacter agglomerans	99		3012730 -7 Enterobacter agglomerans	99	
3012230 -7 Enterobacter agglomerans	99		3012470 -7 Enterobacter agglomerans	99		3012740 -7 Enterobacter agglomerans	99	
3012240 -7 Enterobacter agglomerans	99		3012500 Proteus vulgaris	99		3012750 -7 Enterobacter agglomerans	99	
3012250 -7 Enterobacter agglomerans	99		3012510 Proteus vulgaris	92	7	3012760 -7 Enterobacter agglomerans	99	
3012260 -7 Enterobacter agglomerans	99		3012511 -4 Enterobacter agglomerans	99		3012770 -7 Enterobacter agglomerans	99	
3012270 -7 Enterobacter agglomerans	99		3012530 -7 Enterobacter agglomerans	99		3013000 Proteus mirabilis	81	
3012300 -7 Enterobacter agglomerans	99		3012540 -7 Enterobacter agglomerans	99		Proteus vulgaris	18	
3012310 -7 Enterobacter agglomerans	99		3012550 -7 Enterobacter agglomerans	99		3013054 -4 Enterobacter agglomerans	99	
3012320 -7 Enterobacter agglomerans	99		3012560 -7 Enterobacter agglomerans	99		3013100 -2 Proteus vulgaris	99	
3012330 -7 Enterobacter agglomerans	99		3012570 -7 Enterobacter agglomerans	99		3013110 Proteus vulgaris	99	
3012340 -7 Enterobacter agglomerans	99		3012610 -4 Enterobacter agglomerans	99		3013150 -6 Enterobacter agglomerans	99	
3012350 -7 Enterobacter agglomerans	99		3012630 -7 Enterobacter agglomerans	99		3013170 -7 Enterobacter agglomerans	99	
3012130						3013250 -6 Enterobacter agglomerans	99	
						3013250		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3014430 -7 Enterobacter agglomerans	99		3014640 -7 Enterobacter agglomerans	99		3015013 Providencia rettgeri	99	
3014440 -7 Enterobacter agglomerans	99		3014650 -7 Enterobacter agglomerans	99		3015020 Proteus mirabilis	99	
3014450 -7 Enterobacter agglomerans	99		3014660 -7 Enterobacter agglomerans	99		3015050 -7 Enterobacter agglomerans	99	
3014460 -7 Enterobacter agglomerans	99		3014670 -7 Enterobacter agglomerans	99		3015064 Enterobacter agglomerans	99	
3014470 -7 Enterobacter agglomerans	99		3014700 -7 Enterobacter agglomerans	99		3015070 -7 Enterobacter agglomerans	99	
3014500 -2 Proteus vulgaris	99		3014710 -7 Enterobacter agglomerans	99		3015100 Proteus vulgaris Proteus mirabilis	96 3	
3014504 -6 Enterobacter agglomerans	99		3014720 -7 Enterobacter agglomerans	99		3015101 -4 Proteus vulgaris	99	
3014510 Proteus vulgaris Enterobacter agglomerans	91 8		3014730 -7 Enterobacter agglomerans	99		3015110 Proteus vulgaris	99	
3014511 -7 Enterobacter agglomerans	99		3014740 -7 Enterobacter agglomerans	99		3015114 Enterobacter agglomerans	99	
3014520 -7 Enterobacter agglomerans	99		3014750 -7 Enterobacter agglomerans	99		3015130 -6 Enterobacter agglomerans	99	
3014530 -7 Enterobacter agglomerans	99		3014760 -7 Enterobacter agglomerans	99		3015140 Proteus vulgaris Enterobacter agglomerans	94 5	
3014540 -7 Enterobacter agglomerans	99		3014770 -6 Enterobacter agglomerans	99		3015144 -6 Enterobacter agglomerans	99	
3014550 -7 Enterobacter agglomerans	99		3014777 Klebsiella oxytoca Enterobacter agglomerans	54 45	MOTL/ MOTL/	0 85	3015150 -7 Enterobacter agglomerans	99
3014560 -7 Enterobacter agglomerans	99		3015000 Proteus mirabilis Proteus vulgaris	96 3		3015160 -6 Enterobacter agglomerans	99	
3014570 -7 Enterobacter agglomerans	99		3015002 Proteus mirabilis Providencia rettgeri	78 21	AROL/ AROL/	0 100	3015170 -7 Enterobacter agglomerans	99
3014600 -6 Enterobacter agglomerans	99		3015003 Providencia rettgeri	99		3015200 Proteus mirabilis	99	
3014610 -7 Enterobacter agglomerans	99		3015004 Proteus mirabilis	99		3015214 Enterobacter agglomerans	99	
3014620 -7 Enterobacter agglomerans	99		3015010 Proteus mirabilis Proteus vulgaris	78 21	HALT/ HALT/	0 97	3015230 -6 Enterobacter agglomerans	99
3014430							3015230	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3016240 -7 Enterobacter agglomerans	99		3016420 -7 Enterobacter agglomerans	99		3016640 -7 Enterobacter agglomerans	99	
3016250 -7 Enterobacter agglomerans	99		3016430 -7 Enterobacter agglomerans	99		3016650 -7 Enterobacter agglomerans	99	
3016260 -7 Enterobacter agglomerans	99		3016440 -7 Enterobacter agglomerans	99		3016660 -7 Enterobacter agglomerans	99	
3016270 -7 Enterobacter agglomerans	99		3016450 -7 Enterobacter agglomerans	99		3016670 -7 Enterobacter agglomerans	99	
3016300 -7 Enterobacter agglomerans	99		3016460 -7 Enterobacter agglomerans	99		3016700 -7 Enterobacter agglomerans	99	
3016310 -7 Enterobacter agglomerans	99		3016470 -7 Enterobacter agglomerans	99		3016710 -7 Enterobacter agglomerans	99	
3016320 -7 Enterobacter agglomerans	99		3016500 Proteus vulgaris	97		3016720 -7 Enterobacter agglomerans	99	
3016330 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	2		3016730 -7 Enterobacter agglomerans	99	
3016340 -7 Enterobacter agglomerans	99		3016501 -7 Enterobacter agglomerans	99		3016740 -7 Enterobacter agglomerans	99	
3016350 -7 Enterobacter agglomerans	99		3016510 -7 Enterobacter agglomerans	99		3016750 -7 Enterobacter agglomerans	99	
3016360 Enterobacter agglomerans	99		3016520 -7 Enterobacter agglomerans	99		3016760 -7 Enterobacter agglomerans	99	
3016361 Enterobacter agglomerans	72	GELT/ 2	3016530 -7 Enterobacter agglomerans	99		3016770 -6 Enterobacter agglomerans	99	
Serratia rubidaea	27	GELT/ 90	3016540 -7 Enterobacter agglomerans	99		3016777 Klebsiella oxytoca	81	
3016362 -7 Enterobacter agglomerans	99		3016550 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	18	
3016370 -6 Enterobacter agglomerans	99		3016560 -7 Enterobacter agglomerans	99		3017000 Proteus mirabilis	95	
3016377 Klebsiella pneumoniae	57	MOTL/ 0	3016570 -7 Enterobacter agglomerans	99		Proteus vulgaris	4	
Enterobacter agglomerans	42	MOTL/ 85	3016600 -7 Enterobacter agglomerans	99		3017014 Enterobacter agglomerans	99	
3016400 Proteus vulgaris	90		3016610 -7 Enterobacter agglomerans	99		3017030 -4 Enterobacter agglomerans	99	
Enterobacter agglomerans	9		3016620 -7 Enterobacter agglomerans	99		3017040 -6 Enterobacter agglomerans	99	
3016404 -6 Enterobacter agglomerans	99		3016630 -7 Enterobacter agglomerans	99		3017050 -7 Enterobacter agglomerans	99	
3016410 -7 Enterobacter agglomerans	99							

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3020250 Hafnia alvei	99		3022250 Hafnia alvei	99		3024146 Serratia liquefaciens		
3020264 -6 Serratia liquefaciens	99		3022264 -6 Serratia liquefaciens	99		Serratia marcescens	59	XYLO/100 40 XYLO/ 7
3020360 -7 Serratia liquefaciens	99		3022360 -7 Serratia liquefaciens	99		3024147 Serratia marcescens	99	
3020370 -7 Serratia liquefaciens	99		3022374 Serratia liquefaciens	99		3024154 -6 Serratia liquefaciens	99	
3020400 Shigella species	99		3022376 Serratia liquefaciens	99	78 CELL/ 5 21 CELL/ 99	3024160 -7 Serratia liquefaciens	99	
3020420 Shigella species	99		3022377 Enterobacter aerogenes	99		3024170 -2 Serratia liquefaciens	99	
3020440 Shigella species	99		3023000 Proteus mirabilis	99		3024174 Serratia liquefaciens	74	CELL/ 5 25 CELL/ 99
3020444 Shigella species	78		3024000 Proteus mirabilis	99		3024175 Enterobacter cloacae	78	CELL/ 99
Escherichia coli inaktiv	21		3024002 -7 Serratia marcescens	99		Serratia liquefaciens	21	CELL/ 5
3020454 Escherichia coli inaktiv	99		3024050 Hafnia alvei	99		3024176 Serratia liquefaciens	99	
3020460 -4 Shigella species	99		3024064 -6 Serratia liquefaciens	99		3024177 Enterobacter aerogenes	90	
3020654 Citrobacter amalonaticus	99		3024070 Hafnia alvei	99		Serratia liquefaciens	9	
3020764 -6 Serratia liquefaciens	99		3024074 -6 Serratia liquefaciens	99		3024200 Ewingella americana	99	
3021000 Proteus mirabilis	99		3024100 -7 Serratia marcescens	99		3024224 -6 Serratia liquefaciens	99	
3021200 Proteus mirabilis	99		3024124 Serratia liquefaciens	99		3024250 Hafnia alvei	99	
3021400 Morganella morganii	99		3024126 Serratia liquefaciens	99	59 XYLO/100 40 XYLO/ 7	3024260 -7 Serratia liquefaciens	99	
3022050 Hafnia alvei	99		3024127 Serratia marcescens	99		3024270 -4 Serratia liquefaciens	99	
3022070 Hafnia alvei	99		3024134 -4 Serratia liquefaciens	99		3024305 Serratia marcescens	99	
3022107 Serratia marcescens	99					3024306 Serratia marcescens	80	
3022164 -6 Serratia liquefaciens	99					Serratia liquefaciens	19	
						3024307 Serratia marcescens	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3024320 -7 Serratia liquefaciens	99		3024774 Citrobacter amalonaticus	68	CELL/100	3026165 Enterobacter cloacae	67	CELL/ 99
3024330 -6 Serratia liquefaciens	99		Serratia liquefaciens	31	CELL/ 5	Serratia liquefaciens	32	CELL/ 5
3024340 -7 Serratia liquefaciens	99		3024776 Serratia liquefaciens	99		3026166 -7 Serratia liquefaciens	99	
3024350 Serratia liquefaciens	99		3025000 -4 Proteus mirabilis	99		3026170 Enterobacter cloacae	99	
3024354 Serratia liquefaciens	97		3025010 Proteus mirabilis	99		3026174 Enterobacter cloacae	83	
Enterobacter agglomerans	2		3025020 Proteus mirabilis	99		Serratia liquefaciens	15	
3024356 Serratia liquefaciens	99		3025100 Proteus mirabilis	99		3026175 Enterobacter cloacae	74	AROL/ 15
3024360 -7 Serratia liquefaciens	99		3025200 Proteus mirabilis	99		Enterobacter aerogenes	25	AROL/100
3024370 -2 Serratia liquefaciens	99		3025400 Providencia alcalifaciens	78	GELT/ 0	3026176 Enterobacter cloacae	57	CELL/ 99
3024374 Serratia liquefaciens	91		Proteus mirabilis	21	GELT/ 90	Serratia liquefaciens	42	CELL/ 5
Enterobacter cloacae	8		3025401 Providencia alcalifaciens	99		3026177 Enterobacter aerogenes	91	
3024375 Serratia liquefaciens	33	AROL/ 0	3026000 Proteus mirabilis	99		Enterobacter cloacae	8	
Enterobacter aerogenes	33	AROL/100	3026006 -7 Serratia marcescens	99		3026250 Hafnia alvei	99	
Enterobacter cloacae	33		3026050 Hafnia alvei	99		3026260 -6 Serratia liquefaciens	99	
3024376 Serratia liquefaciens	93		3026064 -6 Serratia liquefaciens	99		3026274 Enterobacter cloacae	50	CELL/ 99
Enterobacter aerogenes	3		3026070 Hafnia alvei	99		Serratia liquefaciens	45	CELL/ 5
Enterobacter cloacae	2		3026100 -7 Serratia marcescens	99		Enterobacter agglomerans	4	
3024377 Enterobacter aerogenes	86		3026124 -6 Serratia liquefaciens	99		3026276 Serratia liquefaciens	99	
Serratia liquefaciens	11		3026160 -2 Serratia liquefaciens	99		3026277 Enterobacter aerogenes	99	
Enterobacter cloacae	3		3026164 Serratia liquefaciens	84		3026307 Serratia marcescens	99	
3024454 -4 Citrobacter amalonaticus	99		Enterobacter cloacae	15		3026320 -6 Serratia liquefaciens	99	
3024506 -7 Serratia marcescens	99		3026350 Enterobacter agglomerans	99		3026340 -6 Serratia liquefaciens	99	
3024564 Serratia liquefaciens	99					3026350 Enterobacter agglomerans	99	
3024654 -4 Citrobacter amalonaticus	99							
3024764 -6 Serratia liquefaciens	99							

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3026354								
Enterobacter cloacae	54	CELL/ 99	3026375			3031010		
Serratia liquefaciens	38	CELL/ 5	Enterobacter aerogenes	50	AROL/100	Proteus mirabilis	99	
Enterobacter agglomerans	6		Enterobacter cloacae	46	AROL/ 15			
			Serratia liquefaciens	3		3031020		
3026355						Proteus mirabilis	99	
Enterobacter cloacae	89		3026376					
Enterobacter agglomerans	10		Serratia liquefaciens	50	AROL/ 0	3031040		
			Enterobacter aerogenes	32	AROL/100	Proteus mirabilis	99	
3026356			Enterobacter cloacae	17				
Serratia liquefaciens	90		3026377			3031100		
Enterobacter agglomerans	9		Enterobacter aerogenes	96		Proteus mirabilis	99	
			Enterobacter cloacae	3		3031200		
3026357						Proteus mirabilis	99	
Enterobacter aerogenes	99		3026455 -5					
			Citrobacter diversus	99		3031400		
3026360 -3						Morganella morganii	91	
Serratia liquefaciens	99		3026764 -6			Proteus mirabilis	8	
			Serratia liquefaciens	99				
3026364						3031404		
Serratia liquefaciens	94		3027000			Morganella morganii	99	
Enterobacter cloacae	5		Proteus mirabilis	99				
						3031500		
3026365			3030000 -4			Morganella morganii	73	GELT/ 0
Serratia liquefaciens	64	CELL/ 5	Proteus mirabilis	99		Proteus vulgaris	26	GELT/ 91
Enterobacter cloacae	35	CELL/ 99						
			3030010			3032000		
3026366			Proteus mirabilis	99		Proteus mirabilis	99	
Serratia liquefaciens	98							
Enterobacter cloacae	1		3030020			3032170 -2		
			Proteus mirabilis	99		Enterobacter gergoviae	99	
3026367								
Enterobacter aerogenes	59	AROL/100	3030100			3032364 -6		
Serratia liquefaciens	40	AROL/ 0	Proteus mirabilis	99		Serratia liquefaciens	99	
3026370			3030106 -7			3032370 -2		
Enterobacter cloacae	59	CELL/ 99	Serratia marcescens	99		Enterobacter gergoviae	99	
Serratia liquefaciens	37	CELL/ 5						
Enterobacter agglomerans	3		3030164 -6			3032374		
			Serratia liquefaciens	99		Serratia liquefaciens	99	
3026371								
Enterobacter cloacae	94		3030200			3033000 -4		
Enterobacter agglomerans	5		Proteus mirabilis	99		Proteus mirabilis	99	
3026372			3030360 -6			3033010		
Serratia liquefaciens	94		Serratia liquefaciens	99		Proteus mirabilis	99	
Enterobacter agglomerans	5							
			3030400			3033020		
3026373			Morganella morganii	94		Proteus mirabilis	99	
Enterobacter aerogenes	99		Proteus mirabilis	5				
						3033100		
3026374			3030544			Proteus mirabilis	99	
Enterobacter cloacae	58	CELL/ 99	Yersinia enterocolitica	99				
Serratia liquefaciens	41	CELL/ 5				3033200		
			3031000 -4			Proteus mirabilis	99	
			Proteus mirabilis	99				

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3033400								
Morganella morganii	59	GELT/ 0	3034506 -7			3036200		
Proteus mirabilis	40	GELT/ 90	Serratia marcescens	99		Proteus mirabilis	99	
3034000 -4			3035000 -4			3036264 -6		
Proteus mirabilis	99		Proteus mirabilis	99		Serratia liquefaciens	99	
3034005 -7			3035010			3036354		
Serratia marcescens	99		Proteus mirabilis	99		Enterobacter agglomerans	99	
3034010			3035020			3036360 -7		
Proteus mirabilis	99		Proteus mirabilis	99		Serratia liquefaciens	99	
3034020			3035040			3036370 -2		
Proteus mirabilis	99		Proteus mirabilis	99		Enterobacter gergoviae	99	
3034100			3035100			3036374		
Proteus mirabilis	83		Proteus mirabilis	99		Enterobacter cloacae	71 CELL/ 99	
Serratia marcescens	16					Serratia liquefaciens	26 CELL/ 5	
3034101 -7			3035200 -4			Enterobacter agglomerans	2	
Serratia marcescens	99		Proteus mirabilis	99		3036375		
3034164 -6			3035210			Enterobacter cloacae	99	
Serratia liquefaciens	99		Proteus mirabilis	99		3036376		
3034200			3035220			Serratia liquefaciens	59 CELL/ 5	
Proteus mirabilis	99		Proteus mirabilis	99		Enterobacter cloacae	40 CELL/ 99	
3034264 -6			3035300			3036377		
Serratia liquefaciens	99		Proteus mirabilis	99		Enterobacter aerogenes	99	
3034306 -7			3035400			3037000 -4		
Serratia marcescens	99		Proteus mirabilis	91		Proteus mirabilis	99	
3034324 -6			Morganella morganii	8		3037010		
Serratia liquefaciens	99		3035500			Proteus mirabilis	99	
3034360 -7			Proteus vulgaris	99		3037020		
Serratia liquefaciens	99		3035600			Proteus mirabilis	99	
3034370 -2			Proteus mirabilis	99		3037040		
Serratia liquefaciens	99		3036000			Proteus mirabilis	99	
3034374			Proteus mirabilis	99		3037100		
Serratia liquefaciens	85		3036006 -7			Proteus mirabilis	99	
Enterobacter cloacae	14		Serratia marcescens	99		3037200		
3034375 -6			3036102 -7			Proteus mirabilis	99	
Serratia liquefaciens	99		Serratia marcescens	99		3037400		
3034377			3036164 -6			Proteus mirabilis	99	
Enterobacter aerogenes	89		Serratia liquefaciens	99		3040040		
Serratia liquefaciens	10		3036170 -2			Shigella species	99	
3034400			Enterobacter gergoviae	99		3040060		
Proteus mirabilis	99		3036174 -5			Shigella species	99	
3033400			Enterobacter cloacae	99		3040060		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3042374 -5 Enterobacter cloacae	99		3046372 Enterobacter agglomerans	99		3064154 -5 Enterobacter cloacae	99	
3044200 Ewingella americana	99		3046374 -5 Enterobacter cloacae	98		3064170 -7 Enterobacter cloacae	99	
3044354 Enterobacter agglomerans	99		3046376 Enterobacter cloacae	97		3064274 -7 Enterobacter cloacae	99	
3044364 Enterobacter cloacae	99		Enterobacter agglomerans	2		3064370 Enterobacter sakazakii	68	
3044374 Enterobacter cloacae	98		3046377 Enterobacter cloacae	97		Enterobacter cloacae	31	
Enterobacter agglomerans	1		Enterobacter agglomerans	2		3064371 Enterobacter cloacae	99	
3044375 -7 Enterobacter cloacae	99		3051000 Proteus mirabilis	99		3064372 Enterobacter sakazakii	95	
3046164 -5 Enterobacter cloacae	99		3051500 Proteus vulgaris	99		Enterobacter cloacae	4	
3046170 -7 Enterobacter cloacae	99		3054374 -5 Enterobacter cloacae	99		3064374 -7 Enterobacter cloacae	99	
3046274 Enterobacter cloacae	93		3055000 Proteus mirabilis	99		3064654 -4 Citrobacter amalonaticus	99	
Enterobacter agglomerans	6		3055500 Proteus vulgaris	99		3064772 Enterobacter sakazakii	99	
3046275 Enterobacter cloacae	99		3056174 -5 Enterobacter cloacae	99		3065000 Proteus mirabilis	99	
3046350 Enterobacter agglomerans	99		3056354 Enterobacter agglomerans	99		3066074 -6 Enterobacter cloacae	99	
3046354 Enterobacter cloacae	90		3056364 Enterobacter cloacae	99		3066106 -7 Serratia marcescens	99	
Enterobacter agglomerans	9		3056374 Enterobacter cloacae	97		3066134 -6 Enterobacter cloacae	99	
3046355 Enterobacter cloacae	91		Enterobacter agglomerans	2		3066160 -7 Enterobacter cloacae	99	
Enterobacter agglomerans	8		3056375 -7 Enterobacter cloacae	99		3066170 -7 Enterobacter cloacae	99	
3046356 Enterobacter agglomerans	99		3060054 Salmonella paratyphi A	99	SEROLOGY	3066264 -5 Enterobacter cloacae	99	
3046364 -6 Enterobacter cloacae	99		3060374 -7 Enterobacter cloacae	99		3066270 -7 Enterobacter cloacae	99	
3046370 Enterobacter cloacae	95		3064074 Enterobacter cloacae	99		3066334 -5 Enterobacter cloacae	99	
Enterobacter agglomerans	4		3064102 -7 Serratia marcescens	99				
3046371 Enterobacter cloacae	95							
Enterobacter agglomerans	4							

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3104275 -7 Enterobacter agglomerans	99		3104373 Enterobacter agglomerans	99		3106040 Hafnia alvei	98	
3104307 Serratia marcescens	99		3104374 Enterobacter agglomerans	99		Enterobacter agglomerans	1	
3104310 -7 Enterobacter agglomerans	99		3104375 Enterobacter agglomerans	99		3106044 -6 Enterobacter agglomerans	99	
3104324 Enterobacter agglomerans	99		3104376 Enterobacter agglomerans	99		3106050 Hafnia alvei	98	
3104330 -7 Enterobacter agglomerans	99		3104377 Klebsiella pneumoniae	91		3106051 -3 Enterobacter agglomerans	99	
3104340 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	8		3106054 Hafnia alvei	84	
3104350 -7 Enterobacter agglomerans	99		3104454 -4 Enterobacter agglomerans	99		Enterobacter agglomerans	15	
3104360 Serratia rubidaea	95		3104550 -6 Enterobacter agglomerans	99		3106055 -7 Enterobacter agglomerans	99	
Enterobacter agglomerans	4		3104570 -7 Enterobacter agglomerans	99		3106060 Hafnia alvei	94	
3104361 Serratia rubidaea	98		3104650 -6 Enterobacter agglomerans	99		Enterobacter agglomerans	5	
3104362 Serratia rubidaea	93		3104670 -7 Enterobacter agglomerans	99		3106061 -6 Enterobacter agglomerans	99	
Enterobacter agglomerans	6		3104714 -4 Enterobacter agglomerans	99		3106070 Hafnia alvei	97	
3104363 Serratia rubidaea	98		3104740 -4 Enterobacter agglomerans	99		Enterobacter agglomerans	2	
Enterobacter agglomerans	1		3104750 -7 Enterobacter agglomerans	99		3106071 -3 Enterobacter agglomerans	99	
3104364 Enterobacter agglomerans	59	GELT/ 2	3104760 -6 Enterobacter agglomerans	99		3106074 Hafnia alvei	58	
Serratia liquefaciens	40	GELT/ 90	3104770 -6 Enterobacter agglomerans	99		Enterobacter agglomerans	42	
3104365 Enterobacter agglomerans	99		3104777 Klebsiella oxytoca	97		3106075 -7 Enterobacter agglomerans	99	
3104366 Serratia liquefaciens	51	GELT/ 90	Enterobacter agglomerans	2		3106102 -7 Serratia marcescens	99	
Enterobacter agglomerans	48	GELT/ 2	3105374 -4 Enterobacter agglomerans	99		3106110 -7 Enterobacter agglomerans	99	
3104367 Enterobacter agglomerans	99		3106030 -6 Enterobacter agglomerans	99		3106124 Enterobacter agglomerans	99	
3104370 Enterobacter agglomerans	99		3106130 -7 Enterobacter agglomerans	99		3106140 -7 Enterobacter agglomerans	99	
3104371 Enterobacter agglomerans	99							
3104372 Enterobacter agglomerans	99							
3104275								

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3106150 -7 Enterobacter agglomerans	99		3106271 -3 Enterobacter agglomerans	99		3106364 Serratia rubidaea	68	GELT/ 90
3106160 Enterobacter agglomerans	99		3106274 Enterobacter agglomerans	83		Enterobacter agglomerans	31	GELT/ 2
3106161 Serratia rubidaea	89		Hafnia alvei	16		3106365 Serratia rubidaea	91	
Enterobacter agglomerans	10		3106275 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	8	
3106162 -7 Enterobacter agglomerans	99		3106300 -6 Enterobacter agglomerans	99		3106366 Enterobacter agglomerans	99	
3106170 -7 Enterobacter agglomerans	99		3106310 -7 Enterobacter agglomerans	99		3106367 Serratia rubidaea	86	
3106210 -7 Enterobacter agglomerans	99		3106320 Enterobacter agglomerans	99		Enterobacter agglomerans	13	
3106224 Enterobacter agglomerans	99		3106321 Serratia rubidaea	95		3106370 Enterobacter agglomerans	99	
3106230 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	4		3106371 Enterobacter agglomerans	51	GELT/ 2
3106240 Hafnia alvei	90		3106322 -6 Enterobacter agglomerans	99		Serratia rubidaea	48	GELT/ 90
Enterobacter agglomerans	9		3106330 -7 Enterobacter agglomerans	99		3106372 -6 Enterobacter agglomerans	99	
3106241 -7 Enterobacter agglomerans	99		3106340 Enterobacter agglomerans	99		3106377 Klebsiella pneumoniae	86	
3106250 Hafnia alvei	95		3106341 Serratia rubidaea	82		Enterobacter aerogenes	12	
Enterobacter agglomerans	4		Enterobacter agglomerans	17		Enterobacter agglomerans	2	
3106251 -3 Enterobacter agglomerans	99		3106342 -7 Enterobacter agglomerans	99		3106450 -6 Enterobacter agglomerans	99	
3106254 Enterobacter agglomerans	56		3106350 -7 Enterobacter agglomerans	99		3106470 -7 Enterobacter agglomerans	99	
Hafnia alvei	43		3106360 Serratia rubidaea	98		3106534 -4 Enterobacter agglomerans	99	
3106255 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	1		3106550 -7 Enterobacter agglomerans	99	
3106260 Hafnia alvei	71		3106361 Serratia rubidaea	99		3106560 -6 Enterobacter agglomerans	99	
Enterobacter agglomerans	28		3106362 Serratia rubidaea	97		3106570 -7 Enterobacter agglomerans	99	
3106261 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	2		3106634 -4 Enterobacter agglomerans	99	
3106270 Hafnia alvei	85		3106363 Serratia rubidaea	99		3106650 -7 Enterobacter agglomerans	99	
Enterobacter agglomerans	14					Enterobacter agglomerans	99	
						3106660 -6 Enterobacter agglomerans	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3114330 -6 Enterobacter agglomerans	99		3116070 Hafnia alvei	82		3116270 Enterobacter agglomerans	59	
3114340 -6 Enterobacter agglomerans	99		Enterobacter agglomerans	17		Hafnia alvei	40	
3114350 -7 Enterobacter agglomerans	99		3116071 -7 Enterobacter agglomerans	99		3116271 -6 Enterobacter agglomerans	99	
3114360 -7 Enterobacter agglomerans	99		3116106 -7 Serratia marcescens	99		3116277 Klebsiella pneumoniae	94	
3114370 -6 Enterobacter agglomerans	99		3116114 Enterobacter agglomerans	99		Enterobacter agglomerans	5	
3114377 Klebsiella pneumoniae	98		3116130 -6 Enterobacter agglomerans	99		3116304 Enterobacter agglomerans	99	
3114554 Enterobacter agglomerans	99		3116140 -6 Enterobacter agglomerans	99		3116310 -7 Enterobacter agglomerans	99	
3114570 -4 Enterobacter agglomerans	99		3116150 -7 Enterobacter agglomerans	99		3116324 Enterobacter agglomerans	99	
3114670 -4 Enterobacter agglomerans	99		3116160 -7 Enterobacter agglomerans	99		3116330 -6 Enterobacter agglomerans	99	
3114750 -4 Enterobacter agglomerans	99		3116170 Enterobacter gergoviae	73		3116337 Klebsiella pneumoniae	95	
3114770 -6 Enterobacter agglomerans	99		Enterobacter agglomerans	26		Enterobacter agglomerans	4	
3114777 Klebsiella oxytoca	99		3116171 -6 Enterobacter agglomerans	99		3116340 -7 Enterobacter agglomerans	99	
3115000 Proteus mirabilis	99		3116177 Klebsiella pneumoniae	94		3116350 -6 Enterobacter agglomerans	99	
3115500 Proteus vulgaris	99		Enterobacter agglomerans	5		3116357 Klebsiella pneumoniae	84	
3116044 Enterobacter agglomerans	99		3116214 Enterobacter agglomerans	99		Enterobacter agglomerans	15	
3116050 Hafnia alvei	94		3116230 -6 Enterobacter agglomerans	99		3116360 Serratia rubidaea	87	
Enterobacter agglomerans	5		3116240 -6 Enterobacter agglomerans	99		Enterobacter agglomerans	12	
3116051 -7 Enterobacter agglomerans	99		3116250 Hafnia alvei	72		3116361 Serratia rubidaea	97	
3116064 Enterobacter agglomerans	99		Enterobacter agglomerans	27		Enterobacter agglomerans	2	
			3116251 -7 Enterobacter agglomerans	99		3116362 Enterobacter agglomerans	99	
			3116260 -7 Enterobacter agglomerans	99		3116363 Serratia rubidaea	95	
						Enterobacter agglomerans	4	
						3116364 -6 Enterobacter agglomerans	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3120276 -6 Serratia liquefaciens	99		3122060 -4 Hafnia alvei	99		3122360 -6 Serratia liquefaciens	99	
3120350 Hafnia alvei	99		3122070 -4 Hafnia alvei	99		3122367 Enterobacter aerogenes	76	AROL/100
3120354 -6 Serratia liquefaciens	99		3122103 -7 Serratia marcescens	99		Serratia liquefaciens	23	AROL/ 0
3120360 -7 Serratia liquefaciens	99		3122150 Hafnia alvei	99		3122370 -2 Serratia liquefaciens	99	
3120370 -6 Serratia liquefaciens	99		3122164 -6 Serratia liquefaciens	99		3122373 Enterobacter aerogenes	99	
3120377 Enterobacter aerogenes	94		3122170 Hafnia alvei	99		3122374 Serratia liquefaciens	99	
Serratia liquefaciens	5		3122174 -6 Serratia liquefaciens	99		3122375 Enterobacter aerogenes	99	
3120400 Edwardsiella tarda	99		3122177 Enterobacter aerogenes	99		3122376 Enterobacter aerogenes	59	AROL/100
3120450 Hafnia alvei	99		3122210 Hafnia alvei	99		Serratia liquefaciens	40	AROL/ 0
3120454 Escherichia coli inaktiv	99		3122240 -4 Hafnia alvei	99		3122377 Enterobacter aerogenes	99	
3120470 Hafnia alvei	99		3122250 -4 Hafnia alvei	99		3122450 Hafnia alvei	99	
3120650 Hafnia alvei	99		3122260 Hafnia alvei	99		3122470 Hafnia alvei	99	
3120764 -6 Serratia liquefaciens	99		3122264 Hafnia alvei	52	ESCL/ 7	3122650 Hafnia alvei	99	
3121050 Hafnia alvei	99		Serratia liquefaciens	47	ESCL/ 97	3122650 Hafnia alvei	99	
3121070 Hafnia alvei	99		3122266 Serratia liquefaciens	99		3123050 Hafnia alvei	99	
3121250 Hafnia alvei	99		3122270 -4 Hafnia alvei	99		3123070 Hafnia alvei	99	
3122010 Hafnia alvei	99		3122277 Enterobacter aerogenes	99		3123250 Hafnia alvei	99	
3122030 Hafnia alvei	99		3122324 -6 Serratia liquefaciens	99		3124000 -7 Serratia marcescens	99	
3122040 -4 Hafnia alvei	99		3122350 Hafnia alvei	99		3124010 Hafnia alvei	99	
3122050 -4 Hafnia alvei	99		3122357 Enterobacter aerogenes	99		3124014 Salmonella choleraesuis	99	SEROLOGY

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3124350 Hafnia alvei	56	ESCL/ 7	3124760 -6 Serratia liquefaciens	99		3126126 Serratia marcescens		73 XYLO/ 7
Serratia liquefaciens	43	ESCL/ 97				Serratia liquefaciens		26 XYLO/100
3124352 -6 Serratia liquefaciens	99		3125000 Proteus mirabilis	99		3126127 Serratia marcescens		99
3124357 Enterobacter aerogenes	99		3125050 Hafnia alvei	99		3126144 Serratia liquefaciens		99
3124360 -6 Serratia liquefaciens	99		3125070 Hafnia alvei	99		3126146 Serratia marcescens		73 XYL/ 7
3124367 Serratia liquefaciens	66	AROL/ 0	3125250 Hafnia alvei	99		Serratia liquefaciens		26 XYLO/100
Enterobacter aerogenes	33	AROL/100				3126147 Serratia marcescens		99
3124370 -2 Serratia liquefaciens	99		3125364 -6 Serratia liquefaciens	99		3126150 Hafnia alvei		99
3124373 Enterobacter aerogenes	96		3126002 -7 Serratia marcescens	99		3126157 Enterobacter aerogenes		99
Serratia liquefaciens	3					3126160 -6 Serratia liquefaciens		99
3124374 Serratia liquefaciens	96		3126010 Hafnia alvei	99		3126167 Enterobacter aerogenes		90
Enterobacter aerogenes	3					Serratia liquefaciens		9
3124375 Enterobacter aerogenes	84		3126030 Hafnia alvei	99		3126170 Enterobacter gergoviae		79 AROL/ 97
Serratia liquefaciens	15					Hafnia alvei		20 AROL/ 0
3124376 Serratia liquefaciens	81		3126040 -4 Hafnia alvei	99		3126173 Enterobacter aerogenes		99
Enterobacter aerogenes	18					3126174 Serratia liquefaciens		99
3124377 Enterobacter aerogenes	97		3126050 -4 Hafnia alvei	99		3126175 Enterobacter aerogenes		99
Serratia liquefaciens	2					3126176 Enterobacter aerogenes		81
3124407 Serratia marcescens	99		3126064 Hafnia alvei	79	ESCL/ 7	Serratia liquefaciens		18
3124450 Hafnia alvei	99		Serratia liquefaciens	20	ESCL/ 97			
3124470 Hafnia alvei	99		3126066 Serratia liquefaciens	99		3126177 Enterobacter aerogenes		99
3124502 -7 Serratia marcescens	99					3126210 Hafnia alvei		99
3124564 -6 Serratia liquefaciens	99		3126070 -4 Hafnia alvei	99		3126240 -4 Hafnia alvei		99
3124650 Hafnia alvei	99							
			3126077 Enterobacter aerogenes	99				
			3126100 -7 Serratia marcescens	99				
			3126124 Serratia liquefaciens	99				

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3136003 -7 Serratia marcescens	99		3136270 -4 Hafnia alvei	99		3140040 Hafnia alvei	99	
3136040 Hafnia alvei	99		3136277 Enterobacter aerogenes	99		3140050 -4 Hafnia alvei	99	
3136050 -4 Hafnia alvei	99		3136306 -7 Serratia marcescens	99		3140060 Hafnia alvei	99	
3136060 Hafnia alvei	99		3136324 -6 Serratia liquefaciens	99		3140070 Hafnia alvei	99	
3136070 Hafnia alvei	92		3136357 Enterobacter aerogenes	99		3140240 Hafnia alvei	99	
Enterobacter gergoviae	7							
3136074 Hafnia alvei	99		3136360 -6 Serratia liquefaciens	99		3140250 Hafnia alvei	99	
3136100 -7 Serratia marcescens	99		3136367 Enterobacter aerogenes	90		3140270 Hafnia alvei	99	
			Serratia liquefaciens	9				
3136127 Serratia marcescens	99		3136370 -2 Enterobacter gergoviae	99		3142040 Hafnia alvei	99	
3136130 Enterobacter gergoviae	99		3136373 Enterobacter aerogenes	99		3142050 -4 Hafnia alvei	99	
3136147 Serratia marcescens	99		3136374 Serratia liquefaciens	99		3142060 Hafnia alvei	99	
3136150 Enterobacter gergoviae	99		3136375 Enterobacter aerogenes	99		3142070 Hafnia alvei	99	
3136160 Enterobacter gergoviae	99		3136376 Enterobacter aerogenes	80	AROL/100	3142240 Hafnia alvei	99	
3136164 -6 Serratia liquefaciens	99		Serratia liquefaciens	19	AROL/ 0	3142250 Hafnia alvei	99	
3136170 -4 Enterobacter gergoviae	99		3136377 Enterobacter aerogenes	99		3142270 Hafnia alvei	99	
3136177 Enterobacter aerogenes	99		3136506 -7 Serratia marcescens	99		3144040 Hafnia alvei	99	
3136240 Hafnia alvei	99		3136570 Enterobacter gergoviae	99		3144050 -4 Hafnia alvei	99	
3136250 -4 Hafnia alvei	99		3136777 Klebsiella ornithinolyt.	99		3144060 Hafnia alvei	99	
3136260 Hafnia alvei	99		3137000 Proteus mirabilis	99		3144070 Hafnia alvei	99	
3136264 -6 Serratia liquefaciens	99		3137170 Enterobacter gergoviae	99		3144106 -7 Serratia marcescens	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3144240 Hafnia alvei	99		3160060 -4 Hafnia alvei	99		3164060 -4 Hafnia alvei	99	
3144250 Hafnia alvei	99		3160070 -4 Hafnia alvei	99		3164070 -4 Hafnia alvei	99	
3144270 Hafnia alvei	99		3160106 -7 Serratia marcescens	99		3164100 -7 Serratia marcescens	99	
3146040 Hafnia alvei	99		3160240 -4 Hafnia alvei	99		3164240 -4 Hafnia alvei	99	
3146050 -4 Hafnia alvei	99		3160250 -4 Hafnia alvei	99		3164250 -4 Hafnia alvei	99	
3146060 Hafnia alvei	99		3160260 Hafnia alvei	99		3164260 Hafnia alvei	99	
3146070 Hafnia alvei	99		3160270 -4 Hafnia alvei	99		3164270 -4 Hafnia alvei	99	
3146240 Hafnia alvei	99		3162040 -4 Hafnia alvei	99		3164306 -7 Serratia marcescens	99	
3146250 Hafnia alvei	99		3162050 -4 Hafnia alvei	99		3164364 -6 Serratia liquefaciens	99	
3146270 Hafnia alvei	99		3162060 -4 Hafnia alvei	99		3164507 Serratia marcescens	99	
3150050 Hafnia alvei	99		3162070 -4 Hafnia alvei	99		3166007 Serratia marcescens	99	
3152050 Hafnia alvei	99		3162240 -4 Hafnia alvei	99		3166040 -4 Hafnia alvei	99	
3154050 Hafnia alvei	99		3162250 -4 Hafnia alvei	99		3166050 -4 Hafnia alvei	99	
3156050 Hafnia alvei	99		3162260 Hafnia alvei	99		3166060 -4 Hafnia alvei	99	
3156377 Klebsiella pneumoniae	99		3162270 -4 Hafnia alvei	99		3166070 -4 Hafnia alvei	99	
3156777 Klebsiella oxytoca	99		3164006 -7 Serratia marcescens	99		3166102 -7 Serratia marcescens	99	
3160010 -4 Salmonella choleraesuis	99	SEROLOGY	3164014 Salmonella choleraesuis	99	SEROLOGY	3166240 -4 Hafnia alvei	99	
3160040 -4 Hafnia alvei	99		3164040 -4 Hafnia alvei	99		3166250 -4 Hafnia alvei	99	
3160050 -4 Hafnia alvei	99		3164050 -4 Hafnia alvei	99		3166260 Hafnia alvei	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3244274 -4 Citrobacter freundii	99		3300004 Salmonella typhi	99	SEROLOGY	3324050 Hafnia alvei	99	
3251000 Proteus mirabilis	99		3300040 Salmonella gallinarum	99	SEROLOGY	3324054 Salmonella subgroup 4 Salmonella subgroup 1	65	SEROLOGY 34 SEROLOGY
3251500 Proteus vulgaris	99		3300050 Salmonella gallinarum	99	SEROLOGY	3324056 Salmonella subgroup 1	99	SEROLOGY
3253500 Proteus vulgaris	99		3300060 Salmonella gallinarum	99	SEROLOGY	3324070 Hafnia alvei	99	
3255000 Proteus mirabilis	99		3311000 Proteus mirabilis	99		3324250 Hafnia alvei	99	
3255500 Proteus vulgaris	99		3311500 Proteus vulgaris	99		3324364 -6 Serratia liquefaciens	99	
3265000 Proteus mirabilis	99		3313500 Proteus vulgaris	99		3325000 Proteus mirabilis	99	
3270000 Proteus mirabilis	99		3315000 Proteus mirabilis	99		3326050 Hafnia alvei	99	
3271000 Proteus mirabilis	99		3315500 Proteus vulgaris	99		3326054 Salmonella subgroup 2	99	SEROLOGY
3273000 Proteus mirabilis	99		3320014 Salmonella choleraesuis	99	SEROLOGY	3326070 Hafnia alvei	99	
3274000 Proteus mirabilis	99		3320050 Salmonella pullorum Hafnia alvei	95	SEROLOGY 4	3326250 Hafnia alvei	99	
3275000 -4 Proteus mirabilis	99		3320054 Salmonella pullorum	99	SEROLOGY	3330000 Proteus mirabilis	99	
3275010 Proteus mirabilis	99		3320070 Hafnia alvei	99		3331000 Proteus mirabilis	99	
3275020 Proteus mirabilis	99		3320250 Hafnia alvei	99		3331400 Morgan. morganii biogr.1	99	
3275100 Proteus mirabilis	99		3320400 Edwardsiella tarda	99		3333000 Proteus mirabilis	99	
3275200 Proteus mirabilis	99		3322050 Hafnia alvei	99		3334000 Proteus mirabilis	99	
3275400 Proteus mirabilis	99		3322070 Hafnia alvei	99		3335000 -4 Proteus mirabilis	99	
3277000 Proteus mirabilis	99		3322250 Hafnia alvei	99		3335010 Proteus mirabilis	99	
3300000 Salmonella gallinarum	99	SEROLOGY	3324014 Salmonella choleraesuis	99	SEROLOGY	3335020 Proteus mirabilis	99	
3244274						3335020		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3400075 Enterobacter agglomerans	91		3400177 Enterobacter agglomerans	94		3400272 Enterobacter agglomerans	96	
Klebsiella ozaenae	8		Klebsiella ozaenae	5		Klebsiella ozaenae	3	
3400076 Enterobacter agglomerans	95		3400200 Ewingella americana	99		3400273 Enterobacter agglomerans	68	MOTL/ 85
Klebsiella ozaenae	4					Klebsiella ozaenae	31	MOTL/ 0
3400077 Enterobacter agglomerans	63	MOTL/ 85	3400201 -7 Enterobacter agglomerans	99		3400274 Enterobacter agglomerans	99	
Klebsiella ozaenae	36	MOTL/ 0						
3400100 -5 Enterobacter agglomerans	99		3400210 Ewingella americana	99		3400275 Enterobacter agglomerans	96	
						Klebsiella ozaenae	3	
3400106 Serratia marcescens	82		3400211 -7 Enterobacter agglomerans	99		3400276 Enterobacter agglomerans	97	
Enterobacter agglomerans	17					Klebsiella ozaenae	2	
3400107 Serratia marcescens	92		3400220 -7 Enterobacter agglomerans	99		3400277 Enterobacter agglomerans	81	
Enterobacter agglomerans	7					Klebsiella ozaenae	18	
3400110 -7 Enterobacter agglomerans	99		3400230 -7 Enterobacter agglomerans	99		3400278 Enterobacter agglomerans	99	
3400120 -7 Enterobacter agglomerans	99		3400240 Enterobacter agglomerans	62	MOTL/ 85	3400300 -7 Enterobacter agglomerans	99	
			Shigella species	37	MOTL/ 0			
3400130 -7 Enterobacter agglomerans	99		3400241 -7 Enterobacter agglomerans	99		3400310 -7 Enterobacter agglomerans	99	
3400140 Enterobacter agglomerans	62	MOTL/ 85	3400250 -7 Enterobacter agglomerans	99		3400320 -7 Enterobacter agglomerans	99	
Shigella species	37	MOTL/ 0						
3400141 -7 Enterobacter agglomerans	99		3400260 Enterobacter agglomerans	69	MOTL/ 85	3400330 -7 Enterobacter agglomerans	99	
			Shigella species	30	MOTL/ 0			
3400150 -7 Enterobacter agglomerans	99		3400261 -2 Enterobacter agglomerans	99		3400340 -7 Enterobacter agglomerans	99	
3400160 Enterobacter agglomerans	69	MOTL/ 85	3400263 Enterobacter agglomerans	79	MOTL/ 85	3400350 -7 Enterobacter agglomerans	99	
Shigella species	30	MOTL/ 0	Klebsiella ozaenae	20	MOTL/ 0			
3400161 -7 Enterobacter agglomerans	99		3400264 -6 Enterobacter agglomerans	99		3400360 Enterobacter agglomerans	62	GELT/ 2
						Serratia rubidaea	37	GELT/ 90
3400170 -2 Enterobacter agglomerans	99		3400267 Enterobacter agglomerans	88		3400361 Serratia rubidaea	75	GELT/ 90
			Klebsiella ozaenae	11		Enterobacter agglomerans	24	GELT/ 2
3400173 Enterobacter agglomerans	89		3400270 Enterobacter agglomerans	99		3400362 Enterobacter agglomerans	99	
Klebsiella ozaenae	10							
3400174 -6 Enterobacter agglomerans	99		3400271 Enterobacter agglomerans	93		3400363 Serratia rubidaea	65	GELT/ 90
			Klebsiella ozaenae	6		Enterobacter agglomerans	34	GELT/ 2

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3402130 -7 Enterobacter agglomerans	99		3402300 -7 Enterobacter agglomerans	99		3402376 Enterobacter agglomerans	98	
3402140 -7 Enterobacter agglomerans	99		3402310 -7 Enterobacter agglomerans	99		Enterobacter cloacae	1	
3402150 -7 Enterobacter agglomerans	99		3402320 -7 Enterobacter agglomerans	99		3402377 Enterobacter agglomerans	56	MOTL/ 85
3402160 -7 Enterobacter agglomerans	99		3402330 -7 Enterobacter agglomerans	99		Klebsiella pneumoniae	43	MOTL/ 0
3402170 -3 Enterobacter agglomerans	99		3402340 -7 Enterobacter agglomerans	99		3402400 -6 Enterobacter agglomerans	99	
3402174 Enterobacter agglomerans	95		3402340 -7 Enterobacter agglomerans	99		3402410 -7 Enterobacter agglomerans	99	
Enterobacter cloacae	4		3402350 -7 Enterobacter agglomerans	99		3402420 -6 Enterobacter agglomerans	99	
3402175 Enterobacter agglomerans	95		3402360 Serratia rubidaea	60	GELT/ 90	3402430 -7 Enterobacter agglomerans	99	
Enterobacter cloacae	4		Enterobacter agglomerans	39	GELT/ 2	3402440 -7 Enterobacter agglomerans	99	
3402176 -7 Enterobacter agglomerans	99		3402361 Serratia rubidaea	88		3402450 -7 Enterobacter agglomerans	99	
3402200 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	11		3402460 -7 Enterobacter agglomerans	99	
3402210 -7 Enterobacter agglomerans	99		3402362 Enterobacter agglomerans	51	GELT/ 2	3402470 -7 Enterobacter agglomerans	99	
3402220 -7 Enterobacter agglomerans	99		Serratia rubidaea	48	GELT/ 90	3402500 -7 Enterobacter agglomerans	99	
3402230 -7 Enterobacter agglomerans	99		3402363 Serratia rubidaea	82		3402510 -7 Enterobacter agglomerans	99	
3402240 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	17		3402520 -7 Enterobacter agglomerans	99	
3402250 -7 Enterobacter agglomerans	99		3402364 Enterobacter agglomerans	99		3402530 -7 Enterobacter agglomerans	99	
3402260 -7 Enterobacter agglomerans	99		3402365 Enterobacter agglomerans	78	GELT/ 2	3402540 -7 Enterobacter agglomerans	99	
3402270 -2 Enterobacter agglomerans	99		Serratia rubidaea	21	GELT/ 90	3402550 -7 Enterobacter agglomerans	99	
3402273 Enterobacter agglomerans	97		3402366 -7 Enterobacter agglomerans	99		3402560 -7 Enterobacter agglomerans	99	
Klebsiella ozaenae	2		3402370 -3 Enterobacter agglomerans	99		3402570 -7 Enterobacter agglomerans	99	
3402274 -7 Enterobacter agglomerans	99		3402374 Enterobacter agglomerans	95				
			Enterobacter cloacae	4				
			3402375 Enterobacter agglomerans	95				
			Enterobacter cloacae	4				

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3402600 -7 Enterobacter agglomerans	99		3403050 -7 Enterobacter agglomerans	99		3403340 -7 Enterobacter agglomerans	99	
3402610 -7 Enterobacter agglomerans	99		3403060 -6 Enterobacter agglomerans	99		3403350 -7 Enterobacter agglomerans	99	
3402620 -7 Enterobacter agglomerans	99		3403070 -7 Enterobacter agglomerans	99		3403360 -7 Enterobacter agglomerans	99	
3402630 -7 Enterobacter agglomerans	99		3403110 -4 Enterobacter agglomerans	99		3403370 -7 Enterobacter agglomerans	99	
3402640 -7 Enterobacter agglomerans	99		3403130 -7 Enterobacter agglomerans	99		3403450 -6 Enterobacter agglomerans	99	
3402650 -7 Enterobacter agglomerans	99		3403140 -7 Enterobacter agglomerans	99		3403470 -4 Enterobacter agglomerans	99	
3402660 -7 Enterobacter agglomerans	99		3403150 -7 Enterobacter agglomerans	99		3403550 -7 Enterobacter agglomerans	99	
3402670 -7 Enterobacter agglomerans	99		3403160 -7 Enterobacter agglomerans	99		3403560 -6 Enterobacter agglomerans	99	
3402700 -7 Enterobacter agglomerans	99		3403170 -7 Enterobacter agglomerans	99		3403570 -7 Enterobacter agglomerans	99	
3402710 -7 Enterobacter agglomerans	99		3403210 -4 Enterobacter agglomerans	99		3403634 -4 Enterobacter agglomerans	99	
3402720 -7 Enterobacter agglomerans	99		3403230 -7 Enterobacter agglomerans	99		3403650 -7 Enterobacter agglomerans	99	
3402730 -7 Enterobacter agglomerans	99		3403240 -7 Enterobacter agglomerans	99		3403660 -6 Enterobacter agglomerans	99	
3402740 -7 Enterobacter agglomerans	99		3403250 -7 Enterobacter agglomerans	99		3403670 -7 Enterobacter agglomerans	99	
3402750 -7 Enterobacter agglomerans	99		3403260 -7 Enterobacter agglomerans	99		3403710 -6 Enterobacter agglomerans	99	
3402760 -7 Enterobacter agglomerans	99		3403270 -7 Enterobacter agglomerans	99		3403730 -6 Enterobacter agglomerans	99	
3402770 -6 Enterobacter agglomerans	99		3403300 -6 Enterobacter agglomerans	99		3403740 -7 Enterobacter agglomerans	99	
3402777 Klebsiella oxytoca	71	MOTL/ 0	3403310 -7 Enterobacter agglomerans	99		3403750 -7 Enterobacter agglomerans	99	
Enterobacter agglomerans	28	MOTL/ 85	3403320 -6 Enterobacter agglomerans	99		3403760 -7 Enterobacter agglomerans	99	
3403014 -4 Enterobacter agglomerans	99		3403330 -7 Enterobacter agglomerans	99		3403770 -7 Enterobacter agglomerans	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3404000 Ewingella americana	99		3404030 Enterobacter agglomerans	58	XYLO/ 93	3404076 Enterobacter agglomerans	98	
3404001 Ewingella americana	95		Ewingella americana	41	XYLO/ 13	3404077 Enterobacter agglomerans	92	
Enterobacter agglomerans	4		3404031 -7 Enterobacter agglomerans	99		Klebsiella ozaenae	7	
3404002 Ewingella americana	95		3404040 Ewingella americana	66	XYLO/ 13	3404100 Ewingella americana	75	XYLO/ 13
Enterobacter agglomerans	4		Enterobacter agglomerans	33	XYLO/ 93	Enterobacter agglomerans	15	XYLO/ 93
3404003 Enterobacter agglomerans	99		3404041 -7 Enterobacter agglomerans	99		Serratia marcescens	10	
3404004 Ewingella americana	88		3404050 Enterobacter agglomerans	85		3404101 Serratia marcescens	62	GELT/ 90
Enterobacter agglomerans	11		Ewingella americana	14		Enterobacter agglomerans	37	GELT/ 2
3404005 Serratia marcescens	53	GELT/ 90	3404051 -7 Enterobacter agglomerans	99		3404102 Serratia marcescens	85	
Enterobacter agglomerans	46	GELT/ 2	3404060 -2 Enterobacter agglomerans	99		Enterobacter agglomerans	14	
3404006 Serratia marcescens	79	GELT/ 90	3404063 Enterobacter agglomerans	91		3404103 Serratia marcescens	93	
Enterobacter agglomerans	20	GELT/ 2	Klebsiella ozaenae	8		Enterobacter agglomerans	6	
3404007 Serratia marcescens	90		3404064 -6 Enterobacter agglomerans	99		3404104 Serratia marcescens	71	GELT/ 90
Enterobacter agglomerans	9		3404067 Enterobacter agglomerans	95		Enterobacter agglomerans	28	GELT/ 2
3404010 Ewingella americana	99		Klebsiella ozaenae	4		3404105 Serratia marcescens	85	
3404011 Ewingella americana	63	XYLO/ 13	3404070 Enterobacter agglomerans	99		Enterobacter agglomerans	14	
Enterobacter agglomerans	36	XYLO/ 93	3404071 Enterobacter agglomerans	97		3404106 Serratia marcescens	95	
3404012 Ewingella americana	63	XYLO/ 13	Klebsiella ozaenae	2		Enterobacter agglomerans	4	
Enterobacter agglomerans	36	XYLO/ 93	3404072 Enterobacter agglomerans	98		3404107 Serratia marcescens	97	
3404013 Enterobacter agglomerans	99		Klebsiella ozaenae	1		Enterobacter agglomerans	2	
3404014 Enterobacter agglomerans	62	XYLO/ 93	3404073 Enterobacter agglomerans	86		3404110 Enterobacter agglomerans	70	XYLO/ 93
Ewingella americana	37	XYLO/ 13	Klebsiella ozaenae	13		Ewingella americana	29	XYLO/ 13
3404015 -7 Enterobacter agglomerans	99		3404074 Enterobacter agglomerans	99		3404111 -7 Enterobacter agglomerans	99	
3404020 Ewingella americana	89		3404075 Enterobacter agglomerans	98		3404120 -7 Enterobacter agglomerans	99	
Enterobacter agglomerans	10		Klebsiella ozaenae	1		3404130 -7 Enterobacter agglomerans	99	
3404021 -7 Enterobacter agglomerans	99		3404076 Enterobacter agglomerans	98		3404140 -7 Enterobacter agglomerans	99	

BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %	BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %	BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %
3406006					
Enterobacter agglomerans	68 GELT/ 2	3406104		3406162	
Serratia marcescens	31 GELT/ 90	Enterobacter agglomerans	77 GELT/ 2	Enterobacter agglomerans	76 GELT/ 2
		Serratia marcescens	22 GELT/ 90	Serratia rubidaea	23 GELT/ 90
3406007					
Serratia marcescens	53 GELT/ 90	3406105		3406163	
Enterobacter agglomerans	46 GELT/ 2	Enterobacter agglomerans	58 GELT/ 2	Serratia rubidaea	62 GELT/ 90
		Serratia marcescens	41 GELT/ 90	Enterobacter agglomerans	37 GELT/ 2
3406010					
Enterobacter agglomerans	69 XYLO/ 93	3406106		3406164	
Ewingella americana	30 XYLO/ 13	Serratia marcescens	71 GELT/ 90	Enterobacter agglomerans	91
		Enterobacter agglomerans	28 GELT/ 2	Enterobacter cloacae	8
3406011 -7					
Enterobacter agglomerans	99	3406107		3406165	
		Serratia marcescens	85	Enterobacter agglomerans	83
3406020 -7		Enterobacter agglomerans	14	Enterobacter cloacae	8
Enterobacter agglomerans	99			Serratia rubidaea	8
		3406110 -7			
3406030 -7		Enterobacter agglomerans	99	3406166	
Enterobacter agglomerans	99			Enterobacter agglomerans	96
		3406120 -7		Enterobacter cloacae	3
3406040 -7		Enterobacter agglomerans	99		
Enterobacter agglomerans	99			3406167	
		3406130 -3		Enterobacter agglomerans	95
3406050 -7		Enterobacter agglomerans	99	Enterobacter cloacae	4
Enterobacter agglomerans	99				
		3406134		3406170	
3406060 -7		Enterobacter agglomerans	97	Enterobacter agglomerans	96
Enterobacter agglomerans	99	Enterobacter cloacae	2	Enterobacter cloacae	3
3406070 -2		3406135 -7		3406171	
Enterobacter agglomerans	99	Enterobacter agglomerans	99	Enterobacter agglomerans	96
				Enterobacter cloacae	3
3406073					
Enterobacter agglomerans	98	3406140 -7		3406172	
		Enterobacter agglomerans	99	Enterobacter agglomerans	98
				Enterobacter cloacae	1
3406074		3406150 -3			
Enterobacter agglomerans	97	Enterobacter agglomerans	99	3406173	
Enterobacter cloacae	2			Enterobacter agglomerans	99
		3406154			
3406075		Enterobacter agglomerans	98	3406174	
Enterobacter agglomerans	97	Enterobacter cloacae	1	Enterobacter agglomerans	84
Enterobacter cloacae	2			Enterobacter cloacae	15
		3406155			
3406076 -7		Enterobacter agglomerans	98	3406175	
Enterobacter agglomerans	99	Enterobacter cloacae	1	Enterobacter agglomerans	83
				Enterobacter cloacae	16
3406100 -1		3406156 -7			
Enterobacter agglomerans	99	Enterobacter agglomerans	99	3406176	
				Enterobacter agglomerans	93
3406102		3406160		Enterobacter cloacae	6
Enterobacter agglomerans	59 GELT/ 2	Enterobacter agglomerans	65 GELT/ 2		
Serratia marcescens	40 GELT/ 90	Serratia rubidaea	34 GELT/ 90	3406177	
				Enterobacter agglomerans	72 HOTL/ 85
3406103		3406161		Klebsiella pneumoniae	18 HOTL/ 0
Serratia marcescens	61 GELT/ 90	Serratia rubidaea	73 GELT/ 90	Enterobacter cloacae	5
Enterobacter agglomerans	38 GELT/ 2	Enterobacter agglomerans	26 GELT/ 2		
3406006					3406177

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3406200								
Ewingella americana	83		3406276			3406340		
Enterobacter agglomerans	16		Enterobacter agglomerans	98		Enterobacter agglomerans	76	GELT/ 2
						Serratia rubidaea	23	GELT/ 90
3406201 -7			3406277					
Enterobacter agglomerans	99		Enterobacter agglomerans	79	MOTL/ 85	3406341		
			Klebsiella pneumoniae	20	MOTL/ 0	Serratia rubidaea	61	GELT/ 90
3406210						Enterobacter agglomerans	38	GELT/ 2
Enterobacter agglomerans	70	XYLO/ 93	3406300 -7					
Ewingella americana	29	XYLO/ 13	Enterobacter agglomerans	99		3406342		
						Enterobacter agglomerans	84	
3406211 -7			3406310 -7			Serratia rubidaea	15	
Enterobacter agglomerans	99		Enterobacter agglomerans	99				
3406220 -7			3406320			3406343		
Enterobacter agglomerans	99		Serratia rubidaea	56	GELT/ 90	Enterobacter agglomerans	51	GELT/ 2
			Enterobacter agglomerans	43	GELT/ 2	Serratia rubidaea	48	GELT/ 90
3406230 -7								
Enterobacter agglomerans	99		3406321			3406344		
			Serratia rubidaea	87		Enterobacter agglomerans	98	
3406240 -7			Enterobacter agglomerans	12				
Enterobacter agglomerans	99					3406345		
						Enterobacter agglomerans	94	
3406250 -7			3406322			Serratia rubidaea	5	
Enterobacter agglomerans	99		Enterobacter agglomerans	56	GELT/ 2			
			Serratia rubidaea	43	GELT/ 90	3406346 -7		
3406260						Enterobacter agglomerans	99	
Enterobacter agglomerans	94		3406323					
Serratia rubidaea	5		Serratia rubidaea	80	GELT/ 90	3406350 -3		
			Enterobacter agglomerans	19	GELT/ 2	Enterobacter agglomerans	99	
3406261								
Enterobacter agglomerans	78	GELT/ 2	3406324			3406354		
Serratia rubidaea	21	GELT/ 90	Enterobacter agglomerans	99		Enterobacter agglomerans	98	
						Enterobacter cloacae	1	
3406262			3406325					
Enterobacter agglomerans	99		Enterobacter agglomerans	81		3406355		
			Serratia rubidaea	18		Enterobacter agglomerans	98	
3406263						Enterobacter cloacae	1	
Enterobacter agglomerans	85		3406326 -7					
Serratia rubidaea	14		Enterobacter agglomerans	99		3406356		
						Enterobacter agglomerans	99	
3406264 -7			3406330 -3					
Enterobacter agglomerans	99		Enterobacter agglomerans	99		3406357		
						Enterobacter agglomerans	92	
3406270 -2			3406334			Klebsiella pneumoniae	6	
Enterobacter agglomerans	99		Enterobacter agglomerans	97				
			Enterobacter cloacae	2		3406360		
3406273						Serratia rubidaea	94	
Enterobacter agglomerans	99		3406335			Enterobacter agglomerans	4	
			Enterobacter agglomerans	97				
3406274			Enterobacter cloacae	2		3406361		
Enterobacter agglomerans	97					Serratia rubidaea	98	
Enterobacter cloacae	2		3406336					
			Enterobacter agglomerans	99		3406362		
3406275						Serratia rubidaea	92	
Enterobacter agglomerans	97		3406337			Enterobacter agglomerans	7	
Enterobacter cloacae	2		Enterobacter agglomerans	75	MOTL/ 85			
			Klebsiella pneumoniae	24	MOTL/ 0			

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3406363			3406377			3406570 -6		
Serratia rubidaea	98		Klebsiella pneumoniae	80	MOTL/ 0	Enterobacter agglomerans	99	
Enterobacter agglomerans	1		Enterobacter agglomerans	16	MOTL/ 85			
3406364			Enterobacter aerogenes	2		3406577		
Enterobacter agglomerans	56	GELT/ 2				Enterobacter agglomerans	54	MOTL/ 85
Serratia rubidaea	38	GELT/ 90	3406400 -7			Klebsiella oxytoca	45	MOTL/ 0
Enterobacter cloacae	5		Enterobacter agglomerans	99				
3406365			3406410 -7			3406600 -7		
Serratia rubidaea	76	GELT/ 90	Enterobacter agglomerans	99		Enterobacter agglomerans	99	
Enterobacter agglomerans	21	GELT/ 2				3406610 -7		
Enterobacter cloacae	2		3406420 -7			Enterobacter agglomerans	99	
3406366			Enterobacter agglomerans	99		3406620 -7		
Enterobacter agglomerans	69	GELT/ 2	3406430 -7			Enterobacter agglomerans	99	
Serratia rubidaea	28	GELT/ 90	Enterobacter agglomerans	99		3406630 -7		
Enterobacter cloacae	2					Enterobacter agglomerans	99	
3406367			3406440 -7			3406640 -7		
Serratia rubidaea	63	GELT/ 90	Enterobacter agglomerans	99		Enterobacter agglomerans	99	
Enterobacter agglomerans	29	GELT/ 2				3406650 -7		
Klebsiella pneumoniae	6		3406450 -4			Enterobacter agglomerans	99	
3406370			Enterobacter agglomerans	99		3406655 -7		
Enterobacter agglomerans	91		3406455			Enterobacter agglomerans	99	
Serratia rubidaea	5		Enterobacter agglomerans	92		3406660 -7		
Enterobacter cloacae	3		Citrobacter diversus	7		Enterobacter agglomerans	99	
3406371			3406456 -7			3406670 -6		
Enterobacter agglomerans	74	GELT/ 2	Enterobacter agglomerans	99		Enterobacter agglomerans	99	
Serratia rubidaea	23	GELT/ 90				3406677		
Enterobacter cloacae	2		3406460 -7			Enterobacter agglomerans	54	MOTL/ 85
3406372			Enterobacter agglomerans	99		Klebsiella oxytoca	45	MOTL/ 0
Enterobacter agglomerans	95		3406470 -7					
Serratia rubidaea	3		Enterobacter agglomerans	99		3406700 -7		
Enterobacter cloacae	1					Enterobacter agglomerans	99	
3406373			3406500 -7			3406710 -7		
Enterobacter agglomerans	73	GELT/ 2	Enterobacter agglomerans	99		Enterobacter agglomerans	99	
Serratia rubidaea	13	GELT/ 90				3406720 -7		
Klebsiella pneumoniae	12		3406510 -7			Enterobacter agglomerans	99	
3406374			Enterobacter agglomerans	99		3406730 -6		
Enterobacter agglomerans	84		3406520 -7			Enterobacter agglomerans	99	
Enterobacter cloacae	14		Enterobacter agglomerans	99		3406737		
3406375			3406530 -7			Klebsiella oxytoca	51	MOTL/ 0
Enterobacter agglomerans	75		Enterobacter agglomerans	99		Enterobacter agglomerans	48	MOTL/ 85
Enterobacter cloacae	13					3406740 -7		
Klebsiella pneumoniae	9		3406540 -7			Enterobacter agglomerans	99	
3406376			Enterobacter agglomerans	99		3406750 -6		
Enterobacter agglomerans	85		3406550 -7			Enterobacter agglomerans	99	
Klebsiella pneumoniae	7		Enterobacter agglomerans	99				
Enterobacter cloacae	6		3406560 -7			Enterobacter agglomerans	99	
			Enterobacter agglomerans	99				

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3406757								
Enterobacter agglomerans	80		3407010 -7			3407240 -7		
Klebsiella oxytoca	19		Enterobacter agglomerans	99		Enterobacter agglomerans	99	
3406760			3407020 -6			3407250 -7		
Enterobacter agglomerans	58 GELT/ 2		Enterobacter agglomerans	99		Enterobacter agglomerans	99	
Serratia rubidaea	41 GELT/ 90		3407030 -7			3407260 -7		
3406761			Enterobacter agglomerans	99		Enterobacter agglomerans	99	
Serratia rubidaea	78 GELT/ 90		3407040 -7			3407270 -7		
Enterobacter agglomerans	21 GELT/ 2		Enterobacter agglomerans	99		Enterobacter agglomerans	99	
3406762			3407050 -7			3407300 -7		
Enterobacter agglomerans	70 GELT/ 2		Enterobacter agglomerans	99		Enterobacter agglomerans	99	
Serratia rubidaea	29 GELT/ 90		3407060 -7			3407310 -7		
3406763			Enterobacter agglomerans	99		Enterobacter agglomerans	99	
Serratia rubidaea	69 GELT/ 90		3407070 -7			3407320 -7		
Enterobacter agglomerans	30 GELT/ 2		Enterobacter agglomerans	99		Enterobacter agglomerans	99	
3406764			3407100 -7			3407330 -7		
Enterobacter agglomerans	99		Enterobacter agglomerans	99		Enterobacter agglomerans	99	
3406765			3407110 -7			3407340 -7		
Enterobacter agglomerans	88		Enterobacter agglomerans	99		Enterobacter agglomerans	99	
Serratia rubidaea	11		3407120 -7			3407350 -7		
3406766			Enterobacter agglomerans	99		Enterobacter agglomerans	99	
Enterobacter agglomerans	99		3407130 -7			3407360		
3406767			Enterobacter agglomerans	99		Enterobacter agglomerans	66 GELT/ 2	
Enterobacter agglomerans	60 MOTL/ 85		3407140 -7			Serratia rubidaea	33 GELT/ 90	
Klebsiella oxytoca	39 MOTL/ 0		Enterobacter agglomerans	99		3407361		
3406770 -2			3407150 -7			Serratia rubidaea	72 GELT/ 90	
Enterobacter agglomerans	99		Enterobacter agglomerans	99		Enterobacter agglomerans	27 GELT/ 2	
3406773			3407160 -7			3407362		
Enterobacter agglomerans	63 MOTL/ 85		Enterobacter agglomerans	99		Enterobacter agglomerans	99	
Klebsiella oxytoca	36 MOTL/ 0		3407170 -7			3407363		
3406774			Enterobacter agglomerans	99		Serratia rubidaea	61 GELT/ 90	
Enterobacter agglomerans	99		3407200 -7			Enterobacter agglomerans	38 GELT/ 2	
3406775			Enterobacter agglomerans	99		3407364 -7		
Enterobacter agglomerans	71 MOTL/ 85		3407210 -7			Enterobacter agglomerans	99	
Klebsiella oxytoca	28 MOTL/ 0		Enterobacter agglomerans	99		3407370 -7		
3406776			3407220 -7			Enterobacter agglomerans	99	
Enterobacter agglomerans	77 MOTL/ 85		Enterobacter agglomerans	99		3407410 -6		
Klebsiella oxytoca	22 MOTL/ 0		3407230 -7			Enterobacter agglomerans	99	
3406777			Enterobacter agglomerans	99		3407430 -7		
Klebsiella oxytoca	93		3407000 -6			Enterobacter agglomerans	99	
Enterobacter agglomerans	6		Enterobacter agglomerans	99				
3407000 -6								
Enterobacter agglomerans	99							
3406757								

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3413550 -6 Enterobacter agglomerans	99		3414050 Enterobacter agglomerans	82		3414130 -7 Enterobacter agglomerans	99	
3413570 -7 Enterobacter agglomerans	99		Yersinia pseudotubercul.	17		3414140 -7 Enterobacter agglomerans	99	
3413650 -6 Enterobacter agglomerans	99		3414051 -7 Enterobacter agglomerans	99		3414150 -7 Enterobacter agglomerans	99	
3413670 -7 Enterobacter agglomerans	99		3414060 -7 Enterobacter agglomerans	99		3414160 -7 Enterobacter agglomerans	99	
3413734 -4 Enterobacter agglomerans	99		3414070 -2 Enterobacter agglomerans	99		3414170 -3 Enterobacter agglomerans	99	
3413750 -7 Enterobacter agglomerans	99		3414073 Enterobacter agglomerans	96		3414174 Enterobacter agglomerans	95	
3413760 -6 Enterobacter agglomerans	99		Klebsiella ozaenae	3		Enterobacter cloacae	4	
3413770 -7 Enterobacter agglomerans	99		3414074 -7 Enterobacter agglomerans	99		3414175 Enterobacter agglomerans	95	
3414000 Ewingella americana	95		3414100 Enterobacter agglomerans	99		Enterobacter cloacae	4	
Enterobacter agglomerans	4		3414101 Serratia marcescens	60	GELT/ 90	3414176 Enterobacter agglomerans	99	
3414001 -5 Enterobacter agglomerans	99		Enterobacter agglomerans	40	GELT/ 2	3414177 Enterobacter agglomerans	62	MOTL/ 85
3414006 Serratia marcescens	78	GELT/ 90	3414102 Serratia marcescens	83		Klebsiella pneumoniae	37	MOTL/ 0
Enterobacter agglomerans	21	GELT/ 2	Enterobacter agglomerans	16		3414200 Ewingella americana	95	
3414007 Serratia marcescens	89		3414103 Serratia marcescens	92		Enterobacter agglomerans	4	
Enterobacter agglomerans	10		Enterobacter agglomerans	7		3414201 -7 Enterobacter agglomerans	99	
3414010 Ewingella americana	48	AROL/ 99	3414104 Serratia marcescens	69	GELT/ 90	3414210 Ewingella americana	66	XYLO/ 13
Yersinia pseudotubercul.	28	AROL/ 0	Enterobacter agglomerans	30	GELT/ 2	Enterobacter agglomerans	33	XYLO/ 93
Enterobacter agglomerans	23		3414105 Serratia marcescens	84		3414211 -7 Enterobacter agglomerans	99	
3414011 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	15		3414220 -7 Enterobacter agglomerans	99	
3414020 -7 Enterobacter agglomerans	99		3414106 Serratia marcescens	94		3414230 -7 Enterobacter agglomerans	99	
3414030 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	5		3414240 -7 Enterobacter agglomerans	99	
3414040 -7 Enterobacter agglomerans	99		3414107 Serratia marcescens	97		3414250 -7 Enterobacter agglomerans	99	
			Enterobacter agglomerans	2				
			3414110 -7 Enterobacter agglomerans	99				
			Enterobacter agglomerans	99				
			3414120 -7 Enterobacter agglomerans	99				

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3414260 -7 Enterobacter agglomerans	99		3414363 Serratia rubidaea	73	GELT/ 90	3414450 -7 Enterobacter agglomerans	99	
3414270 -2 Enterobacter agglomerans	99		Enterobacter agglomerans	26	GELT/ 2	3414460 -7 Enterobacter agglomerans	99	
3414273 Enterobacter agglomerans	98		3414364 Enterobacter agglomerans	99		3414470 -7 Enterobacter agglomerans	99	
Klebsiella ozaenae	1		3414365 Enterobacter agglomerans	86		3414500 Proteus vulgaris	76	GLT/ 91
3414274 -6 Enterobacter agglomerans	99		Serratia rubidaea	13		Enterobacter agglomerans	23	GELT/ 2
3414277 Enterobacter agglomerans	62	MOTL/ 85	3414366 Enterobacter agglomerans	99		3414501 -7 Enterobacter agglomerans	99	
Klebsiella pneumoniae	37	MOTL/ 0	3414367 Enterobacter agglomerans	68	MOTL/ 85	3414510 -7 Enterobacter agglomerans	99	
3414300 -6 Enterobacter agglomerans	99		Klebsiella pneumoniae	31	MOTL/ 0	3414520 -7 Enterobacter agglomerans	99	
3414307 Enterobacter agglomerans	61	GELT/ 2	3414370 -2 Enterobacter agglomerans	99		3414530 -7 Enterobacter agglomerans	99	
Serratia marcescens	38	GELT/ 90	3414373 Enterobacter agglomerans	71	MOTL/ 85	3414540 -7 Enterobacter agglomerans	99	
3414310 -7 Enterobacter agglomerans	99		Klebsiella pneumoniae	28	MOTL/ 0	3414550 -7 Enterobacter agglomerans	99	
3414320 -7 Enterobacter agglomerans	99		3414374 Enterobacter agglomerans	95		3414560 -7 Enterobacter agglomerans	99	
3414330 -6 Enterobacter agglomerans	99		Enterobacter cloacae	4		3414570 -6 Enterobacter agglomerans	99	
3414337 Enterobacter agglomerans	57	MOTL/ 85	3414375 Enterobacter agglomerans	75	MOTL/ 85	3414577 Klebsiella oxytoca	66	MOTL/ 0
Klebsiella pneumoniae	42	MOTL/ 0	Klebsiella pneumoniae	21	MOTL/ 0	Enterobacter agglomerans	33	MOTL/ 85
3414340 -7 Enterobacter agglomerans	99		Enterobacter cloacae	3		3414600 -7 Enterobacter agglomerans	99	
3414350 -6 Enterobacter agglomerans	99		3414376 Enterobacter agglomerans	81		3414610 -7 Enterobacter agglomerans	99	
3414357 Enterobacter agglomerans	85		Klebsiella pneumoniae	17		3414620 -7 Enterobacter agglomerans	99	
Klebsiella pneumoniae	14		Enterobacter cloacae	1		3414630 -7 Enterobacter agglomerans	99	
3414360 Enterobacter agglomerans	52	GELT/ 2	3414377 Klebsiella pneumoniae	91		3414640 -7 Enterobacter agglomerans	99	
Serratia rubidaea	47	GELT/ 90	Enterobacter agglomerans	8		3414640 -7 Enterobacter agglomerans	99	
3414361 Serratia rubidaea	82		3414400 -6 Enterobacter agglomerans	99		3414660 -7 Enterobacter agglomerans	99	
Enterobacter agglomerans	17		3414410 -7 Enterobacter agglomerans	99		3414680 -7 Enterobacter agglomerans	99	
3414362 Enterobacter agglomerans	65	GELT/ 2	3414420 -7 Enterobacter agglomerans	99		3414690 -7 Enterobacter agglomerans	99	
Serratia rubidaea	34	GELT/ 90	3414430 -7 Enterobacter agglomerans	99		3414700 -7 Enterobacter agglomerans	99	
3414260			3414440 -7 Enterobacter agglomerans	99				

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3414650 -7 Enterobacter agglomerans	99		3414775 Enterobacter agglomerans	51	MOTL/ 85	3415170 -7 Enterobacter agglomerans	99	
3414660 -7 Enterobacter agglomerans	99		Klebsiella oxytoca	48	MOTL/ 0	3415204 Enterobacter agglomerans	99	
3414670 -6 Enterobacter agglomerans	99		3414776 Enterobacter agglomerans	59	MOTL/ 85	3415210 -7 Enterobacter agglomerans	99	
3414677 Klebsiella oxytoca	66	MOTL/ 0	Klebsiella oxytoca	97		3415224 Enterobacter agglomerans	99	
Enterobacter agglomerans	33	MOTL/ 85	Enterobacter agglomerans	2		3415230 -7 Enterobacter agglomerans	99	
3414700 -7 Enterobacter agglomerans	99		3415000 Proteus mirabilis	99		3415240 -7 Enterobacter agglomerans	99	
3414710 -7 Enterobacter agglomerans	99		3415010 -4 Enterobacter agglomerans	99		3415250 -7 Enterobacter agglomerans	99	
3414720 -7 Enterobacter agglomerans	99		3415030 -6 Enterobacter agglomerans	99		3415260 -7 Enterobacter agglomerans	99	
3414730 -6 Enterobacter agglomerans	99		3415040 -6 Enterobacter agglomerans	99		3415270 -7 Enterobacter agglomerans	99	
3414737 Klebsiella oxytoca	70	MOTL/ 0	3415050 -7 Enterobacter agglomerans	99		3415300 -6 Enterobacter agglomerans	99	
Enterobacter agglomerans	29	MOTL/ 85	3415060 -7 Enterobacter agglomerans	99		3415310 -7 Enterobacter agglomerans	99	
3414740 -7 Enterobacter agglomerans	99		3415070 -7 Enterobacter agglomerans	99		3415320 -7 Enterobacter agglomerans	99	
3414750 -6 Enterobacter agglomerans	99		3415100 Proteus vulgaris	99		3415330 -7 Enterobacter agglomerans	99	
3414757 Enterobacter agglomerans	64	MOTL/ 85	3415104 Enterobacter agglomerans	99		3415340 -7 Enterobacter agglomerans	99	
Klebsiella oxytoca	35	MOTL/ 0	3415110 -7 Enterobacter agglomerans	99		3415350 -7 Enterobacter agglomerans	99	
3414760 -6 Enterobacter agglomerans	99		3415124 Enterobacter agglomerans	99		3415360 -7 Enterobacter agglomerans	99	
3414767 Klebsiella oxytoca	59	MOTL/ 0	3415130 -7 Enterobacter agglomerans	99		3415370 -7 Enterobacter agglomerans	99	
Enterobacter agglomerans	40	MOTL/ 85	3415140 -7 Enterobacter agglomerans	99		3415400 Proteus vulgaris	99	
3414770 -2 Enterobacter agglomerans	99		3415150 -7 Enterobacter agglomerans	99		3415450 -4 Enterobacter agglomerans	99	
3414773 Klebsiella oxytoca	56	MOTL/ 0	3415160 -7 Enterobacter agglomerans	99				
Enterobacter agglomerans	43	MOTL/ 85						
3414774 Enterobacter agglomerans	99							

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3415470 -7 Enterobacter agglomerans	99		3415760 -7 Enterobacter agglomerans	99		3416110 -7 Enterobacter agglomerans	99	
3415500 -2 Proteus vulgaris	99		3415770 -7 Enterobacter agglomerans	99		3416120 -7 Enterobacter agglomerans	99	
3415510 Proteus vulgaris	99		3416000 -7 Enterobacter agglomerans	99		3416130 -7 Enterobacter agglomerans	99	
3415514 Enterobacter agglomerans	99		3416010 -7 Enterobacter agglomerans	99		3416140 -7 Enterobacter agglomerans	99	
3415530 -4 Enterobacter agglomerans	99		3416020 -7			3416150 -7		
3415540 -6 Enterobacter agglomerans	99		3416030 -7 Enterobacter agglomerans	99		3416160 Enterobacter agglomerans	99	
3415550 -7 Enterobacter agglomerans	99		3416040 -7 Enterobacter agglomerans	99		3416161 Enterobacter agglomerans	75 GELT/ 2	
3415560 -6 Enterobacter agglomerans	99		3416050 -7 Enterobacter agglomerans	99		Serratia rubidaea	24 GELT/ 90	
3415570 -7 Enterobacter agglomerans	99		3416060 -7 Enterobacter agglomerans	99		3416162 -3 Enterobacter agglomerans	99	
3415614 Enterobacter agglomerans	99		3416070 -7 Enterobacter agglomerans	99		3416164 Enterobacter agglomerans	95	
3415630 -4 Enterobacter agglomerans	99		3416100 -1 Enterobacter agglomerans	99		Enterobacter cloacae	4	
3415640 -6 Enterobacter agglomerans	99		3416102 Enterobacter agglomerans	62 GELT/ 2		3416165 -7 Enterobacter agglomerans	99	
3415650 -7 Enterobacter agglomerans	99		Serratia marcescens	37 GELT/ 90		3416170 Enterobacter agglomerans	89	
3415660 -6 Enterobacter agglomerans	99		3416103 Serratia marcescens	59 GELT/ 90		Enterobacter gergoviae	8	
3415670 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	40 GELT/ 2		Enterobacter cloacae	1	
3415710 -4 Enterobacter agglomerans	99		3416104 Enterobacter agglomerans	78 GELT/ 2		3416171 Enterobacter agglomerans	99	
3415730 -7 Enterobacter agglomerans	99		Serratia marcescens	21 GELT/ 90		3416172 Enterobacter agglomerans	92	
3415740 -7 Enterobacter agglomerans	99		3416105 Enterobacter agglomerans	60 GELT/ 2		Enterobacter gergoviae	7	
3415750 -7 Enterobacter agglomerans	99		Serratia marcescens	39 GELT/ 90		3416173 Enterobacter agglomerans	99	
			3416106 Serratia marcescens	69 GELT/ 90		3416174 Enterobacter agglomerans	92	
			Enterobacter agglomerans	30 GELT/ 2		Enterobacter cloacae	7	
			3416107 Serratia marcescens	84		3416175 Enterobacter agglomerans	88	
			Enterobacter agglomerans	15		Enterobacter cloacae	7	
						Klebsiella pneumoniae	5	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3417024 Enterobacter agglomerans	99		3417250 -7 Enterobacter agglomerans	99		3417530 -7 Enterobacter agglomerans	99	
3417030 -7 Enterobacter agglomerans	99		3417260 -7 Enterobacter agglomerans	99		3417540 -7 Enterobacter agglomerans	99	
3417040 -7 Enterobacter agglomerans	99		3417270 -7 Enterobacter agglomerans	99		3417550 -7 Enterobacter agglomerans	99	
3417050 -7 Enterobacter agglomerans	99		3417300 -7 Enterobacter agglomerans	99		3417560 -7 Enterobacter agglomerans	99	
3417060 -7 Enterobacter agglomerans	99		341710 -7 Enterobacter agglomerans	99		3417570 -7 Enterobacter agglomerans	99	
3417070 -7 Enterobacter agglomerans	99		3417320 -7 Enterobacter agglomerans	99		3417610 -4 Enterobacter agglomerans	99	
3417100 -6 Enterobacter agglomerans	99		3417330 -7 Enterobacter agglomerans	99		3417630 -7 Enterobacter agglomerans	99	
3417110 -7 Enterobacter agglomerans	99		3417340 -7 Enterobacter agglomerans	99		3417640 -7 Enterobacter agglomerans	99	
3417120 -7 Enterobacter agglomerans	99		3417350 -7 Enterobacter agglomerans	99		3417650 -7 Enterobacter agglomerans	99	
3417130 -7 Enterobacter agglomerans	99		3417360 -7 Enterobacter agglomerans	99		3417660 -7 Enterobacter agglomerans	99	
3417140 -7 Enterobacter agglomerans	99		3417370 -6 Enterobacter agglomerans	99		3417670 -7 Enterobacter agglomerans	99	
3417150 -7 Enterobacter agglomerans	99		3417377 Klebsiella pneumoniae	52	MOTL/ 0	3417700 -6 Enterobacter agglomerans	99	
3417160 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	47	MOTL/ 85	3417710 -7 Enterobacter agglomerans	99	
3417170 -7 Enterobacter agglomerans	99		3417414 -4 Enterobacter agglomerans	99		3417720 -6 Enterobacter agglomerans	99	
3417200 -6 Enterobacter agglomerans	99		3417440 -6 Enterobacter agglomerans	99		3417730 -7 Enterobacter agglomerans	99	
3417210 -7 Enterobacter agglomerans	99		3417450 -7 Enterobacter agglomerans	99		3417740 -7 Enterobacter agglomerans	99	
3417220 -7 Enterobacter agglomerans	99		3417460 -6 Enterobacter agglomerans	99		3417750 -7 Enterobacter agglomerans	99	
3417230 -7 Enterobacter agglomerans	99		3417470 -7 Enterobacter agglomerans	99		3417760 -7 Enterobacter agglomerans	99	
3417240 -7 Enterobacter agglomerans	99		3417500 Proteus vulgaris	99		3417770 -6 Enterobacter agglomerans	99	
			3417510 -4 Enterobacter agglomerans	99				
3417024						3417770		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3417777								
Klebsiella oxytoca	78	MOTL/ 0	3420250			3420376		
Enterobacter agglomerans	21	MOTL/ 85	Hafnia alvei	99		Serratia liquefaciens	93	
3420006 -7			3420260 -7			Enterobacter aerogenes	4	
Serratia marcescens	99		Serratia liquefaciens	99		Enterobacter cloacae	2	
3420040			3420270			3420377		
Shigella sonnei	83		Hafnia alvei	51	ESCL/ 7	Enterobacter aerogenes	88	
Shigella species	12		Serratia liquefaciens	48	ESCL/ 97	Serratia liquefaciens	9	
Hafnia alvei	4					Enterobacter cloacae	2	
3420050			3420272 -4			3420444		
Shigella sonnei	92		Serratia liquefaciens	99		Escherichia coli inaktiv	99	
Hafnia alvei	7							
3420060			3420306			3420454		
Shigella species	99		Serratia marcescens	53	XYLO/ 7	Citrobacter amalonaticus	89	
			Serratia liquefaciens	46	XYLO/100	Escherichia coli inaktiv	10	
3420064 -6			3420307			3420474		
Serratia liquefaciens	99		Serratia marcescens	99		Citrobacter amalonaticus	99	
3420070			3420320 -7			3420507		
Hafnia alvei	99		Serratia liquefaciens	99		Serratia marcescens	99	
3420074 -6			3420330 -6			3420544		
Serratia liquefaciens	99		Serratia liquefaciens	99		Yersinia enterocolitica	99	
3420100 -7			3420340 -7			3420554		
Serratia marcescens	99		Serratia liquefaciens	99		Citrobacter amalonaticus	99	
3420124 -6			3420350 -6			3420564 -6		
Serratia liquefaciens	99		Serratia liquefaciens	99		Serratia liquefaciens	99	
3420160 -7			3420357			3420644		
Serratia liquefaciens	99		Enterobacter aerogenes	99		Citrobacter amalonaticus	99	
3420170 -2			3420360 -7			3420650		
Serratia liquefaciens	99		Serratia liquefaciens	99		Citrobacter amalonaticus	99	
3420174			3420370 -2			3420654		
Serratia liquefaciens	79	CELL/ 5	Serratia liquefaciens	99		Citrobacter amalonaticus	92	
Enterobacter cloacae	20	CELL/ 99	3420373			Escherichia coli	7	
3420175			Enterobacter aerogenes	87		3420664 -4		
Enterobacter cloacae	73	CELL/ 99	Serratia liquefaciens	12		Citrobacter amalonaticus	99	
Serratia liquefaciens	26	CELL/ 5	3420374			3420754		
3420176			Serratia liquefaciens	93		Citrobacter amalonaticus	80	
Serratia liquefaciens	99		Enterobacter cloacae	6		Escherichia coli	19	
3420177			3420375			3420760 -6		
Enterobacter aerogenes	91		Enterobacter aerogenes	39	AROL/100	Serratia liquefaciens	99	
Serratia liquefaciens	8		Serratia liquefaciens	34	AROL/ 0			
3420224 -6			Enterobacter cloacae	25		3420774		
Serratia liquefaciens	99					Citrobacter amalonaticus	76	CELL/100
						Serratia liquefaciens	23	CELL/ 5

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3424027 Serratia marcescens	99		3424125 Serratia marcescens	81		3424156 Serratia liquefaciens		74 CELL/ 5
3424044 Serratia liquefaciens	99		Serratia liquefaciens	18		Enterobacter cloacae		25 CELL/ 99
3424046 Serratia marcescens	57	XYLO/ 7	3424126 Serratia liquefaciens	51	XYLO/100	3424157 Enterobacter aerogenes		99
Serratia liquefaciens	42	XYLO/100	Serratia marcescens	48	XYLO/ 7	3424160 Serratia liquefaciens		94
3424047 Serratia marcescens	99		3424127 Serratia marcescens	95		Enterobacter cloacae		5
3424050 Hafnia alvei	99		Serratia liquefaciens	4		3424161 -3 Serratia liquefaciens		99
3424060 -7 Serratia liquefaciens	99		3424134 Serratia liquefaciens	65	CELL/ 5	3424164 Serratia liquefaciens		95
3424070 Hafnia alvei	99		Enterobacter cloacae	34	CELL/ 99	Enterobacter cloacae		4
3424074 Enterobacter cloacae	51	CELL/ 99	3424135 Enterobacter cloacae	99		3424165 Serratia liquefaciens		65 CELL/ 5
Serratia liquefaciens	48	CELL/ 5	3424136 Serratia liquefaciens	99		Enterobacter cloacae		34 CELL/ 99
3424075 Enterobacter cloacae	99		3424140 Serratia liquefaciens	99		3424166 Serratia liquefaciens		98
3424076 Serratia liquefaciens	78	CELL/ 5	3424142 Serratia marcescens	63	XYLO/ 7	Enterobacter cloacae		1
Enterobacter cloacae	21	CELL/ 99	Serratia liquefaciens	36	XYLO/100	3424167 Serratia liquefaciens		69 AROL/ 0
3424077 Enterobacter aerogenes	99		3424143 Serratia marcescens	99		Enterobacter aerogenes		21 AROL/100
3424100 -7 Serratia marcescens	99		3424144 Serratia liquefaciens	84		Enterobacter cloacae		9
3424116 -7 Serratia marcescens	99		Serratia marcescens	15		3424170 Enterobacter cloacae		60 CELL/ 99
3424120 Serratia liquefaciens	99		3424145 Serratia marcescens	81		Serratia liquefaciens		39 CELL/ 5
3424122 Serratia marcescens	63	XYLO/ 7	Serratia liquefaciens	18		3424171 Enterobacter cloacae		99
Serratia liquefaciens	36	XYLO/100	3424146 Serratia liquefaciens	51	XYLO/100	3424172 Serratia liquefaciens		72 CELL/ 5
3424123 Serratia marcescens	99		Serratia marcescens	48	XYLO/ 7	Enterobacter cloacae		27 CELL/ 99
3424124 Serratia liquefaciens	84		3424147 Serratia marcescens	95		3424173 Enterobacter aerogenes		83
Serratia marcescens	15		Serratia liquefaciens	4		Enterobacter cloacae		16
3424127 Serratia marcescens	99		3424154 Enterobacter cloacae	56	CELL/ 99	3424174 Enterobacter cloacae		57 CELL/ 99
3424128 Serratia marcescens	99		Serratia liquefaciens	42	CELL/ 5	Serratia liquefaciens		42 CELL/ 5
3424129 Serratia marcescens	99		Enterobacter agglomerans	1		3424175 Enterobacter cloacae		76 AROL/ 15
3424130 Serratia marcescens	99		3424155 Enterobacter cloacae	99		Enterobacter aerogenes		18 AROL/100
3424131 Serratia marcescens	99					Serratia liquefaciens		5

BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %	BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %	BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %
3424176				3424334	
Serratia liquefaciens	67 CELL/ 5	3424274		Serratia liquefaciens	87
Enterobacter cloacae	23 CELL/ 99	Serratia liquefaciens	77 CELL/ 5	Enterobacter cloacae	12
Enterobacter aerogenes	9	Enterobacter cloacae	21 CELL/ 99		
3424177		3424275		3424335	
Enterobacter aerogenes	85	Enterobacter cloacae	48 AROL/ 15	Enterobacter cloacae	60 CELL/ 99
Enterobacter cloacae	11	Enterobacter aerogenes	36 AROL/100	Serratia liquefaciens	39 CELL/ 5
Serratia liquefaciens	3	Serratia liquefaciens	15		
3424200		3424276		3424336	
Ewingella americana	99	Serratia liquefaciens	86	Serratia liquefaciens	96
		Enterobacter aerogenes	8	Enterobacter cloacae	3
		Enterobacter cloacae	6		
3424204		3424277		3424337	
Serratia liquefaciens	99	Enterobacter aerogenes	91	Enterobacter aerogenes	81
		Serratia liquefaciens	5	Serratia liquefaciens	18
3424206		Enterobacter cloacae	4	3424340 -3	
Serratia marcescens	89			Serratia liquefaciens	99
Serratia liquefaciens	10				
3424207		3424300		3424344	
Serratia marcescens	99	Serratia marcescens	67 XYLO/ 7	Serratia liquefaciens	98
		Serratia liquefaciens	32 XYLO/100	Enterobacter cloacae	1
3424210		3424301		3424345	
Ewingella americana	99	Serratia marcescens	99	Serratia liquefaciens	87
				Enterobacter cloacae	12
3424220 -7		3424302		3424346 -7	
Serratia liquefaciens	99	Serratia marcescens	91	Serratia liquefaciens	99
		Serratia liquefaciens	8		
3424234 -6		3424303		3424350	
Serratia liquefaciens	99	Serratia marcescens	99	Serratia liquefaciens	68 CELL/ 5
				Enterobacter cloacae	27 CELL/ 99
3424240 -7		3424304		Enterobacter agglomerans	3
Serratia liquefaciens	99	Serratia marcescens	52 XYLO/ 7		
		Serratia liquefaciens	47 XYLO/100	3424351	
3424250				Enterobacter agglomerans	99
Hafnia alvei	99				
3424254		3424305		3424352	
Serratia liquefaciens	55 CELL/ 5	Serratia marcescens	96	Serratia liquefaciens	96
Citrobacter amalonaticus	41 CELL/100	Serratia liquefaciens	3	Enterobacter agglomerans	3
Enterobacter agglomerans	3	3424306			
3424256		Serratia marcescens	85	3424353	
Serratia liquefaciens	99	Serratia liquefaciens	14	Enterobacter aerogenes	99
3424260 -7		3424307		3424354	
Serratia liquefaciens	99	Serratia marcescens	98	Serratia liquefaciens	73 CELL/ 5
				Enterobacter cloacae	26 CELL/ 99
3424270		3424314 -6		3424355	
Serratia liquefaciens	66 CELL/ 5	Serratia liquefaciens	99	Enterobacter cloacae	49 AROL/ 15
Enterobacter cloacae	21 CELL/ 99			Enterobacter aerogenes	37 AROL/100
Hafnia alvei	10	3424320 -7		Serratia liquefaciens	12
		Serratia liquefaciens	99		
3424272		3424330 -2			
Serratia liquefaciens	99	Serratia liquefaciens	99		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3424356								
Serratia liquefaciens	82		3424373			3424614 -4		
Enterobacter aerogenes	9		Enterobacter aerogenes	89		Citrobacter amalonaticus	99	
Enterobacter cloacae	7		Enterobacter cloacae	5				
			Serratia liquefaciens	5		3424650 -4		
3424357						Citrobacter amalonaticus	99	
Enterobacter aerogenes	92		3424374					
Serratia liquefaciens	4		Serratia liquefaciens	72 CELL/	5	3424664		
Enterobacter cloacae	4		Enterobacter cloacae	25 CELL/	99	Citrobacter amalonaticus	84	
						Serratia liquefaciens	15	
3424360			3424375					
Serratia liquefaciens	96		Enterobacter cloacae	49 AROL/	15	3424666		
Serratia rubidaea	2		Enterobacter aerogenes	37 AROL/	100	Serratia liquefaciens	99	
Enterobacter cloacae	1		Serratia liquefaciens	13				
						3424670 -4		
3424361			3424376			Citrobacter amalonaticus	99	
Serratia rubidaea	49 AOL/	85	Serratia liquefaciens	83				
Serratia liquefaciens	44 AROL/	0	Enterobacter aerogenes	9		3424707		
Enterobacter cloacae	7		Enterobacter cloacae	7		Serratia marcescens	99	
3424362			3424377			3424714		
Serratia liquefaciens	99		Enterobacter aerogenes	92		Citrobacter amalonaticus	99	
			Serratia liquefaciens	4				
3424363			Enterobacter cloacae	4		3424724 -4		
Serratia liquefaciens	70 AROL/	0				Serratia liquefaciens	99	
Serratia rubidaea	29 AROL/	85	3424406 -7					
			Serratia marcescens	99		3424744		
3424364						Citrobacter amalonaticus	69 CELL/	100
Serratia liquefaciens	98		3424414 -4			Serratia liquefaciens	30 CELL/	5
Enterobacter cloacae	1		Citrobacter amalonaticus	99				
						3424746		
3424365			3424450 -4			Serratia liquefaciens	99	
Serratia liquefaciens	85		Citrobacter amalonaticus	99				
Enterobacter cloacae	11					3424750		
Enterobacter aerogenes	3		3424455			Citrobacter amalonaticus	99	
			Citrobacter diversus	99				
3424366						3424754		
Serratia liquefaciens	99		3424464			Citrobacter amalonaticus	98	
			Citrobacter amalonaticus	99		Serratia liquefaciens	1	
3424367								
Serratia liquefaciens	77 AROL/	0	3424470 -4			3424760 -2		
Enterobacter aerogenes	19 AROL/	100	Citrobacter amalonaticus	99		Serratia liquefaciens	99	
Enterobacter cloacae	3							
			3424501 -7					
3424370			Serratia marcescens	99		3424764		
Serratia liquefaciens	60 CELL/	5				Serratia liquefaciens	92	
Enterobacter cloacae	24 CELL/	99	3424554			Citrobacter amalonaticus	7	
Enterobacter sakazakii	15		Citrobacter amalonaticus	99				
						3424765 -7		
3424371			3424564 -6			Serratia liquefaciens	99	
Enterobacter cloacae	58 AROL/	15	Serratia liquefaciens	99				
Enterobacter aerogenes	28 AROL/	100				3424770		
Serratia liquefaciens	13		3424574			Kluyvera cryocrescens	51 ESCL/	100
			Citrobacter amalonaticus	87		Citrobacter amalonaticus	36 ESCL/	10
3424372			Serratia liquefaciens	12		Serratia liquefaciens	12	
Serratia liquefaciens	52 CELL/	5						
Enterobacter sakazakii	37 CELL/	100	3424576			3424772		
Enterobacter cloacae	5		Serratia liquefaciens	99		Enterobacter sakazakii	92	
						Serratia liquefaciens	7	
3424356								
						3424772		

BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %	BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %	BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %
3424774 Citrobacter amalonaticus Serratia liquefaciens Kluyvera cryocrescens	69 CELL/100 23 CELL/ 5 7	3426076 Enterobacter cloacae Serratia liquefaciens	80 19	3426147 Serratia marcescens	99
3424775 -6 Serratia liquefaciens	99	3426077 Enterobacter aerogenes Enterobacter cloacae	89 10	3426150 Enterobacter cloacae Enterobacter agglomerans	93 6
3425000 Proteus mirabilis	99	3426100 -7 Serratia marcescens	99	3426151 Enterobacter cloacae Enterobacter agglomerans	93 6
3425106 -7 Serratia marcescens	99	3426120 Serratia liquefaciens	99	3426152 Enterobacter agglomerans	99
3425364 -6 Serratia liquefaciens	99	3426124 Serratia liquefaciens Enterobacter cloacae	77 CELL/ 5 22 CELL/ 99	3426154 Enterobacter cloacae Serratia liquefaciens	94 4
3426001 -7 Serratia marcescens	99	3426125 Serratia marcescens	99	3426155 Enterobacter cloacae Enterobacter aerogenes	78 AROL/ 15 20 AROL/100
3426050 Hafnia alvei	99	3426126 Serratia liquefaciens Serratia marcescens	51 XYLO/100 48 XYLO/ 7	3426156 Enterobacter cloacae Serratia liquefaciens Enterobacter agglomerans	82 15 2
3426054 Enterobacter cloacae Enterobacter agglomerans	90 9	3426127 Serratia marcescens	99	3426157 Enterobacter aerogenes Enterobacter cloacae	89 10
3426055 Citrobacter diversus	99	3426134 Enterobacter cloacae Serratia liquefaciens	89 11	3426160 Serratia liquefaciens Enterobacter cloacae	54 CELL/ 5 45 CELL/ 99
3426060 Serratia liquefaciens	99	3426135 Enterobacter cloacae	99	3426161 Enterobacter cloacae Serratia liquefaciens	89 10
3426064 Serratia liquefaciens Enterobacter cloacae	62 CELL/ 5 37 CELL/ 99	3426136 Enterobacter cloacae Serratia liquefaciens	67 CELL/ 99 32 CELL/ 5	3426162 Serratia liquefaciens Enterobacter cloacae	82 17
3426065 Enterobacter cloacae	99	3426137 Enterobacter aerogenes Enterobacter cloacae	80 AROL/100 19 AROL/ 15	3426164 Serratia liquefaciens Enterobacter cloacae	57 CELL/ 5 42 CELL/ 99
3426066 Serratia liquefaciens	99	3426140 Serratia liquefaciens	99	3426165 Enterobacter cloacae Serratia liquefaciens Enterobacter aerogenes	83 10 6
3426070 Enterobacter cloacae Hafnia alvei Enterobacter agglomerans	62 CELL/ 99 34 CELL/ 15 3	3426144 Serratia liquefaciens Enterobacter cloacae	57 CELL/ 5 42 CELL/ 99	3426166 Serratia liquefaciens Enterobacter cloacae	84 15
3426071 Enterobacter cloacae	99	3426145 Enterobacter cloacae Serratia marcescens	64 CELL/ 99 35 CELL/ 5	3426166 Serratia liquefaciens Enterobacter cloacae	84 15
3426074 Enterobacter cloacae Serratia liquefaciens	93 5	3426146 Serratia liquefaciens Serratia marcescens	51 XYLO/100 48 XYLO/ 7		
3426075 Enterobacter cloacae Enterobacter aerogenes	79 AROL/ 15 20 AROL/100				
3424774				3426166	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3426167								
Enterobacter aerogenes	63	AROL/100	3426255			3426276		
Enterobacter cloacae	24	AROL/ 15	Enterobacter cloacae	90		Enterobacter aerogenes	41	AROL/100
Serratia liquefaciens	12		Enterobacter agglomerans	9		Enterobacter cloacae	30	AROL/ 15
						Serratia liquefaciens	27	
3426170			3426256					
Enterobacter cloacae	79	AROL/ 15	Serratia liquefaciens	76	GELT/ 90	3426277		
Enterobacter gergoviae	17	AROL/ 97	Enterobacter agglomerans	23	GELT/ 2	Enterobacter aerogenes	95	
Serratia liquefaciens	3					Enterobacter cloacae	4	
			3426257					
3426171			Enterobacter aerogenes	99		3426302 -3		
Enterobacter cloacae	85					Serratia marcescens	99	
Enterobacter aerogenes	14		3426260 -2					
			Serratia liquefaciens	99		3426304		
3426172						Serratia marcescens	52	XYLO/ 7
Enterobacter cloacae	53	AROL/ 15	3426264			Serratia liquefaciens	47	XYLO/100
Enterobacter gergoviae	22	AROL/ 97	Serratia liquefaciens	86				
Enterobacter aerogenes	15		Enterobacter cloacae	13		3426305		
						Serratia marcescens	99	
3426173			3426265					
Enterobacter aerogenes	84		Enterobacter cloacae	62	CELL/ 99	3426306		
Enterobacter cloacae	15		Serratia liquefaciens	37	CELL/ 5	Serratia marcescens	85	
						Serratia liquefaciens	14	
3426174			3426266					
Enterobacter cloacae	93		Serratia liquefaciens	95		3426307		
Serratia liquefaciens	4		Enterobacter cloacae	4		Serratia marcescens	99	
Enterobacter aerogenes	1							
			3426267			3426314		
3426175			Enterobacter aerogenes	77	AROL/100	Enterobacter cloacae	64	CELL/ 99
Enterobacter cloacae	79	AROL/ 15	Serratia liquefaciens	22	AROL/ 0	Serratia liquefaciens	29	CELL/ 5
Enterobacter aerogenes	20	AROL/100				Enterobacter agglomerans	6	
			3426270					
3426176			Enterobacter cloacae	79	CELL/ 99	3426316		
Enterobacter cloacae	61	AROL/ 15	Serratia liquefaciens	16	CELL/ 5	Serratia liquefaciens	99	
Enterobacter aerogenes	26	AROL/100	Enterobacter agglomerans	4				
Serratia liquefaciens	11					3426320 -2		
			3426271			Serratia liquefaciens	99	
3426177			Enterobacter cloacae	95				
Enterobacter aerogenes	88		Enterobacter agglomerans	4		3426324		
Enterobacter cloacae	10					Serratia liquefaciens	92	
			3426272			Enterobacter cloacae	7	
3426207			Enterobacter cloacae	52	CELL/ 99			
Serratia marcescens	99		Serratia liquefaciens	41	CELL/ 5	3426325		
			Enterobacter agglomerans	6		Serratia liquefaciens	54	CELL/ 5
3426224 -6						Enterobacter cloacae	45	CELL/ 99
Serratia liquefaciens	99		3426273					
			Enterobacter aerogenes	99		3426326 -7		
3426250						Serratia liquefaciens	99	
Hafnia alvei	86		3426274					
Enterobacter agglomerans	13		Enterobacter cloacae	80		3426330		
			Serratia liquefaciens	18		Enterobacter cloacae	71	CELL/ 99
3426251 -2						Serratia liquefaciens	28	CELL/ 5
Enterobacter agglomerans	99		3426275					
			Enterobacter cloacae	54	AROL/ 15	3426331		
3426254			Enterobacter aerogenes	43	AROL/100	Enterobacter cloacae	99	
Enterobacter cloacae	74	CELL/ 99	Serratia liquefaciens	1				
Serratia liquefaciens	17	CELL/ 5				3426332		
Enterobacter agglomerans	8					Serratia liquefaciens	99	
3426167								
						3426332		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3426334								
Enterobacter cloacae	68	CELL/ 99	3426353	98		3426367		
Serratia liquefaciens	31	CELL/ 5	Enterobacter aerogenes			Enterobacter aerogenes	73	AROL/100
						Serratia liquefaciens	17	AROL/ 0
						Enterobacter cloacae	9	
3426335			3426354					
Enterobacter cloacae	69	AROL/ 15	Enterobacter cloacae	83		3426370		
Enterobacter aerogenes	27	AROL/100	Serratia liquefaciens	15		Enterobacter cloacae	78	CELL/ 99
Serratia liquefaciens	3					Serratia liquefaciens	12	CELL/ 5
			3426355			Enterobacter gergoviae	5	
3426336			Enterobacter cloacae	54	AROL/ 15			
Serratia liquefaciens	52	CELL/ 5	Enterobacter aerogenes	43	AROL/100	3426371		
Enterobacter cloacae	28	CELL/ 99				Enterobacter cloacae	65	ROL/ 15
Enterobacter aerogenes	19		3426356			Enterobacter aerogenes	33	AROL/100
			Enterobacter aerogenes	43	AROL/100			
3426337			Enterobacter cloacae	32	AROL/ 15	3426372		
Enterobacter aerogenes	91		Serratia liquefaciens	23		Enterobacter cloacae	33	AROL/ 15
Enterobacter cloacae	7					Enterobacter aerogenes	29	AROL/100
Serratia liquefaciens	1		3426357			Serratia liquefaciens	21	
			Enterobacter aerogenes	95				
3426340			Enterobacter cloacae	4		3426373		
Serratia liquefaciens	96					Enterobacter aerogenes	93	
Enterobacter agglomerans	3		3426360			Enterobacter cloacae	5	
			Serratia liquefaciens	59	AROL/ 0			
3426342			Serratia rubidaea	27	AROL/ 85	3426374		
Serratia liquefaciens	99		Enterobacter cloacae	13		Enterobacter cloacae	81	
						Serratia liquefaciens	14	
3426344			3426361			Enterobacter aerogenes	3	
Serratia liquefaciens	83		Serratia rubidaea	87				
Enterobacter cloacae	16		Enterobacter cloacae	8		3426375		
			Serratia liquefaciens	3		Enterobacter cloacae	55	AROL/ 15
3426345						Enterobacter aerogenes	44	AROL/100
Enterobacter cloacae	67	CELL/ 99	3426362					
Serratia liquefaciens	31	CELL/ 5	Serratia liquefaciens	81		3426376		
Enterobacter agglomerans	1		Serratia rubidaea	14		Enterobacter aerogenes	44	AROL/100
			Enterobacter cloacae	4		Enterobacter cloacae	32	AROL/ 15
3426346						Serratia liquefaciens	23	
Serratia liquefaciens	94		3426363					
Enterobacter cloacae	4		Serratia rubidaea	66	GELT/ 90	3426377		
			Enterobacter aerogenes	22	GELT/ 0	Enterobacter aerogenes	95	
3426347			Serratia liquefaciens	7		Enterobacter cloacae	4	
Enterobacter aerogenes	73	AROL/100						
Serratia liquefaciens	17	AROL/ 0	3426364			3426415 -5		
Enterobacter cloacae	9		Serratia liquefaciens	83		Citrobacter diversus	99	
			Enterobacter cloacae	16				
3426350						3426502 -7		
Enterobacter cloacae	81		3426365			Serratia marcescens	99	
Serratia liquefaciens	13		Enterobacter cloacae	56	CELL/ 99			
Enterobacter agglomerans	6		Serratia liquefaciens	26	CELL/ 5	3426515 -5		
			Enterobacter aerogenes	14		Citrobacter diversus	99	
3426351								
Enterobacter cloacae	93		3426366			3426564 -6		
Enterobacter agglomerans	6		Serratia liquefaciens	93		Serratia liquefaciens	99	
			Enterobacter cloacae	4				
3426352			Enterobacter aerogenes	2		3426574		
Enterobacter cloacae	55	CELL/ 99				Enterobacter agglomerans	99	
Serratia liquefaciens	35	CELL/ 5						
Enterobacter agglomerans	10							

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3444350 Enterobacter cloacae	90		3444654 -4 Citrobacter amalonaticus	99		3446157 Enterobacter cloacae	99	
Enterobacter agglomerans	9							
3444351 -2 Enterobacter agglomerans	99		3444770 -2 Enterobacter sakazakii	99		3446160 Enterobacter cloacae	99	
3444354 Enterobacter cloacae	97		3444774 Enterobacter agglomerans	99		3446161 Enterobacter cloacae	99	
Enterobacter agglomerans	2							
3444355 Enterobacter cloacae	97		3444777 Klebsiella oxytoca	99		3446162 Enterobacter cloaca	99	
Enterobacter agglomerans	2							
3444356 Enterobacter cloacae	95		3446054 Enterobacter cloacae	91		3446163 Enterobacter cloacae	99	
Enterobacter agglomerans	4		Enterobacter agglomerans	8				
3444357 Enterobacter cloacae	99		3446064 -5 Enterobacter cloacae	99		3446164 Enterobacter cloacae	99	
3444360 Enterobacter cloacae	99		3446070 Enterobacter cloacae	95		3446165 Enterobacter cloacae	99	
			Enterobacter agglomerans	4				
3444361 Enterobacter cloacae	99		3446071 -7 Enterobacter cloacae	99		3446166 Enterobacter cloacae	99	
3444364 Enterobacter cloacae	99		3446107 Serratia marcescens	99		3446170 Enterobacter cloacae	99	
3444365 Enterobacter cloacae	99		3446124 Enterobacter cloacae	99		3446171 Enterobacter cloacae	99	
3444366 Enterobacter cloacae	99		3446130 -7 Enterobacter cloacae	99		3446172 Enterobacter cloacae	98	
3444367 Enterobacter cloacae	99		3446144 -5 Enterobacter cloacae	99		3446173 Enterobacter cloacae	99	
3444370 Enterobacter sakazakii	93		3446150 Enterobacter cloacae	94		3446174 Enterobacter cloacae	99	
Enterobacter cloacae	6		Enterobacter agglomerans	5				
3444371 Enterobacter cloacae	98		3446151 Enterobacter cloacae	94		3446175 Enterobacter cloacae	99	
			Enterobacter agglomerans	5				
3444372 Enterobacter sakazakii	98		3446152 Enterobacter agglomerans	99		3446176 Enterobacter cloacae	99	
3444373 -6 Enterobacter cloacae	99		3446154 -5 Enterobacter cloacae	98		3446177 Enterobacter cloacae	99	
3444377 Enterobacter cloacae	85		3446156 Enterobacter cloacae	97		3446234 Enterobacter cloacae	99	
Klebsiella pneumoniae	14		Enterobacter agglomerans	2				
						3446250 -2 Enterobacter agglomerans	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3446254								
Enterobacter cloacae	91		3446345			3446370		
Enterobacter agglomerans	8		Enterobacter cloacae	97		Enterobacter cloacae	58	
			Enterobacter agglomerans	2		Enterobacter sakazakii	40	
3446255								
Enterobacter cloacae	91		3446346			3446371		
Enterobacter agglomerans	8		Enterobacter cloacae	95		Enterobacter cloacae	99	
			Enterobacter agglomerans	4				
3446256								
Enterobacter agglomerans	99		3446347			3446372		
			Enterobacter cloacae	99		Enterobacter sakazakii	88	
						Enterobacter cloacae	11	
3446264								
Enterobacter cloacae	98		3446350			3446373 -6		
Enterobacter agglomerans	1		Enterobacter cloacae	94		Enterobacter cloacae	98	
			Enterobacter agglomerans	5				
3446265 -6								
Enterobacter cloacae	99		3446351			3446377		
			Enterobacter cloacae	94		Enterobacter cloacae	74 AROL/ 15	
			Enterobacter agglomerans	5		Klebsiella pneumoniae	25 AROL/ 98	
3446270								
Enterobacter cloacae	95		3446352			3446574 -4		
Enterobacter agglomerans	4		Enterobacter cloacae	87		Enterobacter agglomerans	99	
			Enterobacter agglomerans	12				
3446271								
Enterobacter cloacae	95		3446353			3446750 -6		
Enterobacter agglomerans	4		Enterobacter agglomerans	99		Enterobacter agglomerans	99	
3446272								
Enterobacter cloacae	90		3446354 -5			3446770		
Enterobacter agglomerans	9		Enterobacter cloacae	98		Enterobacter sakazakii	98	
3446274 -5								
Enterobacter cloacae	98		3446356			3446771		
			Enterobacter cloacae	97		Enterobacter agglomerans	99	
			Enterobacter agglomerans	2				
3446276								
Enterobacter cloacae	97		3446357			3446772		
Enterobacter agglomerans	2		Enterobacter cloacae	97		Enterobacter sakazakii	99	
			Enterobacter agglomerans	2				
3446277								
Enterobacter cloacae	97		3446360			3446774		
Enterobacter agglomerans	2		Enterobacter cloacae	84		Enterobacter cloacae	83	
			Serratia rubidaea	14		Enterobacter agglomerans	16	
3446314								
Enterobacter cloacae	92		3446361			3446775 -6		
Enterobacter agglomerans	7		Enterobacter cloacae	54 GELT/ 0		Enterobacter agglomerans	99	
			Serratia rubidaea	45 GELT/ 90				
3446324 -5								
Enterobacter cloacae	99		3446362			3447374 -6		
			Enterobacter cloacae	99		Enterobacter cloacae	99	
3446330 -7								
Enterobacter cloacae	99		3446363			3452377		
			Serratia rubidaea	56 GELT/ 90		Klebsiella pneumoniae	82	
			Enterobacter cloacae	43 GELT/ 0		Enterobacter cloacae	17	
3446340								
Enterobacter agglomerans	99		3446364 -7			3452777		
			Enterobacter cloacae	99		Klebsiella oxytoca	99	
3446344								
Enterobacter cloacae	97					3454074		
Enterobacter agglomerans	2					Enterobacter cloacae	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3460364								
Enterobacter cloacae	94		3462360 -7			3464160 -7		
Serratia liquefaciens	5		Enterobacter cloacae	99		Enterobacter cloacae	99	
3460365			3462370			3464170		
Enterobacter cloacae	99		Enterobacter cloacae	50		Enterobacter cloacae	85	
			Enterobacter sakazakii	49		Enterobacter sakazakii	14	
3460366			3462371			3464171		
Enterobacter cloacae	83		Enterobacter cloacae	99		Enterobacter cloacae	99	
Serratia liquefaciens	16							
3460367			3462372			3464172		
Enterobacter cloacae	99		Enterobacter sakazakii	91		Enterobacter sakazakii	64	
			Enterobacter cloacae	8		Enterobacter cloacae	35	
3460370			3462373 -7			3464173 -7		
Enterobacter sakazakii	95		Enterobacter cloacae	99		Enterobacter cloacae	99	
Enterobacter cloacae	4							
3460371			3462455 -5			3464234 -6		
Enterobacter cloacae	99		Citrobacter diversus	99		Enterobacter cloacae	99	
3460372			3462770 -2			3464260 -7		
Enterobacter sakazakii	99		Enterobacter sakazakii	99		Enterobacter cloacae	99	
3460374 -7			3464002 -7			3464270		
Enterobacter cloacae	99		Serratia marcescens	99		Enterobacter sakazakii	61	
						Enterobacter cloacae	38	
3460454 -4			3464050			3464271		
Citrobacter amalonaticus	99		Hafnia alvei	99		Enterobacter cloacae	99	
3460770 -2			3464054 -6			3464272		
Enterobacter sakazakii	99		Enterobacter cloacae	99		Enterobacter sakazakii	94	
						Enterobacter cloacae	5	
3462050			3464070 -7			3464273 -7		
Hafnia alvei	99		Enterobacter cloacae	99		Enterobacter cloacae	99	
3462074 -6			3464100 -7			3464306 -7		
Enterobacter cloacae	99		Serratia marcescens	99		Serratia marcescens	99	
3462160 -7			3464114 -5			3464314 -6		
Enterobacter cloacae	99		Enterobacter cloacae	99		Enterobacter cloacae	99	
3462170 -7			3464127			3464330		
Enterobacter cloacae	99		Serratia marcescens	99		Enterobacter sakazakii	76	
						Enterobacter cloacae	23	
3462264 -5			3464130 -7			3464331		
Enterobacter cloacae	99		Enterobacter cloacae	99		Enterobacter cloacae	99	
3462270 -7			3464144 -6			3464332		
Enterobacter cloacae	99		Enterobacter cloacae	99		Enterobacter sakazakii	97	
						Enterobacter cloacae	2	
3462334 -5			3464147			3464334 -7		
Enterobacter cloacae	99		Serratia marcescens	99		Enterobacter cloacae	99	
3462350 -7			3464150 -7					
Enterobacter cloacae	99		Enterobacter cloacae	99				

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3500077 Klebsiella ozaenae	99		3500273 Klebsiella ozaenae	99		3500370 Enterobacter agglomerans	99	
3500102 -7 Serratia marcescens	99		3500274 Hafnia alvei	78		3500371 Enterobacter agglomerans	99	
3500150 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	21		3500372 Enterobacter agglomerans	99	
3500164 Enterobacter agglomerans	99		3500275 Klebsiella ozaenae	50	MOTL/ 0	3500373 Klebsiella ozaenae	58	MOTL/ 0
3500170 -2 Enterobacter agglomerans	99		Enterobacter agglomerans	49	MOTL/ 85	Enterobacter agglomerans	41	MOTL/ 85
3500173 Klebsiella ozaenae	99		3500276 Enterobacter agglomerans	99		3500374 -6 Enterobacter agglomerans	99	
3500174 -6 Enterobacter agglomerans	99		3500277 Klebsiella ozaenae	86		3500377 Klebsiella pneumoniae	80	
3500177 Klebsiella ozaenae	63	MOTL/ 0	Enterobacter agglomerans	13		Enterobacter aerogenes	16	
Enterobacter agglomerans	36	MOTL/ 85	3500314 Enterobacter agglomerans	99		Enterobacter agglomerans	1	
3500240 -4 Hafnia alvei	99		3500330 -6 Enterobacter agglomerans	99		3500444 Escherichia coli inaktiv	99	
3500250 Hafnia alvei	99		3500340 -6 Enterobacter agglomerans	99		3500450 -4 Escherichia coli inaktiv	99	
3500251 -2 Enterobacter agglomerans	99		3500350 -7 Enterobacter agglomerans	99		3500554 -4 Enterobacter agglomerans	99	
3500254 Hafnia alvei	93		3500360 Serratia rubidaea	97		3500654 Escherichia coli	86	
Enterobacter agglomerans	6		Enterobacter agglomerans	2		Escherichia coli inaktiv	11	
3500255 -7 Enterobacter agglomerans	99		3500361 Serratia rubidaea	99		Enterobacter agglomerans	2	
3500260 Hafnia alvei	99		3500362 Serratia rubidaea	96		3500674 Escherichia coli	93	
3500264 Enterobacter agglomerans	99		Enterobacter agglomerans	3		Enterobacter agglomerans	6	
3500270 Hafnia alvei	98		3500363 Serratia rubidaea	99		3500750 Enterobacter agglomerans	99	
3500271 Klebsiella ozaenae	67	MOTL/ 0	3500364 Serratia liquefaciens	68	GELT/ 90	3500754 Escherichia coli	94	
Enterobacter agglomerans	33	MOTL/ 85	Enterobacter agglomerans	31	GELT/ 2	Enterobacter agglomerans	5	
3500272 Enterobacter agglomerans	99		3500365 Enterobacter agglomerans	99		3500755 -6 Enterobacter agglomerans	99	
			3500366 Serratia liquefaciens	77	GELT/ 90	3500770 -2 Enterobacter agglomerans	99	
			Enterobacter agglomerans	22	GELT/ 2	3500774 Escherichia coli	82	
			3500367 Enterobacter agglomerans	99		Enterobacter agglomerans	17	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3502372 -4 Enterobacter agglomerans	99		3503374 Enterobacter agglomerans	99		3504072 Enterobacter agglomerans	99	
3502375 Klebsiella pneumoniae	83		3504000 Ewingella americana	99		3504073 Klebsiella ozaenae	81	
Enterobacter agglomerans	16		3504001 -7 Serratia marcescens	99		Enterobacter agglomerans	18	
3502376 Klebsiella pneumoniae	79	MOTL/ 0	3504010 Ewingella americana	97		3504074 Hafnia alvei	75	
Enterobacter agglomerans	20	MOTL/ 85	Enterobacter agglomerans	2		Enterobacter agglomerans	24	
3502377 Klebsiella pneumoniae	89		3504014 Enterobacter agglomerans	99		3504075 Enterobacter agglomerans	74	MOTL/ 85
Enterobacter aerogenes	9		3504030 -6 Enterobacter agglomerans	99		Klebsiella ozaenae	25	MOTL/ 0
3502454 -4 Enterobacter agglomerans	99		3504040 Hafnia alvei	98		3504076 Enterobacter agglomerans	99	
3502550 -6 Enterobacter agglomerans	99		3504044 Hafnia alvei	83		3504077 Klebsiella ozaenae	69	MOTL/ 0
3502570 -7 Enterobacter agglomerans	99		Enterobacter agglomerans	16		Enterobacter agglomerans	30	MOTL/ 85
3502650 -6 Enterobacter agglomerans	99		3504045 -6 Enterobacter agglomerans	99		3504100 -7 Serratia marcescens	99	
3502670 -7 Enterobacter agglomerans	99		3504050 Hafnia alvei	99		3504110 -7 Enterobacter agglomerans	99	
3502734 -4 Enterobacter agglomerans	99		3504051 -3 Enterobacter agglomerans	99		3504124 Enterobacter agglomerans	99	
3502750 -7 Enterobacter agglomerans	99		3504054 Hafnia alvei	92		3504126 -7 Serratia marcescens	99	
3502760 Enterobacter agglomerans	99		Enterobacter agglomerans	7		3504130 -7 Enterobacter agglomerans	99	
3502761 Serratia rubidaea	99		3504055 -7 Enterobacter agglomerans	99		3504140 -5 Enterobacter agglomerans	99	
3502764 -6 Enterobacter agglomerans	99		3504060 Hafnia alvei	97		3504146 Serratia marcescens	89	
3502770 -4 Enterobacter agglomerans	99		Enterobacter agglomerans	2		Enterobacter agglomerans	10	
3502775 Klebsiella oxytoca	94		3504061 -7 Enterobacter agglomerans	99		3504147 Serratia marcescens	95	
Enterobacter agglomerans	5		3504070 Hafnia alvei	98		Enterobacter agglomerans	4	
3502776 Klebsiella oxytoca	92		3504071 Enterobacter agglomerans	59	MOTL/ 85	3504150 -7 Enterobacter agglomerans	99	
Enterobacter agglomerans	7		Klebsiella ozaenae	40	MOTL/ 0	3504160 Serratia rubidaea	94	
3502777 Klebsiella oxytoca	99					Enterobacter agglomerans	5	
3502372						3504161 Serratia rubidaea	98	
						3504161		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3504162								
Serratia rubidaea	92		3504250			3504300		
Enterobacter agglomerans	7		Hafnia alvei	97		Enterobacter agglomerans	99	
			Enterobacter agglomerans	2				
3504163						3504302 -3		
Serratia rubidaea	98		3504251 -3			Serratia marcescens	99	
Enterobacter agglomerans	1		Enterobacter agglomerans	99				
						3504304		
3504164			3504254			Serratia marcescens	88	
Enterobacter agglomerans	66	GELT/ 2	Hafnia alvei	63		Enterobacter agglomerans	11	
Serratia liquefaciens	33	GELT/ 90	Enterobacter agglomerans	36				
						3504305		
3504165 -7			3504255 -7			Serratia marcescens	94	
Enterobacter agglomerans	99		Enterobacter agglomerans	99		Enterobacter agglomerans	5	
3504170 -2			3504260			3504306		
Enterobacter agglomerans	99		Hafnia alvei	85		Serratia marcescens	98	
			Enterobacter agglomerans	14		Enterobacter agglomerans	1	
3504173			3504261			3504307		
Enterobacter agglomerans	68	MOTL/ 85	Serratia rubidaea	91		Serratia marcescens	99	
Klebsiella ozaenae	31	MOTL/ 0	Enterobacter agglomerans	8				
						3504310 -7		
3504174 -6			3504262 -7			Enterobacter agglomerans	99	
Enterobacter agglomerans	99		Enterobacter agglomerans	99				
3504177			3504270			3504320		
Klebsiella pneumoniae	66	MOTL/ 0	Hafnia alvei	92		Serratia rubidaea	97	
Enterobacter aerogenes	29	MOTL/ 97	Enterobacter agglomerans	7		Enterobacter agglomerans	2	
Enterobacter agglomerans	3							
						3504321		
3504200			3504271			Serratia rubidaea	99	
Ewingella americana	99		Enterobacter agglomerans	77	MOTL/ 85			
			Klebsiella ozaenae	22	MOTL/ 0	3504322		
3504204						Serratia rubidaea	96	
Enterobacter agglomerans	99		3504272			Enterobacter agglomerans	3	
			Enterobacter agglomerans	99				
						3504323		
3504210			3504273			Serratia rubidaea	99	
Ewingella americana	97		Klebsiella ozaenae	64	MOTL/ 0			
Enterobacter agglomerans	2		Enterobacter agglomerans	35	MOTL/ 85	3504324 -7		
						Enterobacter agglomerans	99	
3504211 -7			3504274					
Enterobacter agglomerans	99		Enterobacter agglomerans	68		3504330 -6		
			Hafnia alvei	31		Enterobacter agglomerans	99	
3504224			3504275			3504337		
Enterobacter agglomerans	99		Enterobacter agglomerans	87		Klebsiella pneumoniae	95	
			Klebsiella ozaenae	12		Enterobacter agglomerans	4	
3504230 -7			3504276			3504340		
Enterobacter agglomerans	99		Enterobacter agglomerans	99		Serratia rubidaea	91	
						Enterobacter agglomerans	8	
			3504277					
3504240			Klebsiella pneumoniae	90		3504341		
Hafnia alvei	95		Enterobacter agglomerans	5		Serratia rubidaea	97	
Enterobacter agglomerans	4		Klebsiella ozaenae	4		Enterobacter agglomerans	2	
3504241 -7								
Enterobacter agglomerans	99							

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3504773								
Klebsiella oxytoca	97		3506044			3506120		
Enterobacter agglomerans	2		Hafnia alvei	61		Enterobacter agglomerans	99	
			Enterobacter agglomerans	38				
3504774						3506121		
Enterobacter agglomerans	99		3506045 -7			Serratia rubidaea	94	
			Enterobacter agglomerans	99		Enterobacter agglomerans	5	
3504775						3506122 -6		
Klebsiella oxytoca	96		3506050			Enterobacter agglomerans	99	
Enterobacter agglomerans	3		Hafnia alvei	98				
3504776						3506127		
Klebsiella oxytoca	95		3506051 -3			Serratia marcescens	91	
Enterobacter agglomerans	4		Enterobacter agglomerans	99		Enterobacter agglomerans	8	
3504777			3506054			3506130 -7		
Klebsiella oxytoca	99		Hafnia alvei	79		Enterobacter agglomerans	99	
			Enterobacter agglomerans	20				
3505154 -4						3506140		
Enterobacter agglomerans	99		3506055 -7			Enterobacter agglomerans	99	
			Enterobacter agglomerans	99				
3505350 -6						3506141		
Enterobacter agglomerans	99		3506060			Serratia rubidaea	79 GELT/ 90	
			Hafnia alvei	92		Enterobacter agglomerans	20 GELT/ 2	
			Enterobacter agglomerans	7				
3505361						3506142 -6		
Serratia rubidaea	99		3506061 -7			Enterobacter agglomerans	99	
			Enterobacter agglomerans	99				
3505370 -7						3506147		
Enterobacter agglomerans	99		3506070			Serratia marcescens	72 GELT/ 90	
			Hafnia alvei	96		Enterobacter agglomerans	27 GELT/ 2	
3505774			Enterobacter agglomerans	3				
Enterobacter agglomerans	99					3506150 -7		
						Enterobacter agglomerans	99	
3506003			3506071 -3					
Serratia marcescens	99		Enterobacter agglomerans	99		3506160		
						Serratia rubidaea	97	
3506004						Enterobacter agglomerans	2	
Serratia marcescens	92		3506074					
Enterobacter agglomerans	7		Enterobacter agglomerans	50		3506161		
			Hafnia alvei	49		Serratia rubidaea	99	
3506005 -7								
Serratia marcescens	99		3506075 -7			3506162		
			Enterobacter agglomerans	99		Serratia rubidaea	96	
						Enterobacter agglomerans	3	
3506010 -7								
Enterobacter agglomerans	99		3506100			3506163		
			Serratia marcescens	94		Serratia rubidaea	98	
			Enterobacter agglomerans	5		Enterobacter agglomerans		
3506024								
Enterobacter agglomerans	99		3506101 -3			3506164		
			Serratia marcescens	99		Serratia rubidaea	63 GELT/ 90	
						Enterobacter agglomerans	36 GELT/ 2	
3506030 -7								
Enterobacter agglomerans	99		3506104			3506165		
			Serratia marcescens	98		Serratia rubidaea	89	
			Enterobacter agglomerans	1		Enterobacter agglomerans	10	
3506040								
Hafnia alvei	97		3506105 -7			3506165		
Enterobacter agglomerans	2		Serratia marcescens	98		Serratia rubidaea	89	
						Enterobacter agglomerans	10	
3506041 -3								
Enterobacter agglomerans	99		3506110 -7					
			Enterobacter agglomerans	99				
3504773								
						3506165		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3516110 -7 Enterobacter agglomerans	99		3516173 Klebsiella pneumoniae	95		3516270 Enterobacter agglomerans	67	
3516120 -4 Enterobacter agglomerans	99		Enterobacter agglomerans	4		Hafnia alvei	32	
3516130 -6 Enterobacter agglomerans	99		3516174 Enterobacter agglomerans	99		3516271 -2 Enterobacter agglomerans	99	
3516137 Klebsiella pneumoniae	97		3516175 Klebsiella pneumoniae	93		3516273 Klebsiella pneumoniae	95	
Enterobacter agglomerans	2		Enterobacter agglomerans	6		Enterobacter agglomerans	4	
3516140 -7 Enterobacter agglomerans	99		3516176 Klebsiella pneumoniae	91		3516274 Enterobacter agglomerans	99	
3516150 -6 Enterobacter agglomerans	99		Enterobacter agglomerans	8		3516275 Klebsiella pneumoniae	93	
3516157 Klebsiella pneumoniae	90		3516177 Klebsiella pneumoniae	98		Enterobacter agglomerans	6	
Enterobacter agglomerans	9		3516204 Enterobacter agglomerans	99		3516276 Klebsiella pneumoniae	91	
3516160 Serratia rubidaea	85		Enterobacter agglomerans	99		Enterobacter agglomerans	8	
Enterobacter agglomerans	14		3516210 -7 Enterobacter agglomerans	99		3516277 Klebsiella pneumoniae	99	
3516161 Serratia rubidaea	96		3516220 -4 Enterobacter agglomerans	99		3516300 -5 Enterobacter agglomerans	99	
Enterobacter agglomerans	3		3516230 -6 Enterobacter agglomerans	99		3516306 Serratia marcescens	87	
3516162 Enterobacter agglomerans	99		3516237 Klebsiella pneumoniae	97		Enterobacter agglomerans	12	
3516163 Serratia rubidaea	94		Enterobacter agglomerans	2		3516307 Serratia marcescens	99	
Enterobacter agglomerans	5		3516240 -7 Enterobacter agglomerans	99		3516310 -6 Enterobacter agglomerans	99	
3516164 -6 Enterobacter agglomerans	99		3516250 Hafnia alvei	65		3516317 Klebsiella pneumoniae	91	
3516167 Klebsiella pneumoniae	95		Enterobacter agglomerans	34		Enterobacter agglomerans	8	
Enterobacter agglomerans	4		3516251 -6 Enterobacter agglomerans	99		3516320 Serratia rubidaea	93	
3516170 Enterobacter gergoviae	88		3516257 Klebsiella pneumoniae	90		Enterobacter agglomerans	6	
Enterobacter agglomerans	11		Enterobacter agglomerans	9		3516321 Serratia rubidaea	98	
3516171 Enterobacter agglomerans	99		3516260 -6 Enterobacter agglomerans	99		Enterobacter agglomerans	1	
3516172 Enterobacter gergoviae	85		3516267 Klebsiella pneumoniae	95		3516322 Enterobacter agglomerans	99	
Enterobacter agglomerans	14		Enterobacter agglomerans	4		3516323 Serratia rubidaea	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3516324 -6 Enterobacter agglomerans	99		3516354 Enterobacter agglomerans	99		3516374 Klebsiella pneumoniae		84
3516327 Klebsiella pneumoniae	96		3516355 Klebsiella pneumoniae	81		Enterobacter agglomerans		15
Enterobacter agglomerans	3		Enterobacter agglomerans	18		3516375 -7 Klebsiella pneumoniae		98
3516330 -2 Enterobacter agglomerans	99		3516356 Klebsiella pneumoniae	76	HOTL/ 0	3516450 -7 Enterobacter agglomerans		99
3516333 Klebsiella pneumoniae	96		Enterobacter agglomerans	23	HOTL/ 85	3516464 Enterobacter agglomerans		99
Enterobacter agglomerans	3		3516357 Klebsiella pneumoniae	98		3516470 -6 Enterobacter agglomerans		99
3516334 Enterobacter agglomerans	99		3516360 -2 Serratia rubidaea	99		3516477 Klebsiella oxytoca		98
3516335 Klebsiella pneumoniae	94		3516363 Serratia rubidaea	97		3516514 Enterobacter agglomerans		99
Enterobacter agglomerans	5		Klebsiella pneumoniae	2		3516530 -6 Enterobacter agglomerans		99
3516336 Klebsiella pneumoniae	93		3516364 Serratia rubidaea	88		3516537 Klebsiella oxytoca		99
Enterobacter agglomerans	6		Enterobacter agglomerans	11		3516540 -6 Enterobacter agglomerans		99
3516337 Klebsiella pneumoniae	99		3516365 Serratia rubidaea	75	GELT/ 90	3516550 -6 Enterobacter agglomerans		99
3516340 Serratia rubidaea	77	GELT/ 90	Klebsiella pneumoniae	22	GELT/ 0	3516557 Klebsiella oxytoca		96
Enterobacter agglomerans	22	GELT/ 2	Enterobacter agglomerans	2		Enterobacter agglomerans		3
3516341 Serratia rubidaea	94		3516366 Klebsiella pneumoniae	61	GELT/ 0	3516560 -6 Enterobacter agglomerans		99
Enterobacter agglomerans	5		Serratia rubidaea	32	GELT/ 90	3516567 Klebsiella oxytoca		99
3516342 Enterobacter agglomerans	99		Enterobacter agglomerans	7		3516570 -2 Enterobacter agglomerans		99
3516343 Serratia rubidaea	91		3516367 Klebsiella pneumoniae	94		3516575 Klebsiella oxytoca		98
Enterobacter agglomerans	8		Klebsiella pneumoniae	5		Enterobacter agglomerans		1
3516344 -6 Enterobacter agglomerans	99		Serratia rubidaea			3516574 Enterobacter agglomerans		99
3516347 Klebsiella pneumoniae	87		3516370 Enterobacter gergoviae	59		3516575 Klebsiella oxytoca		97
Enterobacter agglomerans	12		Enterobacter agglomerans	24		Enterobacter agglomerans		2
3516350 -2 Enterobacter agglomerans	99		Serratia rubidaea	16		3516575 Klebsiella oxytoca		97
3516353 Klebsiella pneumoniae	86		3516371 Klebsiella pneumoniae	70	GELT/ 0	Enterobacter agglomerans		99
Enterobacter agglomerans	13		Klebsiella pneumoniae	23	GELT/ 90	3516575 Klebsiella oxytoca		98
3516324			Enterobacter agglomerans	6		Enterobacter agglomerans		1
			3516372 Klebsiella pneumoniae	73	KCNT/ 98	3516574 Enterobacter agglomerans		99
			Enterobacter gergoviae	17	KCNT/ 0	3516575 Klebsiella oxytoca		97
			Enterobacter agglomerans	9		Enterobacter agglomerans		2
			3516373 Klebsiella pneumoniae	98		3516575 Klebsiella oxytoca		97
						Enterobacter agglomerans		2
						3516575		3516575

BIOTYPE NUMBER	RP	ADD.TEST	BIOTYPE NUMBER	RP	ADD.TEST	BIOTYPE NUMBER	RP	ADD.TEST
ORGANISM NAME	%	POS. %	ORGANISM NAME	%	POS. %	ORGANISM NAME	%	POS. %
3516576			3516717			3516760		
Klebsiella oxytoca	97		Klebsiella oxytoca	97		Serratia rubidaea	88	
Enterobacter agglomerans	2		Enterobacter agglomerans	2		Enterobacter agglomerans	11	
3516577			3516724			3516761		
Klebsiella oxytoca	99		Enterobacter agglomerans	99		Serratia rubidaea	97	
3516614			3516727			Enterobacter agglomerans	2	
Enterobacter agglomerans	99		Klebsiella oxytoca	99		3516762		
3516630 -6			3516730 -2			Enterobacter agglomerans	99	
Enterobacter agglomerans	99		Enterobacter agglomerans	99		3516763		
3516637			3516733			Klebsiella oxytoca	67	GELT/ 0
Klebsiella oxytoca	99		Klebsiella oxytoca	99		Serratia rubidaea	30	GELT/ 90
3516640 -6			3516734			Enterobacter agglomerans	1	
Enterobacter agglomerans	99		Enterobacter agglomerans	99		3516764		
3516650 -6			3516735			Enterobacter agglomerans	99	
Enterobacter agglomerans	99		Klebsiella oxytoca	98		3516765		
3516657			Enterobacter agglomerans	1		Klebsiella oxytoca	97	
Klebsiella oxytoca	96		3516736			Enterobacter agglomerans	2	
Enterobacter agglomerans	3		Klebsiella oxytoca	97		3516766		
3516660 -6			Enterobacter agglomerans	2		Klebsiella oxytoca	96	
Enterobacter agglomerans	99		3516737			Enterobacter agglomerans	3	
3516667			Klebsiella oxytoca	99		3516767		
Klebsiella oxytoca	99		3516740 -6			Klebsiella oxytoca	99	
3516670 -2			Enterobacter agglomerans	99		3516770		
Enterobacter agglomerans	99		3516747			Enterobacter agglomerans	99	
3516673			Klebsiella oxytoca	95		3516771		
Klebsiella oxytoca	98		Enterobacter agglomerans	4		Klebsiella oxytoca	96	
Enterobacter agglomerans	1		3516750 -2			Enterobacter agglomerans	3	
3516674			Enterobacter agglomerans	99		3516772		
Enterobacter agglomerans	99		3516753			Klebsiella oxytoca	95	
3516675			Klebsiella oxytoca	95		Enterobacter agglomerans	4	
Klebsiella oxytoca	97		Enterobacter agglomerans	4		3516773		
Enterobacter agglomerans	2		3516754			Klebsiella oxytoca	99	
3516676			Enterobacter agglomerans	99		3516774		
Klebsiella oxytoca	97		3516755			Klebsiella oxytoca	94	
Enterobacter agglomerans	2		Klebsiella oxytoca	93		Enterobacter agglomerans	5	
3516677			Enterobacter agglomerans	6		3516775 -7		
Klebsiella oxytoca	99		3516756			Klebsiella oxytoca	99	
3516710 -6			Klebsiella oxytoca	91		3517154		
Enterobacter agglomerans	99		Enterobacter agglomerans	8		Enterobacter agglomerans	99	
3516757			3517170 -4			Enterobacter agglomerans	99	
Klebsiella oxytoca	99		3517170 -4			Enterobacter agglomerans	99	
3516576			3517170			3517170		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3524250 -2 Hafnia alvei	99		3524277 Enterobacter aerogenes	98		3524335 Enterobacter aerogenes Serratia liquefaciens	76	AROL/100 23 AROL/ 0
3524254 Hafnia alvei Serratia liquefaciens	96 3		3524300 Serratia marcescens Serratia liquefaciens	89 10		3524336 Serratia liquefaciens Enterobacter aerogenes	88 11	
3524256 Serratia liquefaciens	99		3524301 Serratia marcescens	99		3524337 Enterobacter aerogenes Serratia liquefaciens	95 4	
3524257 Enterobacter aerogenes	99		3524302 Serratia marcescens Serratia liquefaciens	97 2		3524340 -5 Serratia liquefaciens	99	
3524260 Hafnia alvei Serratia liquefaciens	76 ESCL/ 7 23 ESCL/ 97		3524303 Serratia marcescens	99		3524346 Serratia liquefaciens Serratia marcescens	98 1	
3524261 -3 Serratia liquefaciens	99		3524304 Serratia marcescens Serratia liquefaciens	81 18		3524347 Enterobacter aerogenes Serratia liquefaciens Serratia marcescens	51 AROL/100 35 AROL/ 0 12	
3524264 Serratia liquefaciens Hafnia alvei	98 1		3524305 Serratia marcescens	98		3524350 Serratia liquefaciens Hafnia alvei	74 ESCL/ 97 25 ESCL/ 7	
3524265 -6 Serratia liquefaciens	99		3524306 Serratia marcescens Serratia liquefaciens	95 4		3524351 Enterobacter aerogenes	99	
3524267 Enterobacter aerogenes Serratia liquefaciens	53 AROL/100 46 AROL/ 0		3524307 Serratia marcescens	99		3524352 Serratia liquefaciens	99	
3524270 Hafnia alvei	98		3524314 -6 Serratia liquefaciens	99		3524353 Enterobacter aerogenes	99	
3524271 Hafnia alvei	99		3524317 Enterobacter aerogenes	99		3524354 Serratia liquefaciens Enterobacter aerogenes	92 7	
3524272 Serratia liquefaciens Hafnia alvei	87 12		3524320 -5 Serratia liquefaciens	99		3524355 Enterobacter aerogenes Serratia liquefaciens	93 6	
3524273 Enterobacter aerogenes	99		3524326 Serratia liquefaciens Serratia marcescens	98 1		3524356 Serratia liquefaciens Enterobacter aerogenes	60 AROL/ 0 39 AROL/100	
3524274 Serratia liquefaciens Hafnia alvei Enterobacter aerogenes	61 ESCL/ 97 34 ESCL/ 7 4		3524327 Serratia liquefaciens Serratia marcescens	73 XYLO/100 26 XYLO/ 7		3524357 Enterobacter aerogenes	98	
3524275 Enterobacter aerogenes Serratia liquefaciens	92 7		3524330 -2 Serratia liquefaciens	99		3524360 Serratia liquefaciens Serratia rubidáea	96 3	
3524276 Serratia liquefaciens Enterobacter aerogenes	66 AROL/ 0 33 AROL/100		3524333 Enterobacter aerogenes	99		3524334 Serratia liquefaciens	99	
3524250						3524360		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3526171 Enterobacter aerogenes	99		3526255 Enterobacter aerogenes	99		3526276 Enterobacter aerogenes	89	
3526172 Enterobacter aerogenes	68	KCNT/ 98	3526256 Enterobacter aerogenes	89		Serratia liquefaciens	10	
Enterobacter gergoviae	24	KCNT/ 0	Serratia liquefaciens	10		3526277 Enterobacter aerogenes	99	
Serratia liquefaciens	7		3526257 Enterobacter aerogenes	99		3526300 -3 Serratia marcescens	99	
3526173 Enterobacter aerogenes	99		3526260 Hafnia alvei	89		Serratia marcescens	81	
3526174 Enterobacter aerogenes	56	AROL/100	Serratia liquefaciens	10		Serratia liquefaciens	18	
Serratia liquefaciens	34	AROL/ 0	3526262 Serratia liquefaciens	99		3526305 Serratia marcescens	99	
Kluyvera ascorbata	9		3526263 Enterobacter aerogenes	99		3526306 Serratia marcescens	95	
3526175 Enterobacter aerogenes	99		3526264 Serratia liquefaciens	95		Serratia liquefaciens	4	
3526176 Enterobacter aerogenes	92		Hafnia alvei	4		3526307 Serratia marcescens	99	
Serratia liquefaciens	7		3526265 Enterobacter aerogenes	70	AROL/100	3526314 -6 Serratia liquefaciens	99	
3526177 Enterobacter aerogenes	99		Serratia liquefaciens	29	AROL/ 0	3526317 Enterobacter aerogenes	99	
3526203 -7 Serratia marcescens	99		3526266 Serratia liquefaciens	91		3526320 Serratia liquefaciens	99	
3526210 Hafnia alvei	99		Enterobacter aerogenes	8		3526321 Serratia rubidaea	99	
3526224 -6 Serratia liquefaciens	99		3526267 Enterobacter aerogenes	94		3526322 -6 Serratia liquefaciens	99	
3526237 Enterobacter aerogenes	99		Serratia liquefaciens	5		3526327 Enterobacter aerogenes	82	
3526240 Hafnia alvei	99		3526270 Hafnia alvei	99		Serratia liquefaciens	17	
3526244 Hafnia alvei	69	ESCL/ 7	3526271 Enterobacter aerogenes	99		3526330 Kluyvera ascorbata	74	CELL/100
Serratia liquefaciens	30	ESCL/ 97	3526272 Enterobacter aerogenes	86		Serratia liquefaciens	25	CELL/ 5
3526246 Serratia liquefaciens	99		Serratia liquefaciens	13		3526332 Serratia liquefaciens	99	
3526250 -2 Hafnia alvei	99		3526273 Enterobacter aerogenes	99		3526333 Enterobacter aerogenes	99	
3526253 Enterobacter aerogenes	99		3526274 Hafnia alvei	42	AROL/ 0	3526334 Serratia liquefaciens	65	CELL/ 5
3526254 Hafnia alvei	98		Enterobacter aerogenes	29	AROL/100	Kluyvera ascorbata	17	CELL/100
Serratia liquefaciens	1		Serratia liquefaciens	28		Enterobacter aerogenes	17	
3526171			3526275 Enterobacter aerogenes	99		3526334		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3526335								
Enterobacter aerogenes	97		3526355			3526372		
Serratia liquefaciens	2		Enterobacter aerogenes	99		Enterobacter aerogenes	85	
						Serratia liquefaciens	11	
3526336			3526356			Enterobacter gergoviae	3	
Enterobacter aerogenes	67	AROL/100	Enterobacter aerogenes	91				
Serratia liquefaciens	32	AROL/ 0	Serratia liquefaciens	8		3526373		
						Enterobacter aerogenes	99	
3526337			3526357					
Enterobacter aerogenes	99		Enterobacter aerogenes	99		3526374		
						Enterobacter aerogenes	44	AROL/100
3526340			3526360			Serratia liquefaciens	33	AROL/ 0
Serratia liquefaciens	99		Serratia liquefaciens	57	AROL/ 0	Kluyvera ascorbata	22	
			Serratia rubidaea	42	AROL/ 85			
3526341			3526361			3526375		
Serratia rubidaea	99		Serratia rubidaea	92		Enterobacter aerogenes	99	
			Enterobacter aerogenes	5				
3526342			Serratia liquefaciens	2		3526376		
Serratia liquefaciens	99					Enterobacter aerogenes	91	
						Serratia liquefaciens	8	
3526343			3526362					
Enterobacter aerogenes	99		Serratia liquefaciens	73	AROL/ 0	3526377		
			Serratia rubidaea	20	AROL/ 85	Enterobacter aerogenes	99	
3526344			Enterobacter aerogenes	6				
Serratia liquefaciens	99					3526406 -7		
						Serratia marcescens	99	
3526345			3526363					
Enterobacter aerogenes	74	AROL/100	Enterobacter aerogenes	51	GELT/ 0	3526450		
Serratia liquefaciens	25	AROL/ 0	Serratia rubidaea	45	GELT/ 90	Hafnia alvei	99	
			Serratia liquefaciens	3				
3526346			3526364					
Serratia liquefaciens	89		Serratia liquefaciens	97		3526470		
Enterobacter aerogenes	10		Enterobacter aerogenes	1		Hafnia alvei	99	
3526347			3526365			3526500 -7		
Enterobacter aerogenes	95		Enterobacter aerogenes	71	AROL/100	Serratia marcescens	99	
Serratia liquefaciens	4		Serratia liquefaciens	24	AROL/ 0			
			Serratia rubidaea	5		3526550 -4		
3526350						Kluyvera ascorbata	99	
Kluyvera ascorbata	79	CELL/100	3526366					
Serratia liquefaciens	10	CELL/ 5	Serratia liquefaciens	89		3526564 -6		
Hafnia alvei	9		Enterobacter aerogenes	10		Serratia liquefaciens	99	
3526351			3526367			3526570 -4		
Enterobacter aerogenes	99		Enterobacter aerogenes	95		Kluyvera ascorbata	99	
			Serratia liquefaciens	4				
3526352						3526577		
Enterobacter aerogenes	88		3526370			Enterobacter aerogenes	99	
Serratia liquefaciens	11		Kluyvera ascorbata	75	CELL/100			
			Serratia liquefaciens	10	CELL/ 5	3526650		
3526353			Enterobacter aerogenes	9		Hafnia alvei	99	
Enterobacter aerogenes	99							
						3526664 -6		
3526354			3526371			Serratia liquefaciens	99	
Enterobacter aerogenes	44	AROL/100	Enterobacter aerogenes	98				
Serratia liquefaciens	33	AROL/ 0				3526670 -4		
Kluyvera ascorbata	22					Kluyvera ascorbata	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3532164 -6 Serratia liquefaciens	99		3532374 Serratia liquefaciens	99		3534120 -3 Serratia marcescens	99	
3532170 -2 Enterobacter gergoviae	99		3532375 Enterobacter aerogenes	99		3534124 Serratia marcescens	75	XYLO/ 7
3532174 Serratia liquefaciens	99		3532376 Enterobacter aerogenes	82		Serratia liquefaciens	24	XYLO/100
3532175 Enterobacter aerogenes	99		Serratia liquefaciens	17		3534125 Serratia marcescens	99	
3532176 Serratia liquefaciens	99		3532377 Enterobacter aerogenes	99		3534126 Serratia marcescens	93	
3532177 Enterobacter aerogenes	99		3532777 Klebsiella ornithinolyt.	99		Serratia liquefaciens	6	
3532240 Hafnia alvei	99		3534000 -7 Serratia marcescens	99		3534127 Serratia marcescens	99	
3532250 -4 Hafnia alvei	99		3534026 -7 Serratia marcescens	99		3534134 -6 Serratia liquefaciens	99	
3532260 Hafnia alvei	99		3534040 -4 Hafnia alvei	99		3534140 -3 Serratia marcescens	99	
3532264 -6 Serratia liquefaciens	99		3534046 -7 Serratia marcescens	99		3534144 Serratia marcescens	75	XYLO/ 7
3532270 -4 Hafnia alvei	99		3534050 -4 Hafnia alvei	99		Serratia liquefaciens	24	XYLO/100
3532277 Enterobacter aerogenes	99		3534064 -6 Serratia liquefaciens	99		3534145 Serratia marcescens	99	
3532307 Serratia marcescens	99		3534060 Hafnia alvei	99		3534146 Serratia marcescens	93	
3532324 -6 Serratia liquefaciens	99		3534064 -6 Serratia liquefaciens	99		Serratia liquefaciens	6	
3532357 Enterobacter aerogenes	99		3534070 Hafnia alvei	99		3534147 Serratia marcescens	99	
3532360 -6 Serratia liquefaciens	99		3534074 Hafnia alvei	87		3534154 -6 Serratia liquefaciens	99	
3532367 Enterobacter aerogenes	91		3534076 Serratia liquefaciens	12		3534157 Enterobacter aerogenes	99	
Serratia liquefaciens	8		3534077 Enterobacter aerogenes	99		3534160 -6 Serratia liquefaciens	99	
3532370 -2 Enterobacter gergoviae	99		3534100 -7 Serratia marcescens	99		3534167 Enterobacter aerogenes	66	AROL/100
3532373 Enterobacter aerogenes	99		3534116 -7 Serratia marcescens	99		Serratia liquefaciens	33	AROL/ 0
3532164						3534170 Enterobacter gergoviae	94	
						Serratia liquefaciens	5	
						3534172 Enterobacter gergoviae	89	
						Serratia liquefaciens	10	
						3534172		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3534173 Enterobacter aerogenes	99		3534300 -3 Serratia marcescens	99		3534370 Enterobacter gergoviae Serratia liquefaciens	60 39	AROL/ AROL/ 97 0
3534174 Serratia liquefaciens	99		3534304 Serratia marcescens Serratia liquefaciens	94 5		3534371 Enterobacter aerogenes Serratia liquefaciens	92 7	
3534175 Enterobacter aerogenes Serratia liquefaciens	95 4		3534305 -7 Serratia marcescens	99		3534372 Serratia liquefaciens Enterobacter gergoviae Enterobacter aerogenes	43 33 23	AROL/ AROL/ 97 0
3534176 Serratia liquefaciens Enterobacter aerogenes	53 46	AROL/ AROL/100 0	3534320 -6 Serratia liquefaciens	99		3534373 Enterobacter aerogenes	98	
3534177 Enterobacter aerogenes	98		3534327 Serratia marcescens Serratia liquefaciens	60 39	XYLO/ XYLO/100 7	3534374 Serratia liquefaciens Enterobacter aerogenes	91 8	
3534202 -7 Serratia marcescens	99		3534330 -6 Serratia liquefaciens	99		3534375 Enterobacter aerogenes Serratia liquefaciens	94 5	
3534224 -6 Serratia liquefaciens	99		3534337 Enterobacter aerogenes	99		3534376 Serratia liquefaciens Enterobacter aerogenes	58 41	AROL/ AROL/100 0
3534240 Hafnia alvei	99		3534340 -6 Serratia liquefaciens	99		3534377 Enterobacter aerogenes	98	
3534244 -6 Serratia liquefaciens	99		3534347 Serratia marcescens Serratia liquefaciens	60 39	XYLO/ XYLO/100 7	3534403 -7 Serratia marcescens	99	
3534250 -4 Hafnia alvei	99		3534350 -4 Serratia liquefaciens	99		3534500 -7 Serratia marcescens	99	
3534260 Hafnia alvei Serratia liquefaciens	63 36	ESCL/ ESCL/ 7 97	3534355 Enterobacter aerogenes	99		3534564 Serratia liquefaciens	99	
3534262 -7 Serratia liquefaciens	99		3534356 Serratia liquefaciens Enterobacter aerogenes	58 41	AROL/ AROL/100 0	3534706 -7 Serratia marcescens	99	
3534270 Hafnia alvei Serratia liquefaciens	97 2		3534357 Enterobacter aerogenes	99		3534764 -6 Serratia liquefaciens	99	
3534272 Serratia liquefaciens	99		3534360 -4 Serratia liquefaciens	99		3534777 Klebsiella ornithinolyt. Klebsiella oxytoca	92 7	
3534274 Serratia liquefaciens Hafnia alvei	77 22	ESCL/ ESCL/ 97 7	3534365 Serratia liquefaciens Enterobacter aerogenes	83 16		3535000 Proteus mirabilis	99	
3534275 Enterobacter aerogenes	99		3534366 Serratia liquefaciens	99		3535106 -7 Serratia marcescens	99	
3534276 Serratia liquefaciens	99		3534367 Enterobacter aerogenes Serratia liquefaciens	61 38	AROL/ AROL/ 100 0			
3534277 Enterobacter aerogenes	99							
3534173							3535106	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
3726050 Hafnia alvei	99		3764050 Hafnia alvei	99				
3726070 Hafnia alvei	99		3764054 Salmonella subgroup 5	99	SEROLOGY			
3726106 -7 Serratia marcescens	99		3766050 Hafnia alvei	99				
3726177 Enterobacter aerogenes	99		3766054 Salmonella subgroup 2	99	SEROLOGY			
3726250 Hafnia alvei	99							
3726364 -6 Serratia liquefaciens	99							
3726375 -7 Enterobacter aerogenes	99							
3731000 Proteus mirabilis	99							
3734106 -7 Serratia marcescens	99							
3735000 Proteus mirabilis	99							
3736107 Serratia marcescens	99							
3736170 -2 Enterobacter gergoviae	99							
3736370 -2 Enterobacter gergoviae	99							
3736377 Enterobacter aerogenes	99							
3736777 Klebsiella ornithinolyt.	99							
3760050 Hafnia alvei	99							
3760054 Salmonella subgroup 3a	85	SEROLOGY						
Salmonella subgroup 3b	14	SEROLOGY						
3760254 Salmonella subgroup 3b	83	SEROLOGY						
Salmonella subgroup 3a	16	SEROLOGY						
3762050 Hafnia alvei	99							
3726050								