



IH-500 The Automated Antibody Titration

BIO-RAD

IH-500 Antibody Titration

Save Time
and Eliminate
Manual Errors





- Bio-Rad is bringing **automated antibody titration** to an outstanding level by offering you the latest robotic procedure in this field.
- **Save time**, performing with a simple click of a button, automated titrations using the IH-500 system empowered with a **ID-Titration solution** specifically engineered for ID-Cards*.
- **Eliminate potential errors** from manual techniques while IH-500 performs automatically **two-fold serial dilutions**.
- According to your needs, IH-500 performs dilution from 1:1 to 1:32, from 1:1 to 1:2048 and from 1:64 to 1:2048.
- The flexibility of the titration procedure allows you to use **LISS Coombs**, **Coombs Anti-IgG** or **NaCl ID-Cards** in respect to your practice.




* The ID-Titration Solution can be perfectly used for manual method and integrated with IH-Com.

Enjoy the IH-500 System, with a Complete Automated Titration Solution.



Seamless Integration in your Lab

- Automatically adapted to your laboratory practice, the IH-500 titration procedure allows you to use the Bio-Rad ID-Cells and your customized red blood cells.
- **Fully integrated with IH-Com**, the titration result is added to the overall patient profile allowing you to provide **better care, faster timings** in an highly reliable system.



Main Features

- **Automated antibody titration with two fold dilutions from:**
 - 1:1 to 1:32
 - 1:64 to 1:2048
 - 1:1 to 1:2048
- **Use any ID-Cell for the titration on the IH-500 according to the conclusion of the antibody identification**
- **Use the ID-Card you need:**
 - LISS / Coombs (Id-n°: 50531)
 - Coombs Anti-IgG (Id-n°: 50540)
 - NaCl, Enzyme Test and Cold Agglutinins (Id-n°: 50520)
- **ID-Titration Solution: a titration buffer designed to match with the ID-System**
- **Fully integrated with IH-Com**

Ordering Information

Catalog #	Description
001500	IH-500 , Fully Automated System for ID-Cards
009000	IH-Com Kit Fullversion , Data Management Software
009360	ID-Titration Solution (Id-n°: 05780), 10 x 10 ml
009390	ID-Titration Racks (Id-n°: 05770), 10 racks of 60 wells

Visit www.bio-rad.com/immunohematology for more information.



BIO-RAD

**Bio-Rad
Laboratories**

*Clinical
Diagnostics Group*

Website www.bio-rad.com/diagnostics **Australia** +61(2)9914-2800 **Austria** +43-1-877-8901 **Belgium** +32(3)710-53-00 **Brazil** +55(31)3689-6600
Canada +1 514 334-4372 **China** +86-21-61698500 **Czech Republic** +420-241-430-532 **Denmark** +45-4452-1000 **Finland** +358-9-804-22-00
France +33 1 47 95 60 00 **Germany** +49(0)89-318-840 **Greece** +30-210-7774396 **Hong Kong** +852-2789-3300 **Hungary** +36-1-459-6100
India +1-800-180-1224 **Israel** +972-3-9636050 **Italy** +39-02-216091 **Japan** +81-3-6361-7070 **Korea** +82-2-3473-4460 **Mexico** +52(55)5488-7670
The Netherlands +31-318-540666 **New Zealand** +64-9-415-2280 **Norway** +47-23-38-41-30 **Poland** +48-22-3319999 **Portugal** +351-21-472-7700
Russia +7-495-721-14-04 **Singapore** +65-6415-3170 **South Africa** +27-11-442-85-08 **Spain** +34-91-590-5200 **Sweden** +46-8-555-127-00
Switzerland +41 26 674 55 05/06 **Taiwan** +886-2-2578-7189 **Thailand** +662-651-8311 **United Kingdom** +44(0)20-8328-2000