

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4100455 -7 Pseudomonas cepacia	99		4100572 -3 Pseudomonas cepacia	99		4100705 -7 Pseudomonas cepacia	99	
4100461 -3 Pseudomonas cepacia	99		4100575 -7 Pseudomonas cepacia	99		4100711 -3 Pseudomonas cepacia	99	
4100465 -7 Pseudomonas cepacia	99		4100601 -3 Pseudomonas cepacia	99		4100715 -7 Pseudomonas cepacia	99	
4100473 Pseudomonas cepacia	99		4100605 -7 Pseudomonas cepacia	99		4100721 -3 Pseudomonas cepacia	99	
4100476 -7 Pseudomonas cepacia	99		4100611 -3 Pseudomonas cepacia	99		4100725 -7 Pseudomonas cepacia	99	
4100501 -3 Pseudomonas cepacia	99		4100615 -7 Pseudomonas cepacia	99		4100731 -3 Pseudomonas cepacia	99	
4100505 -7 Pseudomonas cepacia	99		4100621 -3 Pseudomonas cepacia	99		4100735 -7 Pseudomonas cepacia	99	
4100511 -3 Pseudomonas cepacia	99		4100625 -7 Pseudomonas cepacia	99		4100741 -3 Pseudomonas cepacia	99	
4100515 -7 Pseudomonas cepacia	99		4100631 -3 Pseudomonas cepacia	99		4100745 -7 Pseudomonas cepacia	99	
4100521 -3 Pseudomonas cepacia	99		4100635 -7 Pseudomonas cepacia	99		4100751 -3 Pseudomonas cepacia	99	
4100525 -7 Pseudomonas cepacia	99		4100641 -3 Pseudomonas cepacia	99		4100755 -7 Pseudomonas cepacia	99	
4100531 -3 Pseudomonas cepacia	99		4100645 -7 Pseudomonas cepacia	99		4100761 -3 Pseudomonas cepacia	99	
4100535 -7 Pseudomonas cepacia	99		4100651 -3 Pseudomonas cepacia	99		4100765 -7 Pseudomonas cepacia	99	
4100541 -3 Pseudomonas cepacia	99		4100655 -7 Pseudomonas cepacia	99		4100771 -3 Pseudomonas cepacia	99	
4100545 -7 Pseudomonas cepacia	99		4100661 -3 Pseudomonas cepacia	99		4100775 -7 Pseudomonas cepacia	99	
4100551 -3 Pseudomonas cepacia	99		4100665 -7 Pseudomonas cepacia	99		4101003 Pseudomonas cepacia	99	
4100555 -7 Pseudomonas cepacia	99		4100672 -3 Pseudomonas cepacia	99		4101007 Pseudomonas cepacia	99	
4100561 -3 Pseudomonas cepacia	99		4100675 -7 Pseudomonas cepacia	99		4101010 Weeksella virosa	99	
4100565 -7 Pseudomonas cepacia	99		4100701 -3 Pseudomonas cepacia	99		4101013 Pseudomonas cepacia	99	

4100455

4101013

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4101373 Pseudomonas cepacia	99		4101507 Pseudomonas cepacia	99		4101617 Pseudomonas cepacia	99	
4101377 Pseudomonas cepacia	99		4101513 Pseudomonas cepacia	99		4101623 Pseudomonas cepacia	99	
4101403 Pseudomonas cepacia	99		4101517 Pseudomonas cepacia	99		4101627 Pseudomonas cepacia	99	
4101407 Pseudomonas cepacia	99		4101523 Pseudomonas cepacia	99		4101633 Pseudomonas cepacia	99	
4101413 Pseudomonas cepacia	99		4101527 Pseudomonas cepacia	99		4101637 Pseudomonas cepacia	99	
4101417 Pseudomonas cepacia	99		4101533 Pseudomonas cepacia	99		4101643 Pseudomonas cepacia	99	
4101423 Pseudomonas cepacia	99		4101537 Pseudomonas cepacia	99		4101647 Pseudomonas cepacia	99	
4101427 Pseudomonas cepacia	99		4101543 Pseudomonas cepacia	99		4101653 Pseudomonas cepacia	99	
4101433 Pseudomonas cepacia	99		4101547 Pseudomonas cepacia	99		4101657 Pseudomonas cepacia	99	
4101437 Pseudomonas cepacia	99		4101550 Flavobacterium menigosep.	99		4101663 Pseudomonas cepacia	99	
4101443 Pseudomonas cepacia	99		4101553 Pseudomonas cepacia	99		4101667 Pseudomonas cepacia	99	
4101447 Pseudomonas cepacia	99		4101557 Pseudomonas cepacia	99		4101673 Pseudomonas cepacia	99	
4101453 Pseudomonas cepacia	99		4101563 Pseudomonas cepacia	99		4101677 Pseudomonas cepacia	99	
4101457 Pseudomonas cepacia	99		4101567 Pseudomonas cepacia	99		4101703 Pseudomonas cepacia	99	
4101463 Pseudomonas cepacia	99		4101573 Pseudomonas cepacia	99		4101707 Pseudomonas cepacia	99	
4101467 Pseudomonas cepacia	99		4101577 Pseudomonas cepacia	99		4101713 Pseudomonas cepacia	99	
4101473 Pseudomonas cepacia	99		4101603 Pseudomonas cepacia	99		4101717 Pseudomonas cepacia	99	
4101477 Pseudomonas cepacia	99		4101607 Pseudomonas cepacia	99		4101723 Pseudomonas cepacia	99	
4101503 Pseudomonas cepacia	99		4101613 Pseudomonas cepacia	99		4101727 Pseudomonas cepacia	99	
4101373						4101727		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4101733 Pseudomonas cepacia	99		4102047 Pseudomonas cepacia	99		4102163 Pseudomonas cepacia	99	
4101737 Pseudomonas cepacia	99		4102053 Pseudomonas cepacia	99		4102167 Pseudomonas cepacia	99	
4101743 Pseudomonas cepacia	99		4102057 Pseudomonas cepacia	99		4102173 Pseudomonas cepacia	99	
4101747 Pseudomonas cepacia	99		4102063 Pseudomonas cepacia	99		4102177 Pseudomonas cepacia	99	
4101753 Pseudomonas cepacia	99		4102067 Pseudomonas cepacia	99		4102203 Alcaligenes species Pseudomonas cepacia	50 OFFR/ 0 49 OFFR/100	
4101757 Pseudomonas cepacia	99		4102073 Pseudomonas cepacia	99		4102206 Pseudomonas cepacia	99	
4101763 Pseudomonas cepacia	99		4102077 Pseudomonas cepacia	99		4102207 Pseudomonas cepacia Alcaligenes species	64 OFFR/100 35 OFFR/ 0	
4101767 Pseudomonas cepacia	99		4102103 Pseudomonas cepacia	99		4102213 Pseudomonas cepacia	99	
4101773 Pseudomonas cepacia	99		4102106 -7 Pseudomonas cepacia	99		4102217 Pseudomonas cepacia	99	
4101777 Pseudomonas cepacia	99		4102113 Pseudomonas cepacia	99		4102223 Pseudomonas cepacia	99	
4102003 Pseudomonas cepacia	99		4102117 Pseudomonas cepacia	99		4102227 Pseudomonas cepacia	99	
4102007 Pseudomonas cepacia	99		4102123 Pseudomonas cepacia	99		4102233 Pseudomonas cepacia	99	
4102013 Pseudomonas cepacia	99		4102127 Pseudomonas cepacia	99		4102237 Pseudomonas cepacia	99	
4102017 Pseudomonas cepacia	99		4102133 Pseudomonas cepacia	99		4102243 Pseudomonas cepacia	99	
4102023 Pseudomonas cepacia	99		4102137 Pseudomonas cepacia	99		4102247 Pseudomonas cepacia	99	
4102027 Pseudomonas cepacia	99		4102143 Pseudomonas cepacia	99		4102253 Pseudomonas cepacia	99	
4102033 Pseudomonas cepacia	99		4102147 Pseudomonas cepacia	99		4102257 Pseudomonas cepacia	99	
4102037 Pseudomonas cepacia	99		4102153 Pseudomonas cepacia	99		4102263 Pseudomonas cepacia	99	
4102043 Pseudomonas cepacia	99		4102157 Pseudomonas cepacia	99				
4101733						4102263		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4102267 Pseudomonas cepacia	99		4102403 Pseudomonas cepacia	99		4102517 Pseudomonas cepacia	99	
4102273 Pseudomonas cepacia	99		4102406 -7 Pseudomonas cepacia	99		4102523 Pseudomonas cepacia	99	
4102277 Pseudomonas cepacia	99		4102413 Pseudomonas cepacia	99		4102526 -7 Pseudomonas cepacia	99	
4102302 -3 Pseudomonas cepacia	99		4102417 Pseudomonas cepacia	99		4102533 Pseudomonas cepacia	99	
4102305 -7 Pseudomonas cepacia	99		4102423 Pseudomonas cepacia	99		4102537 Pseudomonas cepacia	99	
4102313 Pseudomonas cepacia	99		4102427 Pseudomonas cepacia	99		4102543 Pseudomonas cepacia	99	
4102317 Pseudomonas cepacia	99		4102433 Pseudomonas cepacia	99		4102546 -7 Pseudomonas cepacia	99	
4102323 Pseudomonas cepacia	99		4102437 Pseudomonas cepacia	99		4102553 Pseudomonas cepacia	99	
4102326 -7 Pseudomonas cepacia	99		4102443 Pseudomonas cepacia	99		4102557 Pseudomonas cepacia	99	
4102333 Pseudomonas cepacia	99		4102447 Pseudomonas cepacia	99		4102563 Pseudomonas cepacia	99	
4102337 Pseudomonas cepacia	99		4102453 Pseudomonas cepacia	99		4102566 -7 Pseudomonas cepacia	99	
4102343 Pseudomonas cepacia	99		4102457 Pseudomonas cepacia	99		4102573 Pseudomonas cepacia	99	
4102346 -7 Pseudomonas cepacia	99		4102463 Pseudomonas cepacia	99		4102577 Pseudomonas cepacia	99	
4102353 Pseudomonas cepacia	99		4102467 Pseudomonas cepacia	99		4102602 -3 Pseudomonas cepacia	99	
4102357 Pseudomonas cepacia	99		4102473 Pseudomonas cepacia	99		4102605 -7 Pseudomonas cepacia	99	
4102363 Pseudomonas cepacia	99		4102477 Pseudomonas cepacia	99		4102613 Pseudomonas cepacia	99	
4102366 -7 Pseudomonas cepacia	99		4102502 -3 Pseudomonas cepacia	99		4102616 -7 Pseudomonas cepacia	99	
4102373 Pseudomonas cepacia	99		4102505 -7 Pseudomonas cepacia	99		4102622 -3 Pseudomonas cepacia	99	
4102377 Pseudomonas cepacia	99		4102513 Pseudomonas cepacia	99		4102626 -7 Pseudomonas cepacia	99	
4102267						4102627		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4103527 Pseudomonas cepacia	99		4103707 Pseudomonas cepacia	99		4104013 Pseudomonas cepacia	99	
4103537 Pseudomonas cepacia	99		4103713 Pseudomonas cepacia	99		4104017 Pseudomonas cepacia	99	
4103543 Pseudomonas cepacia	99		4103717 Pseudomonas cepacia	99		4104023 Pseudomonas cepacia	99	
4103547 Pseudomonas cepacia	99		4103723 Pseudomonas cepacia	99		4104026 -7 Pseudomonas cepacia	99	
4103557 Pseudomonas cepacia	99		4103727 Pseudomonas cepacia	99		4104033 Pseudomonas cepacia	99	
4103563 Pseudomonas cepacia	99		4103733 Pseudomonas cepacia	99		4104037 Pseudomonas cepacia	99	
4103567 Pseudomonas cepacia	99		4103737 Pseudomonas cepacia	99		4104043 Pseudomonas cepacia	99	
4103603 Pseudomonas cepacia	99		4103743 Pseudomonas cepacia	99		4104046 -7 Pseudomonas cepacia	99	
4103607 Pseudomonas cepacia	99		4103747 Pseudomonas cepacia	99		4104053 Pseudomonas cepacia	99	
4103617 Pseudomonas cepacia	99		4103753 Pseudomonas cepacia	99		4104057 Pseudomonas cepacia	99	
4103623 Pseudomonas cepacia	99		4103757 Pseudomonas cepacia	99		4104063 Pseudomonas cepacia	99	
4103627 Pseudomonas cepacia	99		4103763 Pseudomonas cepacia	99		4104067 Pseudomonas cepacia	99	
4103637 Pseudomonas cepacia	99		4103767 Pseudomonas cepacia	99		4104073 Pseudomonas cepacia	99	
4103643 Pseudomonas cepacia	99		4103777 Pseudomonas cepacia	99		4104077 Pseudomonas cepacia	99	
4103647 Pseudomonas cepacia	99		4104003 Pseudomonas cepacia	58	OFFR/100	4104102 -3 Pseudomonas cepacia	99	
4103657 Pseudomonas cepacia	99		Alcaligenes species	41	OFFR/ 0	4104105 -7 Pseudomonas cepacia	99	
4103663 Pseudomonas cepacia	99		4104006 Pseudomonas cepacia	99		4104113 Pseudomonas cepacia	99	
4103667 Pseudomonas cepacia	99		4104007 Achromobacter xylosoxid.	57	GRSS/100	4104116 -7 Pseudomonas cepacia	99	
4103703 Pseudomonas cepacia	99		Pseudomonas cepacia	31	GRSS/ 9			
			Alcaligenes species	12				
			4104012 Flavobacterium odoratum	99		4104122 -3 Pseudomonas cepacia	99	
4103527						4104123		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4104126 -7 Pseudomonas cepacia	99		4104217 Pseudomonas cepacia Alcaligenes species	62	OFFR/100 37 OFFR/ 0	4104316 -7 Pseudomonas cepacia	99	
4104133 Pseudomonas cepacia	99		4104222 -3 Pseudomonas cepacia	99		4104321 -3 Pseudomonas cepacia	99	
4104137 Pseudomonas cepacia	99		4104225 -6 Pseudomonas cepacia	99		4104325 -7 Pseudomonas cepacia	99	
4104142 -3 Pseudomonas cepacia	99		4104227 Pseudomonas cepacia	61	GRSS/ 9	4104333 Pseudomonas cepacia	99	
4104146 -7 Pseudomonas cepacia	99		Achromobacter xylosoxid. 4104233 Pseudoonas cepacia	38	GRSS/100	4104336 -7 Pseudomonas cepacia	99	
4104153 Pseudomonas cepacia	99		4104237 Pseudomonas cepacia	99		4104341 -3 Pseudomonas cepacia	99	
4104157 Pseudomonas cepacia	99		4104242 -3 Pseudomonas cepacia	99		4104345 -7 Pseudomonas cepacia	99	
4104163 Pseudomonas cepacia	99		4104245 -7 Pseudomonas cepacia	99		4104353 Pseudomonas cepacia	99	
4104166 -7 Pseudomonas cepacia	99		4104253 Pseudomonas cepacia	99		4104356 -7 Pseudomonas cepacia	99	
4104173 Pseudomonas cepacia	99		4104257 Pseudomonas cepacia	99		4104362 -3 Pseudomonas cepacia	99	
4104177 Pseudomonas cepacia	99		4104263 Pseudomonas cepacia	99		4104365 -7 Pseudomonas cepacia	99	
4104202 Pseudomonas cepacia	99		4104266 -7 Pseudomonas cepacia	99		4104373 Pseudomonas cepacia	99	
4104203 Alcaligenes species Pseudomonas cepacia	54 45	OFFR/ 0 OFFR/100	4104273 Pseudomonas cepacia	99		4104376 -7 Pseudomonas cepacia	99	
4104205 -6 Pseudomonas cepacia	99		4104277 Pseudomonas cepacia	99		4104402 -3 Pseudomonas cepacia	99	
4104207 Achromobacter xylosoxid. Pseudomonas cepacia Alcaligenes species	72 17 11	GRSS/100 GRSS/ 9	4104301 -3 Pseudomonas cepacia	99		4104405 -7 Pseudomonas cepacia	99	
4104213 Alcaligenes species Pseudomonas cepacia	51 48	OFFR/ 0 OFFR/100	4104305 -6 Pseudomonas cepacia	99		4104413 Pseudomonas cepacia	99	
4104216 Pseudomonas cepacia	99		4104307 Pseudomonas cepacia	57	GRSS/ 9	4104416 -7 Pseudomonas cepacia	99	
			Achromobacter xylosoxid.	42	GRSS/100	4104422 -3 Pseudomonas cepacia	99	
			4104312 -3 Pseudomonas cepacia	99		4104425 -7 Pseudomonas cepacia	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4106023 Pseudomonas cepacia	99		4106137 Pseudomonas cepacia	99		4106233 Pseudomonas cepacia	99	
4106027 Pseudomonas cepacia	99		4106143 Pseudomonas cepacia	99		4106237 Pseudomonas cepacia	99	
4106033 Pseudomonas cepacia	99		4106147 Pseudomonas cepacia	99		4106243 Pseudomonas cepacia	99	
4106037 Pseudomonas cepacia	99		4106153 Pseudomonas cepacia	99		4106247 Pseudomonas cepacia	99	
4106043 Pseudomonas cepacia	99		4106157 Pseudomonas cepacia	99		4106253 Pseudomonas cepacia	99	
4106047 Pseudomonas cepacia	99		4106163 Pseudomonas cepacia	99		4106257 Pseudomonas cepacia	99	
4106053 Pseudomonas cepacia	99		4106167 Pseudomonas cepacia	99		4106263 Pseudomonas cepacia	99	
4106057 Pseudomonas cepacia	99		4106173 Pseudomonas cepacia	99		4106267 Pseudomonas cepacia	99	
4106063 Pseudomonas cepacia	99		4106177 Pseudomonas cepacia	99		4106273 Pseudomonas cepacia	99	
4106067 Pseudomonas cepacia	99		4106202 Alcaligenes species	99		4106277 Pseudomonas cepacia	99	
4106073 Pseudomonas cepacia	99		4106203 Alcaligenes species Pseudomonas cepacia	92 7		4106303 Pseudomonas cepacia	99	
4106077 Pseudomonas cepacia	99		4106206 Alcaligenes species	99		4106307 Pseudomonas cepacia	99	
4106103 Pseudomonas cepacia	99		4106207 Alcaligenes species	59		4106313 Pseudomonas cepacia	99	
4106107 Pseudomonas cepacia	99		Achromobacter xylooxid. Pseudomonas cepacia	32 8		4106317 Pseudomonas cepacia	99	
4106113 Pseudomonas cepacia	99		4106213 Alcaligenes species Pseudomonas cepacia	92 7		4106323 Pseudomonas cepacia	99	
4106117 Pseudomonas cepacia	99		4106217 Alcaligenes species Pseudomonas cepacia	86 13		4106327 Pseudomonas cepacia	99	
4106123 Pseudomonas cepacia	99		4106223 Pseudomonas cepacia	99		4106333 Pseudomonas cepacia	99	
4106127 Pseudomonas cepacia	99		4106227 Pseudomonas cepacia	99		4106337 Pseudomonas cepacia	99	
4106133 Pseudomonas cepacia	99					4106343 Pseudomonas cepacia	99	
4106023						4106343		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4106713 Pseudomonas cepacia	99		4107547 Pseudomonas cepacia	99		4120007 Pseudomonas cepacia	99	
4106717 Pseudomonas cepacia	99		4107607 Pseudomonas cepacia	99		4120013 Pseudomonas cepacia	99	
4106723 Pseudomonas cepacia	99		4107627 Pseudomonas cepacia	99		4120017 Pseudomonas cepacia	99	
4106726 -7 Pseudomonas cepacia	99		4107647 Pseudomonas cepacia	99		4120023 Pseudomonas cepacia	99	
4106733 Pseudomonas cepacia	99		4107703 Pseudomonas cepacia	99		4120027 Pseudomonas cepacia	99	
4106737 Pseudomonas cepacia	99		4107707 Pseudomonas cepacia	99		4120033 Pseudomonas cepacia	99	
4106743 Pseudomonas cepacia	99		4107717 Pseudomonas cepacia	99		4120037 Pseudomonas cepacia	99	
4106746 -7 Pseudomonas cepacia	99		4107723 Pseudomonas cepacia	99		4120043 Pseudomonas cepacia	99	
4106753 Pseudomonas cepacia	99		4107727 Pseudomonas cepacia	99		4120047 Pseudomonas cepacia	99	
4106757 Pseudomonas cepacia	99		4107743 Pseudomonas cepacia	99		4120053 Pseudomonas cepacia	99	
4106763 Pseudomonas cepacia	99		4107747 Pseudomonas cepacia	99		4120057 Pseudomonas cepacia	99	
4106766 -7 Pseudomonas cepacia	99		4107767 Pseudomonas cepacia	99		4120063 Pseudomonas cepacia	99	
4106773 Pseudomonas cepacia	99		4110707 Pseudomonas cepacia	99		4120067 Pseudomonas cepacia	99	
4106777 Pseudomonas cepacia	99		4114007 Achromobacter species	99		4120073 Pseudomonas cepacia	99	
4107307 Pseudomonas cepacia	99		4114012 Flavobacterium odoratum	99		4120077 Pseudomonas cepacia	99	
4107327 Pseudomonas cepacia	99		4114207 Bordetella bronchiseptica	99		4120103 Pseudomonas cepacia	99	
4107347 Pseudomonas cepacia	99		4115110 Flavobacterium species	99		4120107 Pseudomonas cepacia	99	
4107507 Pseudomonas cepacia	99		4115510 Flavobacterium species	99		4120113 Pseudomonas cepacia	99	
4107527 Pseudomonas cepacia	99		4120003 Pseudomonas cepacia	99		4120117 Pseudomonas cepacia	99	
4106713						4120117		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4120123 Pseudomonas cepacia	99		4120237 Pseudomonas cepacia	99		4120353 Pseudomonas cepacia	99	
4120127 Pseudomonas cepacia	99		4120243 Pseudomonas cepacia	99		4120357 Pseudomonas cepacia	99	
4120133 Pseudomonas cepacia	99		4120247 Pseudomonas cepacia	99		4120363 Pseudomonas cepacia	99	
4120137 Pseudomonas cepacia	99		4120253 Pseudomonas cepacia	99		4120367 Pseudomonas cepacia	99	
4120143 Pseudomonas cepacia	99		4120257 Pseudomonas cepacia	99		4120373 Pseudomonas cepacia	99	
4120147 Pseudomonas cepacia	99		4120263 Pseudomonas cepacia	99		4120377 Pseudomonas cepacia	99	
4120153 Pseudomonas cepacia	99		4120267 Pseudomonas cepacia	99		4120403 Pseudomonas cepacia	99	
4120157 Pseudomonas cepacia	99		4120273 Pseudomonas cepacia	99		4120407 Pseudomonas cepacia	99	
4120163 Pseudomonas cepacia	99		4120277 Pseudomonas cepacia	99		4120413 Pseudomonas cepacia	99	
4120167 Pseudomonas cepacia	99		4120303 Pseudomonas cepacia	99		4120417 Pseudomonas cepacia	99	
4120173 Pseudomonas cepacia	99		4120306 -7 Pseudomonas cepacia	99		4120423 Pseudomonas cepacia	99	
4120177 Pseudomonas cepacia	99		4120313 Pseudomonas cepacia	99		4120427 Pseudomonas cepacia	99	
4120203 Pseudomonas cepacia	99		4120317 Pseudomonas cepacia	99		4120433 Pseudomonas cepacia	99	
4120207 Pseudomonas cepacia	99		4120323 Pseudomonas cepacia	99		4120437 Pseudomonas cepacia	99	
4120213 Pseudomonas cepacia	99		4120327 Pseudomonas cepacia	99		4120443 Pseudomonas cepacia	99	
4120217 Pseudomonas cepacia	99		4120333 Pseudomonas cepacia	99		4120447 Pseudomonas cepacia	99	
4120223 Pseudomonas cepacia	99		4120337 Pseudomonas cepacia	99		4120453 Pseudomonas cepacia	99	
4120227 Pseudomonas cepacia	99		4120343 Pseudomonas cepacia	99		4120457 Pseudomonas cepacia	99	
4120233 Pseudomonas cepacia	99		4120347 Pseudomonas cepacia	99		4120463 Pseudomonas cepacia	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4122243 Pseudomonas cepacia	99		4122357 Pseudomonas cepacia	99		4122473 Pseudomonas cepacia	99	
4122247 Pseudomonas cepacia	99		4122363 Pseudomonas cepacia	99		4122477 Pseudomonas cepacia	99	
4122253 Pseudomonas cepacia	99		4122367 Pseudomonas cepacia	99		4122503 Pseudomonas cepacia	99	
4122257 Pseudomonas cepacia	99		4122373 Pseudomonas cepacia	99		4122507 Pseudomonas cepacia	99	
4122263 Pseudomonas cepacia	99		4122377 Pseudomonas cepacia	99		4122513 Pseudomonas cepacia	99	
4122267 Pseudomonas cepacia	99		4122403 Pseudomonas cepacia	99		4122517 Pseudomonas cepacia	99	
4122273 Pseudomonas cepacia	99		4122407 Pseudomonas cepacia	99		4122523 Pseudomonas cepacia	99	
4122277 Pseudomonas cepacia	99		4122413 Pseudomonas cepacia	99		4122527 Pseudomonas cepacia	99	
4122303 Pseudomonas cepacia	99		4122417 Pseudomonas cepacia	99		4122533 Pseudomonas cepacia	99	
4122307 Pseudomonas cepacia	99		4122423 Pseudomonas cepacia	99		4122537 Pseudomonas cepacia	99	
4122313 Pseudomonas cepacia	99		4122427 Pseudomonas cepacia	99		4122543 Pseudomonas cepacia	99	
4122317 Pseudomonas cepacia	99		4122433 Pseudomonas cepacia	99		4122547 Pseudomonas cepacia	99	
4122323 Pseudomonas cepacia	99		4122437 Pseudomonas cepacia	99		4122553 Pseudomonas cepacia	99	
4122327 Pseudomonas cepacia	99		4122443 Pseudomonas cepacia	99		4122557 Pseudomonas cepacia	99	
4122333 Pseudomonas cepacia	99		4122447 Pseudomonas cepacia	99		4122563 Pseudomonas cepacia	99	
4122337 Pseudomonas cepacia	99		4122453 Pseudomonas cepacia	99		4122567 Pseudomonas cepacia	99	
4122343 Pseudomonas cepacia	99		4122457 Pseudomonas cepacia	99		4122573 Pseudomonas cepacia	99	
4122347 Pseudomonas cepacia	99		4122463 Pseudomonas cepacia	99		4122577 Pseudomonas cepacia	99	
4122353 Pseudomonas cepacia	99		4122467 Pseudomonas cepacia	99		4122603 Pseudomonas cepacia	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4126127 Pseudomonas cepacia	99		4126337 Pseudomonas cepacia	99		4126517 Pseudomonas cepacia	99	
4126143 Pseudomonas cepacia	99		4126343 Pseudomonas cepacia	99		4126523 Pseudomonas cepacia	99	
4126147 Pseudomonas cepacia	99		4126347 Pseudomonas cepacia	99		4126527 Pseudomonas cepacia	99	
4126167 Pseudomonas cepacia	99		4126357 Pseudomonas cepacia	99		4126537 Pseudomonas cepacia	99	
4126203 Pseudomonas cepacia	99		4126363 Pseudomonas cepacia	99		4126543 Pseudomonas cepacia	99	
4126207 Pseudomonas cepacia	99		4126367 Pseudomonas cepacia	99		4126547 Pseudomonas cepacia	99	
4126217 Pseudomonas cepacia	99		4126377 Pseudomonas cepacia	99		4126557 Pseudomonas cepacia	99	
4126223 Pseudomonas cepacia	99		4126403 Pseudomonas cepacia	99		4126563 Pseudomonas cepacia	99	
4126227 Pseudomonas cepacia	99		4126407 Pseudomonas cepacia	99		4126567 Pseudomonas cepacia	99	
4126243 Pseudomonas cepacia	99		4126417 Pseudomonas cepacia	99		4126577 Pseudomonas cepacia	99	
4126247 Pseudomonas cepacia	99		4126423 Pseudomonas cepacia	99		4126603 Pseudomonas cepacia	99	
4126263 Pseudomonas cepacia	99		4126427 Pseudomonas cepacia	99		4126607 Pseudomonas cepacia	99	
4126267 Pseudomonas cepacia	99		4126443 Pseudomonas cepacia	99		4126613 Pseudomonas cepacia	99	
4126303 Pseudomonas cepacia	99		4126447 Pseudomonas cepacia	99		4126617 Pseudomonas cepacia	99	
4126307 Pseudomonas cepacia	99		4126463 Pseudomonas cepacia	99		4126623 Pseudomonas cepacia	99	
4126313 Pseudomonas cepacia	99		4126467 Pseudomonas cepacia	99		4126627 Pseudomonas cepacia	99	
4126317 Pseudomonas cepacia	99		4126503 Pseudomonas cepacia	99		4126633 Pseudomonas cepacia	99	
4126323 Pseudomonas cepacia	99		4126507 Pseudomonas cepacia	99		4126637 Pseudomonas cepacia	99	
4126327 Pseudomonas cepacia	99		4126513 Pseudomonas cepacia	99		4126643 Pseudomonas cepacia	99	
4126127						4126643		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4126647 Pseudomonas cepacia	99		4126767 Pseudomonas cepacia	99		4146317 Pseudomonas aeruginosa	99	
4126653 Pseudomonas cepacia	99		4126773 Pseudomonas cepacia	99		4146337 Pseudomonas aeruginosa	99	
4126657 Pseudomonas cepacia	99		4126777 Pseudomonas cepacia	99		4147157 Aeromonas hydrophila	99	
4126663 Pseudomonas cepacia	99		4140557 Aeromonas hydrophila	99		4147557 Aeromonas hydrophila	99	
4126667 Pseudomonas cepacia	99		4140707 Pseudomonas cepacia	99		4161007 Plesiomonas shigelloides	99	
4126677 Pseudomonas cepacia	99		4141117 Aeromonas hydrophila	99		4161046 -7 Plesiomons shigelloides	99	
4126703 Pseudomonas cepacia	99		4141157 Aeromonas hydrophila	99		4161106 -7 Plesiomonas shigelloides	99	
4126707 Pseudomonas cepacia	99		4141517 Aeromonas hydrophila	99		4161146 -7 Plesiomonas shigelloides	99	
4126713 Pseudomonas cepacia	99		4141557 Aeromonas hydrophila	99		4161307 Plesiomonas shigelloides	99	
4126717 Pseudomonas cepacia	99		4143157 Aeromonas hydrophila	99		4161347 Plesiomonas shigelloides	99	
4126723 Pseudomonas cepacia	99		4143557 Aeromonas hydrophila	99		4200000 Moraxella species	99	
4126727 Pseudomonas cepacia	99		4144003 Pseudomonas putida	99		4200002 Oligella ureolytica	99	
4126733 Pseudomonas cepacia	99		4144103 Pseudomonas fluorescens	99		4200003 Alteromonas putrefaciens	56	
4126737 Pseudomonas cepacia	99		4144107 Pseudomonas stutzeri	99		Pseudomonas species	43	
4126743 Pseudomonas cepacia	99		4144557 Aeromonas hydrophila	99		4200007 Alteromonas putrefaciens	99	
4126747 Pseudomonas cepacia	99		4145117 Aeromonas hydrophila	99		4200013 Alteromonas putrefaciens	99	
4126753 Pseudomonas cepacia	99		4145157 Aeromonas hydrophila	99		4200017 Alteromonas putrefaciens	99	
4126757 Pseudomonas cepacia	99		4145517 Aeromonas hydrophila	99		4200047 Alteromonas putrefaciens	99	
4126763 Pseudomonas cepacia	99		4145557 Aeromonas hydrophila	99		4200057 Alteromonas putrefaciens	99	
4126647						4200107 Alteromonas putrefaciens	99	
						4200107		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4200117 Alteromonas putrefaciens	99		4204147 Alteromonas putrefaciens	99		4220015 -7 Alteromonas putrefaciens	99	
4201010 Weeksella virosa	99		4204157 Alteromonas putrefaciens	99		4220027 Alteromonas putrefaciens	99	
4201210 Weeksella virosa	99		4204207 Achromobacter xylosoxid.	99		4220037 Alteromonas putrefaciens	99	
4201550 Flavobacterium menigosep.	99		4205010 Weeksella virosa	99		4220043 Alteromonas putrefaciens	99	
4204003 Alteromonas putrefaciens	99		4205210 Weeksella virosa	99		4220047 Alteromonas putrefaciens	99	
4204005 -7 Alteromonas putrefaciens	99		4206203 Alcaligenes species	99		4220053 Alteromonas putrefaciens	99	
4204012 Flavobacterium odoratum	99		4214007 Alteromonas putrefaciens Achromobacter species	92 7		4220057 Alteromonas putrefaciens	99	
4204013 Alteromonas putrefaciens	99		4214012 Flavobacterium odoratum	99		4220103 Alteromonas putrefaciens	99	
4204015 -7 Alteromonas putrefaciens	99		4214017 Alteromonas putrefaciens	99		4220107 Alteromonas putrefaciens	99	
4204027 Alteromonas putrefaciens	99		4214047 Alteromonas putrefaciens	99		4220113 Alteromonas putrefaciens	99	
4204037 Alteromonas putrefaciens	99		4214057 Alteromonas putrefaciens	99		4220117 Alteromonas putrefaciens	99	
4204043 Alteromonas putrefaciens	99		4214107 Alteromonas putrefaciens	99		4220147 Alteromonas putrefaciens	99	
4204047 Alteromonas putrefaciens	99		4214117 Alteromonas putrefaciens	99		4220157 Alteromonas putrefaciens	99	
4204053 Alteromonas putrefaciens	99		4214207 Bordetella bronchiseptica	99		4224001 -3 Alteromonas putrefaciens	99	
4204057 Alteromonas putrefaciens	99		4215110 Flavobacterium species	99		4224005 -7 Alteromonas putrefaciens	99	
4204103 Alteromonas putrefaciens	99		4215510 Flavobacterium species	99		4224011 -3 Alteromonas putrefaciens	99	
4204107 Alteromonas putrefaciens	99		4220003 Alteromonas putrefaciens	99		4224015 -7 Alteromonas putrefaciens	99	
4204113 Alteromonas putrefaciens	99		4220005 -7 Alteromonas putrefaciens	99		4224023 Alteromonas putrefaciens	99	
4204117 Alteromonas putrefaciens	99		4220013 Alteromonas putrefaciens	99		4224027 Alteromonas putrefaciens	99	
4200117						4224027		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4400557 Pseudomonas cepacia	99		4400663 Pseudomonas cepacia	99		4400766 -7 Pseudomonas cepacia	99	
4400560 Flavobacterium meningosep.	99		4400667 Pseudomonas cepacia	99		4400773 Pseudomonas cepacia	99	
4400563 Pseudomonas cepacia	99		4400673 Pseudomonas cepacia	99		4400777 Pseudomonas cepacia	99	
4400567 Pseudomonas cepacia	99		4400677 Pseudomonas cepacia	99		4401000 Flavobacterium meningosep.	99	
4400570 Flavobacterium meningosep.	99		4400700 -1 Sphingomonas paucimobilis	99		4401002 Flavobacterium meningosep.	99	
4400573 Pseudomonas cepacia	99		4400702 -3 Pseudomonas cepacia	99		4401010 Flavobacterium meningosep. Weeksella virosa	78 OFFR/100 21 OFFR/ 0	
4400577 Pseudomonas cepacia	99		4400705 -7 Pseudomonas cepacia	99		4401012 Flavobacterium meningosep.	99	
4400603 Pseudomonas cepacia	99		4400710 Flavobacterium species	99		4401020 Flavobacterium meningosep.	99	
4400606 -7 Pseudomonas cepacia	99		4400713 Pseudomonas cepacia	99		4401022 Flavobacterium meningosep.	99	
4400613 Pseudomonas cepacia	99		4400716 -7 Pseudomonas cepacia	99		4401030 Flavobacterium meningosep.	99	
4400617 Pseudomonas cepacia	99		4400722 -3 Pseudomonas cepacia	99		4401032 Flavobacterium meningosep.	99	
4400623 Pseudomonas cepacia	99		4400725 -7 Pseudomonas cepacia	99		4401040 Flavobacterium meningosep.	99	
4400626 -7 Pseudomonas cepacia	99		4400733 Pseudomonas cepacia	99		4401042 Flavobacterium meningosep.	99	
4400633 Pseudomonas cepacia	99		4400737 Pseudomonas cepacia	99		4401050 Flavobacterium meningosep.	99	
4400637 Pseudomonas cepacia	99		4400742 -3 Pseudomonas cepacia	99		4401052 Flavobacterium meningosep.	99	
4400643 Pseudomonas cepacia	99		4400745 -7 Pseudomonas cepacia	99		4401054 Flavobacterium meningosep.	99	
4400646 -7 Pseudomonas cepacia	99		4400753 Pseudomonas cepacia	99		4401060 Flavobacterium meningosep.	99	
4400653 Pseudomonas cepacia	99		4400757 Pseudomonas cepacia	99		4401062 Flavobacterium meningosep.	99	
4400657 Pseudomonas cepacia	99		4400763 Pseudomonas cepacia	99		4401070 Flavobacterium meningosep.	99	
4400557						4401070		

BIOTYPE NUMBER	RP	ADD.TEST	BIOTYPE NUMBER	RP	ADD.TEST	BIOTYPE NUMBER	RP	ADD.TEST
ORGANISM NAME	%	POS. %	ORGANISM NAME	%	POS. %	ORGANISM NAME	%	POS. %
4401406			4401454			4401512		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	99		Flavobacterium menigosep.	95	
						Flavobacterium species	4	
4401407			4401456			4401514		
Pseudomonas cepacia	99		Flavobacterium menigosep.	99		Flavobacterium menigosep.	91	
						Flavobacterium species	8	
4401410			4401460			4401516		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	99		Flavobacterium menigosep.	99	
4401412			4401462			4401520		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	99		Flavobacterium menigosep.	99	
4401414			4401464			4401522		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	99		Flavobacterium menigosep.	99	
4401416			4401466			4401524		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	99		Flavobacterium menigosep.	99	
4401420			4401470			4401526		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	99		Flavobacterium menigosep.	99	
4401422			4401472			4401527		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	99		Pseudomonas cepacia	99	
4401424			4401474			4401530		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	99		Flavobacterium menigosep.	99	
4401430			4401476			4401532		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	99		Flavobacterium menigosep.	99	
4401432			4401500			4401534		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	95		Flavobacterium menigosep.	99	
			Sphingomonas paucimobilis	4				
4401434			4401501			4401536		
Flavobacterium menigosep.	99		Sphingomonas paucimobilis	99		Flavobacterium menigosep.	99	
4401436			4401502			4401540		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	99		Flavobacterium menigosep.	99	
4401440			4401503			4401542		
Flavobacterium menigosep.	99		Pseudomonas cepacia	99		Flavobacterium menigosep.	99	
4401442			4401504			4401544		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	99		Flavobacterium menigosep.	99	
4401444			4401506			4401546		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	99		Flavobacterium menigosep.	99	
4401446			4401507			4401547		
Flavobacterium menigosep.	99		Pseudomonas cepacia	99		Pseudomonas cepacia	99	
4401450			4401510			4401550		
Flavobacterium menigosep.	99		Flavobacterium menigosep.	96		Flavobacterium menigosep.	99	
			Flavobacterium species	3				
4401452								
Flavobacterium menigosep.	99							

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4402557 Pseudomonas cepacia	99		4402673 Pseudomonas cepacia	99		4403510 Flavobacterium species	99	
4402563 Pseudomonas cepacia	99		4402677 Pseudomonas cepacia	99		4404003 Pseudomonas cepacia	99	
4402567 Pseudomonas cepacia	99		4402703 Pseudomonas cepacia	99		4404007 Pseudomonas cepacia	99	
4402573 Pseudomonas cepacia	99		4402707 Pseudomonas cepacia	99		4404012 Flavobacterium odoratum	99	
4402577 Pseudomonas cepacia	99		4402713 Pseudomonas cepacia	99		4404013 Pseudomonas cepacia	99	
4402603 Pseudomonas cepacia	99		4402717 Pseudomonas cepacia	99		4404017 Pseudomonas cepacia	99	
4402607 Pseudomonas cepacia	99		4402723 Pseudomonas cepacia	99		4404023 Pseudomonas cepacia	99	
4402613 Pseudomonas cepacia	99		4402727 Pseudomonas cepacia	99		4404027 Pseudomonas cepacia	99	
4402617 Pseudomonas cepacia	99		4402733 Pseudomonas cepacia	99		4404033 Pseudomonas cepacia	99	
4402623 Pseudomonas cepacia	99		4402737 Pseudomonas cepacia	99		4404037 Pseudomonas cepacia	99	
4402627 Pseudomonas cepacia	99		4402743 Pseudomonas cepacia	99		4404043 Pseudomonas cepacia	99	
4402633 Pseudomonas cepacia	99		4402747 Pseudomonas cepacia	99		4404047 Pseudomonas cepacia	99	
4402637 Pseudomonas cepacia	99		4402753 Pseudomonas cepacia	99		4404053 Pseudomonas cepacia	99	
4402643 Pseudomonas cepacia	99		4402757 Pseudomonas cepacia	99		4404057 Pseudomonas cepacia	99	
4402647 Pseudomonas cepacia	99		4402763 Pseudomonas cepacia	99		4404063 Pseudomonas cepacia	99	
4402653 Pseudomonas cepacia	99		4402767 Pseudomonas cepacia	99		4404067 Pseudomonas cepacia	99	
4402657 Pseudomonas cepacia	99		4402773 Pseudomonas cepacia	99		4404073 Pseudomonas cepacia	99	
4402663 Pseudomonas cepacia	99		4402777 Pseudomonas cepacia	99		4404077 Pseudomonas cepacia	99	
4402667 Pseudomonas cepacia	99		4403110 Flavobacterium species	99		4404103 Pseudomonas cepacia	99	
4402557						4404103		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4404714 Flavobacterium species	99		4405012 Flavobacterium menigosep. Flavobacterium species	70 29	OFLA/ OFLA/	81 0	4405120 Flavobacterium menigosep.	99
4404717 Pseudomonas cepacia	99		4405020 Flavobacterium menigosep.	99		4405130 Flavobacterium menigosep. Flavobacterium species	68 31	OFLA/ OFLA/
4404723 Pseudomonas cepacia	99		4405030 Flavobacterium menigosep.	99		4405132 Flavobacterium menigosep. Flavobacterium species	63 36	OFLA/ OFLA/
4404727 Pseudomonas cepacia	99		4405032 Flavobacterium menigosep.	99		4405140 Flavobacterium menigosep.	99	
4404733 Pseudomonas cepacia	99		4405040 Flavobacterium menigosep.	99		4405142 Flavobacterium menigosep.	99	
4404737 Pseudomonas cepacia	99		4405042 Flavobacterium menigosep.	99		4405150 Flavobacterium menigosep. Flavobacterium species	83 16	
4404743 Pseudomonas cepacia	99		4405050 Flavobacterium menigosep.	99		4405152 Flavobacterium menigosep. Flavobacterium species	79 20	OFLA/ OFLA/
4404747 Pseudomonas cepacia	99		4405052 Flavobacterium menigosep.	99		4405154 Flavobacterium menigosep.	99	
4404750 Flavobacterium species	99		4405060 Flavobacterium menigosep.	99		4405160 Flavobacterium menigosep.	99	
4404753 Pseudomonas cepacia	99		4405062 Flavobacterium menigosep.	99		4405162 Flavobacterium menigosep.	99	
4404757 Pseudomonas cepacia	99		4405070 Flavobacterium menigosep.	99		4405170 Flavobacterium menigosep.	99	
4404763 Pseudomonas cepacia	99		4405072 Flavobacterium menigosep.	99		4405172 Flavobacterium menigosep.	99	
4404767 Pseudomonas cepacia	99		4405100 Flavobacterium menigosep.	99		4405210 Weeksella virosa Flavobacterium species Flavobacterium menigosep.	59 24 16	OFFR/ OFFR/ OFFR/
4404773 Pseudomonas cepacia	99		4405102 Flavobacterium menigosep.	99		4405240 Flavobacterium menigosep.	99	
4404777 Pseudomonas cepacia	99		4405110 Flavobacterium species Flavobacterium menigosep.	66 33	OFLA/ OFLA/	0 81	4405250 Flavobacterium menigosep.	99
4405000 Flavobacterium menigosep.	99		4405112 Flavobacterium species Flavobacterium menigosep.	71 28	OFLA/ OFLA/	0 81	4405252 Flavobacterium menigosep.	99
4405002 Flavobacterium menigosep.	99		4405114 Flavobacterium species	99		4405270 Flavobacterium menigosep.	99	
4405010 Flavobacterium menigosep. Weeksella virosa Flavobacterium species	50 32 17	OFFR/ OFFR/ 0	4405116 Flavobacterium species	99		4405270 Flavobacterium menigosep.	99	
4404714						4405270	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4405310 Flavobacterium species Flavobacterium menigosep.	89 10		4405422 Flavobacterium menigosep.	99		4405500 Flavobacterium menigosep.	99	
4405312 Flavobacterium species	99		4405430 Flavobacterium menigosep.	99		4405502 Flavobacterium menigosep.	99	
4405314 Flavobacterium species	99		4405432 Flavobacterium menigosep.	99		4405504 Flavobacterium menigosep.	99	
4405316 Flavobacterium species	99		4405434 Flavobacterium menigosep.	99		4405506 Flavobacterium menigosep.	99	
4405330 Flavobacterium species	99		4405440 Flavobacterium menigosep.	99		4405510 Flavobacterium menigosep. Flavobacterium species	81 18	
4405340 Flavobacterium menigosep.	99		4405442 Flavobacterium menigosep.	99		4405512 Flavobacterium menigosep. Flavobacterium species	78 OFLA/ 21 OFLA/	81 0
4405350 Flavobacterium menigosep. Flavobacterium species	51 OFLA/ 48 OFLA/	81 0	4405444 Flavobacterium menigosep.	99		4405514 Flavobacterium menigosep. Flavobacterium species	60 OFLA/ 39 OFLA/	81 0
4405352 Flavobacterium species Flavobacterium menigosep.	53 OFLA/ 46 OFLA/	0 81	4405446 Flavobacterium menigosep.	99		4405516 Flavobacterium menigosep. Flavobacterium species	54 OFLA/ 45 OFLA/	81 0
4405370 Flavobacterium menigosep.	99		4405450 Flavobacterium menigosep.	99		4405520 Flavobacterium menigosep.	99	
4405400 Flavobacterium menigosep.	99		4405452 Flavobacterium menigosep.	99		4405522 Flavobacterium menigosep.	99	
4405402 Flavobacterium menigosep.	99		4405454 Flavobacterium menigosep.	99		4405524 Flavobacterium menigosep.	99	
4405404 Flavobacterium menigosep.	99		4405456 Flavobacterium menigosep.	99		4405524 Flavobacterium menigosep.	99	
4405410 Flavobacterium menigosep. Flavobacterium species	96 4		4405460 Flavobacterium menigosep.	99		4405530 Flavobacterium menigosep. Flavobacterium species	94 5	
4405412 Flavobacterium menigosep. Flavobacterium species	95 4		4405462 Flavobacterium menigosep.	99		4405532 Flavobacterium menigosep. Flavobacterium species	93 6	
4405414 Flavobacterium menigosep.	99		4405464 Flavobacterium menigosep.	99		4405534 Flavobacterium menigosep.	99	
4405416 Flavobacterium menigosep.	99		4405470 Flavobacterium menigosep.	99		4405536 Flavobacterium menigosep.	99	
4405420 Flavobacterium menigosep.	99		4405472 Flavobacterium menigosep.	99		4405540 Flavobacterium menigosep.	99	
			4405474 Flavobacterium menigosep.	99		4405542 Flavobacterium menigosep.	99	
			4405476 Flavobacterium menigosep.	99				
4405310						4405542		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4405767 Pseudomonas cepacia	99		4406163 Pseudomonas cepacia	99		4406333 Pseudomonas cepacia	99	
4405770 Flavobacterium menigosep.	99		4406167 Pseudomonas cepacia	99		4406337 Pseudomonas cepacia	99	
4405772 Flavobacterium menigosep.	99		4406203 Alcaligenes species Pseudomonas cepacia	90 9		4406343 Pseudomonas cepacia	99	
4405774 Flavobacterium menigosep.	99		4406207 Pseudomonas cepacia	99		4406347 Pseudomonas cepacia	99	
4406003 Pseudomonas cepacia	99		4406217 Pseudomonas cepacia	99		4406353 Pseudomonas cepacia	99	
4406007 Pseudomonas cepacia	99		4406223 Pseudomonas cepacia	99		4406357 Pseudomonas cepacia	99	
4406027 Pseudomonas cepacia	99		4406227 Pseudomonas cepacia	99		4406363 Pseudomonas cepacia	99	
4406047 Pseudomonas cepacia	99		4406237 Pseudomonas cepacia	99		4406367 Pseudomonas cepacia	99	
4406067 Pseudomonas cepacia	99		4406243 Pseudomonas cepacia	99		4406377 Pseudomonas cepacia	99	
4406103 Pseudomonas cepacia	99		4406247 Pseudomonas cepacia	99		4406403 Pseudomonas cepacia	99	
4406107 Pseudomonas cepacia	99		4406257 Pseudomonas cepacia	99		4406407 Pseudomonas cepacia	99	
4406110 Flavobacterium species	99		4406263 Pseudomonas cepacia	99		4406417 Pseudomonas cepacia	99	
4406117 Pseudomonas cepacia	99		4406267 Pseudomonas cepacia	99		4406423 Pseudomonas cepacia	99	
4406123 Pseudomonas cepacia	99		4406303 Pseudomonas cepacia	99		4406427 Pseudomonas cepacia	99	
4406127 Pseudomonas cepacia	99		4406307 Pseudomonas cepacia	99		4406437 Pseudomonas cepacia	99	
4406137 Pseudomonas cepacia	99		4406313 Pseudomonas cepacia	99		4406443 Pseudomonas cepacia	99	
4406143 Pseudomonas cepacia	99		4406317 Pseudomonas cepacia	99		4406447 Pseudomonas cepacia	99	
4406147 Pseudomonas cepacia	99		4406323 Pseudomonas cepacia	99		4406457 Pseudomonas cepacia	99	
4406157 Pseudomonas cepacia	99		4406327 Pseudomonas cepacia	99		4406463 Pseudomonas cepacia	99	
4405767						4406463		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4410110 Flavobacterium species	99		4410514 Flavobacterium species	99		4411130 Flavobacterium species	99	
4410112 Flavobacterium species	99		4410516 Flavobacterium species	99		4411132 Flavobacterium species	99	
4410114 Flavobacterium species	99		4410530 Flavobacterium species	99		4411150 Flavobacterium species	99	
4410116 Flavobacterium species	99		4410550 Flavobacterium species	99		4411152 Flavobacterium species	99	
4410130 Flavobacterium species	99		4410552 Flavobacterium species	99		4411154 Flavobacterium species	99	
4410150 Flavobacterium species	99		4410610 Flavobacterium species	99		4411210 Flavobacterium species	99	
4410152 Flavobacterium species	99		4410710 Flavobacterium species	99		4411212 Flavobacterium species	99	
4410210 Flavobacterium species	99		4410712 Flavobacterium species	99		4411310 Flavobacterium species	99	
4410310 Flavobacterium species	99		4410714 Flavobacterium species	99		4411312 Flavobacterium species	99	
4410312 Flavobacterium species	99		4410716 Flavobacterium species	99		4411314 Flavobacterium species	99	
4410314 Flavobacterium species	99		4410730 Flavobacterium species	99		4411316 Flavobacterium species	99	
4410316 Flavobacterium species	99		4410750 Flavobacterium species	99		4411330 Flavobacterium species	99	
4410330 Flavobacterium species	99		4410752 Flavobacterium species	99		4411332 Flavobacterium species	99	
4410350 Flavobacterium species	99		4411010 Flavobacterium species	99		4411350 Flavobacterium species	99	
4410352 Flavobacterium species	99		4411012 Flavobacterium species	99		4411352 Flavobacterium species	99	
4410410 Flavobacterium species	99		4411110 Flavobacterium species	99		4411354 Flavobacterium species	99	
4410500 -1 Sphingomonas paucimobilis	99		4411112 Flavobacterium species	99		4411410 Flavobacterium species	99	
4410510 Flavobacterium species	99		4411114 Flavobacterium species	99		4411412 Flavobacterium species	99	
4410512 Flavobacterium species	99		4411116 Flavobacterium species	99		4411510 Flavobacterium species	99	

4410110

4411510

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4411512 Flavobacterium species	99		4412110 Flavobacterium species	99		4413750 Flavobacterium species	99	
4411514 Flavobacterium species	99		4412112 Flavobacterium species	99		4414007 Achromobacter species	99	
4411516 Flavobacterium species	99		4412310 Flavobacterium species	99		4414010 Flavobacterium species	99	
4411530 Flavobacterium species	99		4412510 Flavobacterium species	99		4414012 Flavobacterium odoratum Flavobacterium species	73 DASE/100 26 DASE/ 0	
4411532 Flavobacterium species	99		4412512 Flavobacterium species	99		4414110 Flavobacterium species	99	
4411550 Flavobacterium species	99		4412710 Flavobacterium species	99		4414112 Flavobacterium species	99	
4411552 Flavobacterium species	99		4413110 Flavobacterium species	99		4414114 Flavobacterium species	99	
4411554 Flavobacterium species	99		4413112 Flavobacterium species	99		4414116 Flavobacterium species	99	
4411610 Flavobacterium species	99		4413114 Flavobacterium species	99		4414130 Flavobacterium species	99	
4411612 Flavobacterium species	99		4413150 Flavobacterium species	99		4414132 Flavobacterium species	99	
4411710 Flavobacterium species	99		4413310 Flavobacterium species	99		4414150 Flavobacterium species	99	
4411712 Flavobacterium species	99		4413312 Flavobacterium species	99		4414152 Flavobacterium species	99	
4411714 Flavobacterium species	99		4413350 Flavobacterium species	99		4414154 Flavobacterium species	99	
4411716 Flavobacterium species	99		4413510 Flavobacterium species	99		4414207 Bordetella bronchiseptica	99	
4411730 Flavobacterium species	99		4413512 Flavobacterium species	99		4414210 Flavobacterium species	99	
4411732 Flavobacterium species	99		4413514 Flavobacterium species	99		4414212 Flavobacterium species	99	
4411750 Flavobacterium species	99		4413550 Flavobacterium species	99		4414310 Flavobacterium species	99	
4411752 Flavobacterium species	99		4413710 Flavobacterium species	99		4414312 Flavobacterium species	99	
4411754 Flavobacterium species	99		4413712 Flavobacterium species	99		4414314 Flavobacterium species	99	
4411512						4414314		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4415252 Flavobacterium species	99		4415416 Flavobacterium species	99		4415572 Flavobacterium species	99	
4415300 Flavobacterium species	99		4415430 Flavobacterium species	99		4415610 Flavobacterium species	99	
4415310 Flavobacterium species	99		4415450 Flavobacterium species	99		4415612 Flavobacterium species	99	
4415312 Flavobacterium species	99		4415452 Flavobacterium species	99		4415614 Flavobacterium species	99	
4415314 Flavobacterium species	99		4415500 Flavobacterium species	99		4415616 Flavobacterium species	99	
4415316 Flavobacterium species	99		4415502 Flavobacterium species	99		4415630 Flavobacterium species	99	
4415330 Flavobacterium species	99		4415510 Flavobacterium species	99		4415650 Flavobacterium species	99	
4415332 Flavobacterium species	99		4415512 Flavobacterium species	99		4415652 Flavobacterium species	99	
4415334 Flavobacterium species	99		4415514 Flavobacterium species	99		4415700 Flavobacterium species	99	
4415336 Flavobacterium species	99		4415516 Flavobacterium species	99		4415710 Flavobacterium species	99	
4415350 Flavobacterium species	99		4415530 Flavobacterium species	99		4415712 Flavobacterium species	99	
4415352 Flavobacterium species	99		4415532 Flavobacterium species	99		4415714 Flavobacterium species	99	
4415354 Flavobacterium species	99		4415534 Flavobacterium species	99		4415716 Flavobacterium species	99	
4415356 Flavobacterium species	99		4415536 Flavobacterium species	99		4415730 Flavobacterium species	99	
4415370 Flavobacterium species	99		4415550 Flavobacterium species	99		4415732 Flavobacterium species	99	
4415372 Flavobacterium species	99		4415552 Flavobacterium species	99		4415734 Flavobacterium species	99	
4415410 Flavobacterium species	99		4415554 Flavobacterium species	99		4415736 Flavobacterium species	99	
4415412 Flavobacterium species	99		4415556 Flavobacterium species	99		4415750 Flavobacterium species	99	
4415414 Flavobacterium species	99		4415570 Flavobacterium species	99		4415752 Flavobacterium species	99	
4415252						4415752		

BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %	BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %	BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %
4415754 Flavobacterium species	99	4417012 Flavobacterium species	99	4417510 Flavobacterium species	99
4415756 Flavobacterium species	99	4417110 Flavobacterium species	99	4417512 Flavobacterium species	99
4415770 Flavobacterium species	99	4417112 Flavobacterium species	99	4417514 Flavobacterium species	99
4415772 Flavobacterium species	99	4417114 Flavobacterium species	99	4417516 Flavobacterium species	99
4416110 Flavobacterium species	99	4417116 Flavobacterium species	99	4417530 Flavobacterium species	99
4416112 Flavobacterium species	99	4417130 Flavobacterium species	99	4417532 Flavobacterium species	99
4416114 Flavobacterium species	99	4417132 Flavobacterium species	99	4417550 Flavobacterium species	99
4416150 Flavobacterium species	99	4417150 Flavobacterium species	99	4417552 Flavobacterium species	99
4416310 Flavobacterium species	99	4417152 Flavobacterium species	99	4417610 Flavobacterium species	99
4416312 Flavobacterium species	99	4417210 Flavobacterium species	99	4417710 Flavobacterium species	99
4416350 Flavobacterium species	99	4417310 Flavobacterium species	99	4417712 Flavobacterium species	99
4416510 Flavobacterium species	99	4417312 Flavobacterium species	99	4417714 Flavobacterium species	99
4416512 Flavobacterium species	99	4417314 Flavobacterium species	99	4417716 Flavobacterium species	99
4416514 Flavobacterium species	99	4417316 Flavobacterium species	99	4417730 Flavobacterium species	99
4416550 Flavobacterium species	99	4417330 Flavobacterium species	99	4417750 Flavobacterium species	99
4416710 Flavobacterium species	99	4417350 Flavobacterium species	99	4417752 Flavobacterium species	99
4416712 Flavobacterium species	99	4417352 Flavobacterium species	99	4420003 Pseudomonas cepacia	99
4416750 Flavobacterium species	99	4417410 Flavobacterium species	99	4420007 Pseudomonas cepacia	99
4417010 Flavobacterium species	99	4417412 Flavobacterium species	99	4420017 Pseudomonas cepacia	99
4415754				4420017	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4420463 Pseudomonas cepacia	99		4420577 Pseudomonas cepacia	99		4420713 Pseudomonas cepacia	99	
4420467 Pseudomonas cepacia	99		4420603 Pseudomonas cepacia	99		4420717 Pseudomonas cepacia	99	
4420477 Pseudomonas cepacia	99		4420607 Pseudomonas cepacia	99		4420723 Pseudomonas cepacia	99	
4420500 -1 Sphingomonas paucimobilis	99		4420613 Pseudomonas cepacia	99		4420727 Pseudomonas cepacia	99	
4420503 Pseudomonas cepacia	99		4420617 Pseudomonas cepacia	99		4420733 Pseudomonas cepacia	99	
4420507 Pseudomonas cepacia	99		4420623 Pseudomonas cepacia	99		4420737 Pseudomonas cepacia	99	
4420513 Pseudomonas cepacia	99		4420627 Pseudomonas cepacia	99		4420743 Pseudomonas cepacia	99	
4420517 Pseudomonas cepacia	99		4420633 Pseudomonas cepacia	99		4420747 Pseudomonas cepacia	99	
4420523 Pseudomonas cepacia	99		4420637 Pseudomonas cepacia	99		4420753 Pseudomonas cepacia	99	
4420527 Pseudomonas cepacia	99		4420643 Pseudomonas cepacia	99		4420757 Pseudomonas cepacia	99	
4420533 Pseudomonas cepacia	99		4420647 Pseudomonas cepacia	99		4420763 Pseudomonas cepacia	99	
4420537 Pseudomonas cepacia	99		4420653 Pseudomonas cepacia	99		4420767 Pseudomonas cepacia	99	
4420543 Pseudomonas cepacia	99		4420657 Pseudomonas cepacia	99		4420773 Pseudomonas cepacia	99	
4420547 Pseudomonas cepacia	99		4420663 Pseudomonas cepacia	99		4420777 Pseudomonas cepacia	99	
4420553 Pseudomonas cepacia	99		4420667 Pseudomonas cepacia	99		4422307 Pseudomonas cepacia	99	
4420557 Pseudomonas cepacia	99		4420673 Pseudomonas cepacia	99		4422327 Pseudomonas cepacia	99	
4420563 Pseudomonas cepacia	99		4420677 Pseudomonas cepacia	99		4422347 Pseudomonas cepacia	99	
4420567 Pseudomonas cepacia	99		4420703 Pseudomonas cepacia	99		4422507 Pseudomonas cepacia	99	
4420573 Pseudomonas cepacia	99		4420707 Pseudomonas cepacia	99		4422527 Pseudomonas cepacia	99	
4420463						4422527		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4422547 Pseudomonas cepacia	99		4424147 Pseudomonas cepacia	99		4424407 Pseudomonas cepacia	99	
4422603 Pseudomonas cepacia	99		4424167 Pseudomonas cepacia	99		4424427 Pseudomonas cepacia	99	
4422607 Pseudomonas cepacia	99		4424203 Pseudomonas cepacia	99		4424447 Pseudomonas cepacia	99	
4422627 Pseudomonas cepacia	99		4424207 Pseudomonas cepacia	99		4424467 Pseudomonas cepacia	99	
4422647 Pseudomonas cepacia	99		4424227 Pseudomonas cepacia	99		4424503 Pseudomonas cepacia	99	
4422667 Pseudomonas cepacia	99		4424247 Pseudomonas cepacia	99		4424507 Pseudomonas cepacia	99	
4422703 Pseudomonas cepacia	99		4424267 Pseudomonas cepacia	99		4424517 Pseudomonas cepacia	99	
4422707 Pseudomonas cepacia	99		4424303 Pseudomonas cepacia	99		4424523 Pseudomonas cepacia	99	
4422717 Pseudomonas cepacia	99		4424307 Pseudomonas cepacia	99		4424527 Pseudomonas cepacia	99	
4422723 Pseudomonas cepacia	99		4424317 Pseudomonas cepacia	99		4424537 Pseudomonas cepacia	99	
4422727 Pseudomonas cepacia	99		4424323 Pseudomonas cepacia	99		4424543 Pseudomonas cepacia	99	
4422743 Pseudomonas cepacia	99		4424327 Pseudomonas cepacia	99		4424547 Pseudomonas cepacia	99	
4422747 Pseudomonas cepacia	99		4424337 Pseudomonas cepacia	99		4424557 Pseudomonas cepacia	99	
4422763 Pseudomonas cepacia	99		4424343 Pseudomonas cepacia	99		4424563 Pseudomonas cepacia	99	
4422767 Pseudomonas cepacia	99		4424347 Pseudomonas cepacia	99		4424567 Pseudomonas cepacia	99	
4424007 Pseudomonas cepacia	99		4424357 Pseudomonas cepacia	99		4424603 Pseudomonas cepacia	99	
4424103 Pseudomonas cepacia	99		4424363 Pseudomonas cepacia	99		4424607 Pseudomonas cepacia	99	
4424107 Pseudomonas cepacia	99		4424367 Pseudomonas cepacia	99		4424617 Pseudomonas cepacia	99	
4424127 Pseudomonas cepacia	99		4424403 Pseudomonas cepacia	99		4424623 Pseudomonas cepacia	99	

4422547

4424623

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS.	%	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS.	%	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS.	%
4424627 Pseudomonas cepacia	99			4424763 Pseudomonas cepacia	99			4444107 Pseudomonas stutzeri	99		
4424637 Pseudomonas cepacia	99			4424767 Pseudomonas cepacia	99			4444157 Aeromonas hydrophila	99		
4424643 Pseudomonas cepacia	99			4424777 Pseudomonas cepacia	99			4444517 Aeromonas hydrophila	99		
4424647 Pseudomonas cepacia	99			4426707 Pseudomonas cepacia	99			4444557 Aeromonas hydrophila	99		
4424657 Pseudomonas cepacia	99			4440157 Aeromonas hydrophila	99			4445117 Aeromonas hydrophila	99		
4424663 Pseudomonas cepacia	99			4440500 -1 Sphingomonas paucimobilis	99			4445157 Aeromonas hydrophila	99		
4424667 Pseudomonas cepacia	99			4440517 Aeromonas hydrophila	99			4445517 Aeromonas hydrophila	99		
4424703 Pseudomonas cepacia	99			4440557 Aeromonas hydrophila	99			4445557 Aeromonas hydrophila	99		
4424707 Pseudomonas cepacia	99			4441117 Aeromonas hydrophila	99			4446317 Pseudomonas aeruginosa	99		
4424713 Pseudomonas cepacia	99			4441157 Aeromonas hydrophila	99			4446337 Pseudomonas aeruginosa	99		
4424717 Pseudomonas cepacia	99			4441517 Aeromonas hydrophila	99			4446557 Aeromonas hydrophila	99		
4424723 Pseudomonas cepacia	99			4441557 Aeromonas hydrophila	99			4447117 Aeromonas hydrophila	99		
4424727 Pseudomonas cepacia	99			4442557 Aeromonas hydrophila	99			4447157 Aeromonas hydrophila	99		
4424733 Pseudomonas cepacia	99			4443117 Aeromonas hydrophila	99			4447517 Aeromonas hydrophila	99		
4424737 Pseudomonas cepacia	99			4443157 Aeromonas hydrophila	99			4447557 Aeromonas hydrophila	99		
4424743 Pseudomonas cepacia	99			4443517 Aeromonas hydrophila	99			4451557 Aeromonas hydrophila	99		
4424747 Pseudomonas cepacia	99			4443557 Aeromonas hydrophila	99			4455157 Aeromonas hydrophila	99		
4424753 Pseudomonas cepacia	99			4444003 Pseudomonas putida	99			4455557 Aeromonas hydrophila	99		
4424757 Pseudomonas cepacia	99			4444103 Pseudomonas fluorescens	99			4461147 Plesiomonas shigelloides	99		
4424627								4461147			

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4500345 -7 Pseudomonas cepacia	99		4500463 Pseudomonas cepacia	99		4500573 Pseudomonas cepacia	99	
4500353 Pseudomonas cepacia	99		4500467 Pseudomonas cepacia	99		4500577 Pseudomonas cepacia	99	
4500357 Pseudomonas cepacia	99		4500473 Pseudomonas cepacia	99		4500602 -3 Pseudomonas cepacia	99	
4500363 Pseudomonas cepacia	99		4500477 Pseudomonas cepacia	99		4500605 -7 Pseudomonas cepacia	99	
4500366 -7 Pseudomonas cepacia	99		4500500 -1 Sphingomonas paucimobilis	99		4500613 Pseudomonas cepacia	99	
4500373 Pseudomonas cepacia	99		4500502 -3 Pseudomonas cepacia	99		4500616 -7 Pseudomonas cepacia	99	
4500377 Pseudomonas cepacia	99		4500505 -7 Pseudomonas cepacia	99		4500622 -3 Pseudomonas cepacia	99	
4500403 Pseudomonas cepacia	99		4500513 Pseudomonas cepacia	99		4500625 -7 Pseudomonas cepacia	99	
4500406 -7 Pseudomonas cepacia	99		4500516 -7 Pseudomonas cepacia	99		4500633 Pseudomonas cepacia	99	
4500413 Pseudomonas cepacia	99		4500522 -3 Pseudomonas cepacia	99		4500637 Pseudomonas cepacia	99	
4500417 Pseudomonas cepacia	99		4500525 -7 Pseudomonas cepacia	99		4500642 -3 Pseudomonas cepacia	99	
4500423 Pseudomonas cepacia	99		4500533 Pseudomonas cepacia	99		4500645 -7 Pseudomonas cepacia	99	
4500426 -7 Pseudomonas cepacia	99		4500537 Pseudomonas cepacia	99		4500653 Pseudomonas cepacia	99	
4500433 Pseudomonas cepacia	99		4500542 -3 Pseudomonas cepacia	99		4500657 Pseudomonas cepacia	99	
4500437 Pseudomonas cepacia	99		4500545 -7 Pseudomonas cepacia	99		4500662 -3 Pseudomonas cepacia	99	
4500443 Pseudomonas cepacia	99		4500553 Pseudomonas cepacia	99		4500666 -7 Pseudomonas cepacia	99	
4500446 -7 Pseudomonas cepacia	99		4500557 Pseudomonas cepacia	99		4500673 Pseudomonas cepacia	99	
4500453 Pseudomonas cepacia	99		4500563 Pseudomonas cepacia	99		4500677 Pseudomonas cepacia	99	
4500457 Pseudomonas cepacia	99		4500566 -7 Pseudomonas cepacia	99		4500701 -3 Pseudomonas cepacia	99	
4500345						4500703		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4500705 -7 Pseudomonas cepacia	99		4501147 Pseudomonas cepacia	99		4501357 Pseudomonas cepacia	99	
4500712 -3 Pseudomonas cepacia	99		4501167 Pseudomonas cepacia	99		4501363 Pseudomonas cepacia	99	
4500715 -7 Pseudomonas cepacia	99		4501203 Pseudomonas cepacia	99		4501367 Pseudomonas cepacia	99	
4500721 -3 Pseudomonas cepacia	99		4501207 Pseudomonas cepacia	99		4501377 Pseudomonas cepacia	99	
4500725 -7 Pseudomonas cepacia	99		4501217 Pseudomonas cepacia	99		4501403 Pseudomonas cepacia	99	
4500733 Pseudomonas cepacia	99		4501223 Pseudomonas cepacia	99		4501407 Pseudomonas cepacia	99	
4500736 -7 Pseudomonas cepacia	99		4501227 Pseudomonas cepacia	99		4501417 Pseudomonas cepacia	99	
4500741 -3 Pseudomonas cepacia	99		4501243 Pseudomonas cepacia	99		4501423 Pseudomonas cepacia	99	
4500745 -7 Pseudomonas cepacia	99		4501247 Pseudomonas cepacia	99		4501427 Pseudomonas cepacia	99	
4500753 Pseudomonas cepacia	99		4501267 Pseudomonas cepacia	99		4501443 Pseudomonas cepacia	99	
4500756 -7 Pseudomonas cepacia	99		4501303 Pseudomonas cepacia	99		4501447 Pseudomonas cepacia	99	
4500762 -3 Pseudomonas cepacia	99		4501307 Pseudomonas cepacia	99		4501467 Pseudomonas cepacia	99	
4500765 -7 Pseudomonas cepacia	99		4501313 Pseudomonas cepacia	99		4501503 Pseudomonas cepacia	99	
4500773 Pseudomonas cepacia	99		4501317 Pseudomonas cepacia	99		4501507 Pseudomonas cepacia	99	
4500776 -7 Pseudomonas cepacia	99		4501323 Pseudomonas cepacia	99		4501513 Pseudomonas cepacia	99	
4501007 Pseudomonas cepacia	99		4501327 Pseudomonas cepacia	99		4501517 Pseudomonas cepacia	99	
4501103 Pseudomonas cepacia	99		4501337 Pseudomonas cepacia	99		4501523 Pseudomonas cepacia	99	
4501107 Pseudomonas cepacia	99		4501343 Pseudomonas cepacia	99		4501527 Pseudomonas cepacia	99	
4501127 Pseudomonas cepacia	99		4501347 Pseudomonas cepacia	99		4501537 Pseudomonas cepacia	99	
4500705						4501537		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4504027 Pseudomonas cepacia	99		4504143 Pseudomonas cepacia	99		4504257 Pseudomonas cepacia	99	
4504033 Pseudomonas cepacia	99		4504147 Pseudomonas cepacia	99		4504263 Pseudomonas cepacia	99	
4504037 Pseudomonas cepacia	99		4504153 Pseudomonas cepacia	99		4504267 Pseudomonas cepacia	99	
4504043 Pseudomonas cepacia	99		4504157 Pseudomonas cepacia	99		4504273 Pseudomonas cepacia	99	
4504047 Pseudomonas cepacia	99		4504163 Pseudomonas cepacia	99		4504277 Pseudomonas cepacia	99	
4504053 Pseudomonas cepacia	99		4504167 Pseudomonas cepacia	99		4504303 Pseudomonas cepacia	99	
4504057 Pseudomonas cepacia	99		4504173 Pseudomonas cepacia	99		4504306 -7 Pseudomonas cepacia	99	
4504063 Pseudomonas cepacia	99		4504177 Pseudomonas cepacia	99		4504313 Pseudomonas cepacia	99	
4504067 Pseudomonas cepacia	99		4504203 Pseudomonas cepacia	99		4504317 Pseudomonas cepacia	99	
4504073 Pseudomonas cepacia	99		4504207 Pseudomonas cepacia	99		4504323 Pseudomonas cepacia	99	
4504077 Pseudomonas cepacia	99		4504213 Pseudomonas cepacia	99		4504327 Pseudomonas cepacia	99	
4504103 Pseudomonas cepacia	99		4504217 Pseudomonas cepacia	99		4504333 Pseudomonas cepacia	99	
4504107 Pseudomonas cepacia	99		4504223 Pseudomonas cepacia	99		4504337 Pseudomonas cepacia	99	
4504113 Pseudomonas cepacia	99		4504227 Pseudomonas cepacia	99		4504343 Pseudomonas cepacia	99	
4504117 Pseudomonas cepacia	99		4504233 Pseudomonas cepacia	99		4504347 Pseudomonas cepacia	99	
4504123 Pseudomonas cepacia	99		4504237 Pseudomonas cepacia	99		4504353 Pseudomonas cepacia	99	
4504127 Pseudomonas cepacia	99		4504243 Pseudomonas cepacia	99		4504357 Pseudomonas cepacia	99	
4504133 Pseudomonas cepacia	99		4504247 Pseudomonas cepacia	99		4504363 Pseudomonas cepacia	99	
4504137 Pseudomonas cepacia	99		4504253 Pseudomonas cepacia	99		4504367 Pseudomonas cepacia	99	
4504027						4504367		

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4504373 Pseudomonas cepacia	99		4504506 -7 Pseudomonas cepacia	99		4504623 Pseudomonas cepacia	99	
4504377 Pseudomonas cepacia	99		4504513 Pseudomonas cepacia	99		4504627 Pseudomonas cepacia	99	
4504403 Pseudomonas cepacia	99		4504517 Pseudomonas cepacia	99		4504633 Pseudomonas cepacia	99	
4504407 Pseudomonas cepacia	99		4504523 Pseudomonas cepacia	99		4504637 Pseudomonas cepacia	99	
4504413 Pseudomonas cepacia	99		4504527 Pseudomonas cepacia	99		4504643 Pseudomonas cepacia	99	
4504417 Pseudomonas cepacia	99		4504533 Pseudomonas cepacia	99		4504647 Pseudomonas cepacia	99	
4504423 Pseudomonas cepacia	99		4504537 Pseudomonas cepacia	99		4504653 Pseudomonas cepacia	99	
4504427 Pseudomonas cepacia	99		4504543 Pseudomonas cepacia	99		4504657 Pseudomonas cepacia	99	
4504433 Pseudomonas cepacia	99		4504547 Pseudomonas cepacia	99		4504663 Pseudomonas cepacia	99	
4504437 Pseudomonas cepacia	99		4504553 Pseudomonas cepacia	99		4504667 Pseudomonas cepacia	99	
4504443 Pseudomonas cepacia	99		4504557 Pseudomonas cepacia	99		4504673 Pseudomonas cepacia	99	
4504447 Pseudomonas cepacia	99		4504563 Pseudomonas cepacia	99		4504677 Pseudomonas cepacia	99	
4504453 Pseudomonas cepacia	99		4504567 Pseudomonas cepacia	99		4504702 -3 Pseudomonas cepacia	99	
4504457 Pseudomonas cepacia	99		4504573 Pseudomonas cepacia	99		4504705 -7 Pseudomonas cepacia	99	
4504463 Pseudomonas cepacia	99		4504577 Pseudomonas cepacia	99		4504713 Pseudomonas cepacia	99	
4504467 Pseudomonas cepacia	99		4504603 Pseudomonas cepacia	99		4504717 Pseudomonas cepacia	99	
4504473 Pseudomonas cepacia	99		4504606 -7 Pseudomonas cepacia	99		4504723 Pseudomonas cepacia	99	
4504477 Pseudomonas cepacia	99		4504613 Pseudomonas cepacia	99		4504726 -7 Pseudomonas cepacia	99	
4504503 Pseudomonas cepacia	99		4504617 Pseudomonas cepacia	99		4504733 Pseudomonas cepacia	99	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
4522127 Pseudomonas cepacia	99		4522447 Pseudomonas cepacia	99		4522657 Pseudomonas cepacia	99	
4522147 Pseudomonas cepacia	99		4522467 Pseudomonas cepacia	99		4522663 Pseudomonas cepacia	99	
4522203 Pseudomonas cepacia	99		4522503 Pseudomonas cepacia	99		4522667 Pseudomonas cepacia	99	
4522207 Pseudomonas cepacia	99		4522507 Pseudomonas cepacia	99		4522703 Pseudomonas cepacia	99	
4522227 Pseudomonas cepacia	99		4522517 Pseudomonas cepacia	99		4522707 Pseudomonas cepacia	99	
4522247 Pseudomonas cepacia	99		4522523 Pseudomonas cepacia	99		4522713 Pseudomonas cepacia	99	
4522267 Pseudomonas cepacia	99		4522527 Pseudomonas cepacia	99		4522717 Pseudomonas cepacia	99	
4522303 Pseudomonas cepacia	99		4522543 Pseudomonas cepacia	99		4522723 Pseudomonas cepacia	99	
4522307 Pseudomonas cepacia	99		4522547 Pseudomonas cepacia	99		4522727 Pseudomonas cepacia	99	
4522317 Pseudomonas cepacia	99		4522563 Pseudomonas cepacia	99		4522733 Pseudomonas cepacia	99	
4522323 Pseudomonas cepacia	99		4522567 Pseudomonas cepacia	99		4522737 Pseudomonas cepacia	99	
4522327 Pseudomonas cepacia	99		4522603 Pseudomonas cepacia	99		4522743 Pseudomonas cepacia	99	
4522343 Pseudomonas cepacia	99		4522607 Pseudomonas cepacia	99		4522747 Pseudomonas cepacia	99	
4522347 Pseudomonas cepacia	99		4522617 Pseudomonas cepacia	99		4522753 Pseudomonas cepacia	99	
4522363 Pseudomonas cepacia	99		4522623 Pseudomonas cepacia	99		4522757 Pseudomonas cepacia	99	
4522367 Pseudomonas cepacia	99		4522627 Pseudomonas cepacia	99		4522763 Pseudomonas cepacia	99	
4522403 Pseudomonas cepacia	99		4522637 Pseudomonas cepacia	99		4522767 Pseudomonas cepacia	99	
4522407 Pseudomonas cepacia	99		4522643 Pseudomonas cepacia	99		4522777 Pseudomonas cepacia	99	
4522427 Pseudomonas cepacia	99		4522647 Pseudomonas cepacia	99		4524003 Pseudomonas cepacia	99	
4522127						4524003		

BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %	BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %	BIOTYPE NUMBER ORGANISM NAME	RP ADD.TEST % POS. %
4524527 Pseudomonas cepacia	99	4524647 Pseudomonas cepacia	99	4524763 Pseudomonas cepacia	99
4524533 Pseudomonas cepacia	99	4524653 Pseudomonas cepacia	99	4524767 Pseudomonas cepacia	99
4524537 Pseudomonas cepacia	99	4524657 Pseudomonas cepacia	99	4524773 Pseudomonas cepacia	99
4524543 Pseudomonas cepacia	99	4524663 Pseudomonas cepacia	99	4524777 Pseudomonas cepacia	99
4524547 Pseudomonas cepacia	99	4524667 Pseudomonas cepacia	99	4525511 Vibrio cholerae	99
4524553 Pseudomonas cepacia	99	4524673 Pseudomonas cepacia	99	4525711 Vibrio cholerae	99
4524557 Pseudomonas cepacia	99	4524677 Pseudomonas cepacia	99	4526307 Pseudomonas cepacia	99
4524563 Pseudomonas cepacia	99	4524703 Pseudomonas cepacia	99	4526507 Pseudomonas cepacia	99
4524567 Pseudomonas cepacia	99	4524707 Pseudomonas cepacia	99	4526607 Pseudomonas cepacia	99
4524577 Pseudomonas cepacia	99	4524713 Pseudomonas cepacia	99	4526627 Pseudomonas cepacia	99
4524603 Pseudomonas cepacia	99	4524717 Pseudomonas cepacia	99	4526647 Pseudomonas cepacia	99
4524607 Pseudomonas cepacia	99	4524723 Pseudomonas cepacia	99	4526703 Pseudomonas cepacia	99
4524613 Pseudomonas cepacia	99	4524727 Pseudomonas cepacia	99	4526707 Pseudomonas cepacia	99
4524617 Pseudomonas cepacia	99	4524733 Pseudomonas cepacia	99	4526717 Pseudomonas cepacia	99
4524623 Pseudomonas cepacia	99	4524737 Pseudomonas cepacia	99	4526723 Pseudomonas cepacia	99
4524627 Pseudomonas cepacia	99	4524743 Pseudomonas cepacia	99	4526727 Pseudomonas cepacia	99
4524633 Pseudomonas cepacia	99	4524747 Pseudomonas cepacia	99	4526743 Pseudomonas cepacia	99
4524637 Pseudomonas cepacia	99	4524753 Pseudomonas cepacia	99	4526747 Pseudomonas cepacia	99
4524643 Pseudomonas cepacia	99	4524757 Pseudomonas cepacia	99	4526767 Pseudomonas cepacia	99
4524527				4526767	

BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %	BIOTYPE NUMBER ORGANISM NAME	RP %	ADD.TEST POS. %
5000107 Pseudomonas cepacia	99		5000213 Pseudomonas cepacia	99		5000327 Pseudomonas cepacia	99	
5000113 Pseudomonas cepacia	99		5000217 Pseudomonas cepacia	99		5000333 Pseudomonas cepacia	99	
5000117 Pseudomonas cepacia	99		5000223 Pseudomonas cepacia	99		5000337 Pseudomonas cepacia	99	
5000123 Pseudomonas cepacia	99		5000227 Pseudomonas cepacia	99		5000343 Pseudomonas cepacia	99	
5000127 Pseudomonas cepacia	99		5000233 Pseudomonas cepacia	99		5000347 Pseudomonas cepacia	99	
5000133 Pseudomonas cepacia	99		5000237 Pseudomonas cepacia	99		5000353 Pseudomonas cepacia	99	
5000137 Pseudomonas cepacia	99		5000243 Pseudomonas cepacia	99		5000357 Pseudomonas cepacia	99	
5000143 Pseudomonas cepacia	99		5000247 Pseudomonas cepacia	99		5000363 Pseudomonas cepacia	99	
5000147 Pseudomonas cepacia	99		5000253 Pseudomonas cepacia	99		5000367 Pseudomonas cepacia	99	
5000153 Pseudomonas cepacia	99		5000257 Pseudomonas cepacia	99		5000373 Pseudomonas cepacia	99	
5000157 Pseudomonas cepacia	99		5000263 Pseudomonas cepacia	99		5000377 Pseudomonas cepacia	99	
5000163 Pseudomonas cepacia	99		5000267 Pseudomonas cepacia	99		5000403 Pseudomonas cepacia	99	
5000167 Pseudomonas cepacia	99		5000273 Pseudomonas cepacia	99		5000407 Pseudomonas cepacia	99	
5000173 Pseudomonas cepacia	99		5000277 Pseudomonas cepacia	99		5000413 Pseudomonas cepacia	99	
5000177 Pseudomonas cepacia	99		5000303 Pseudomonas cepacia	99		5000417 Pseudomonas cepacia	99	
5000200 Pseudomonas species Moraxella species	63 36		5000307 Achromobacter xylosoxid. Pseudomonas cepacia	72 27	GRSS/100 GRSS/ 9	5000423 Pseudomonas cepacia	99	
5000201 -6 Pseudomonas species	99		5000313 Pseudomonas cepacia	99		5000427 Pseudomonas cepacia	99	
5000207 Pseudomonas species Achromobacter xylosoxid. Bordetella bronchiseptica	81 15 2		5000317 Pseudomonas cepacia	99		5000433 Pseudomonas cepacia	99	
			5000323 Pseudomonas cepacia	99		5000437 Pseudomonas cepacia	99	
5000107						5000437		