



Certificate

Product:	IH-Card Group ABO	REF	12012727
Intended Use:	<p>IH-Card Group ABO is intended to be used for the confirmation of the A and B antigen typing, on the human red blood cells, conducted with anti-A (ABO1) and anti-B (ABO2) Blood Grouping Reagents, on a single Gel Card.</p> <p>The IH-Card Group ABO can be used for performing three forward ABO groupings on a single IH-Card confirming the ABO blood group of donors and patients.</p> <p>Each card contains six microtubes.</p>	LOT	9315010
			2024-08-11
			18...25°C

Reagents:

Reagent	Clone	Source	Raw Material Lot.-No.:
Anti-A	15750F7	Murine Monoclonal	64444302
Anti-B	X9	Murine Monoclonal	64451850, 64400920

- **Information, Instructions and Safety:**

- For further information, instructions and safety refer to the included instructions for use of this product.

- **Analytic Results, Quality Control:**

- The product lot was released following evaluation of sensitivity and specificity.

- **Regulation:**

- DIN EN ISO 13485
- MDSAP
- FDA 21 CFR Part 820, 809, 600

This document has been produced electronically and is valid without a signature, according to SOP SV-D:A-042-00