




Certificate

Product:	Anti-C3d	REF	12013874
Intended Use:	Anti-C3d is used for qualitative test of complement bounding to human red blood cells. It is used for the direct antiglobulin-test to determine the presence or absence of complement (C3b and/or C3d) on red blood cells in the tube test in patient or donor samples.	LOT	3333010-00
			2023-08-14
			2025-02-04
		ACT	1:16
		VOL	5 ml
			2°C / 8°C

- **Specifications:**

Clone: /
Anti-C3d (rabbit, polyclonal)

- **Information, Instructions and Safety:**

- For further information, instructions and safety refer to the Instructions for Use and the Material Safety Data Sheet of this product.

- **Analytic Results, Quality Control:**

- The product lot was released following evaluation of sensitivity and specificity.

- **Regulation:**

- The product has been manufactured and tested by Bio-Rad Medical Diagnostics GmbH in accordance with the currently applicable regulations for the manufacture of in-vitro diagnostics.

This document has been produced electronically and is valid without a signature, according to SOP SV-D:A-042-00