

Certificate of Serafol® Products

CERTIFIED QUALITY MANAGEMENT SYSTEM IN
ACCORDANCE WITH DIN EN ISO 9001 AND DIN EN ISO 13485

sifin diagnostics gmbh

Berliner Allee 317-321
13088 Berlin, Germany


T: +49 30 700 144-0
F: +49 30 700 144-30
info@sifin.de
www.sifin.de

Product	Serafol® AB0+D	
Lot No.	6950922	
Storage	2...25 °C	
Expiry date	September 30, 2024	
Materials	Reagents	Clone
	Anti-A monoclonal	A003
	Anti-B monoclonal	B005
	Anti-D monoclonal	BS226
	Serafol® products are cards with dried specific monoclonal reagents for bedside testing.	
Preservative	<0.1 % NaN ₃	
Manufacturer, Distributor	sifin diagnostics gmbh, Berlin, Germany	
Distributor	Bio-Rad Medical Diagnostics GmbH, Dreieich, Germany	

sifin diagnostics gmbh is regularly controlled by Notified Body 0483 and Paul-Ehrlich-Institut.

For practical purposes, the risk of contamination with infectious agents can be excluded by virtue of the use of biotechnology production methods of all monoclonal reagents.

Appropriate safety precautions are recommended.



Johner Medical Schweiz GmbH
Tafelstattstrasse 13a, 6415 Arth, Switzerland



Bio-Rad Laboratories AG
Pra Rond 23, 1785 Cressier FR, Switzerland

This document (Rev.1) has been established electronically and is valid without signature.
The batch release was done in accordance with SOP QMV 7.