



## Certificate / Zertifikat

<b>Product:</b>	Seraclone® Anti- $\bar{s}$ (MNS4)	<b>REF</b>	808068
<b>Intended Use:</b>	Seraclone® Anti- $\bar{s}$ (MNS4) is used for the determination of the red blood cell antigen $\bar{s}$ . Seraclone® Anti- $\bar{s}$ (MNS4) wird zur Bestimmung des erythrozytären Antigens $\bar{s}$ eingesetzt.	<b>LOT</b>	3221031-00
			2024-05-08
		<b>ACT</b>	1:16
			2 ... 8 °C

- **Specifications**

**Clone:** P3YAN3 (IgG), human monoclonal/ Human monoklonal

- **Information, Instructions and Safety:**

- For further information, instructions and safety refer to the included instructions for use of this product.

- **Analytic Results, Quality Control:**

- The product lot was released following evaluation of sensitivity and specificity.

- **Regulation:**

- DIN EN ISO 13485
- MDSAP
- Directive 98/79/EC
- FDA 21 CFR Part 820, 809, 600, relevant only for products for the US-market
- Recommended Methods for Blood Grouping Reagents Evaluation, Docket No.84S-0181, March 1992

***This document has been produced electronically and is valid without a signature, according to SOP SV-D:A-042-00 /***