

DRB SSP Kit

Primer positions/ Primerpositionen

REF

826215

LOT

B115151-db2020-10V1

187837/03-08/2012

| PRIMERMIX | 5' PRIMER | 3' PRIMER |
|-----------|-----------------------------------|---|
| 1 (H1) | Negative | control |
| 2 (G1) | GAA Exon 2, 42 | GAG Exon 2, 199 |
| 3 (F1) | GAA Exon 2, 42 | GAT Exon 2, 199 |
| 4 (E1) | AGG Exon 2, 39 | CGC Exon 2, 211 GAT Exon 2, 199 |
| 5 (D1) | ACT Exon 2, 91 | CAA Intron 2 |
| 6 (C1) | CGT Exon 2, 37 | GCT Exon 2, 212 |
| 7 (B1) | GCG Intron 1 | CTC Exon 2, 84 |
| 8 (A1) | ATC Exon 2, 97 | TCT Exon 2, 260 |
| 9 (H2) | GGG Exon 2, 32 | CAC Exon 2, 257 |
| 10 (G2) | GTT Exon 2, 46 | CAG Exon 2, 220 CAC Exon 2, 257 CCT Exon 2, 173 |
| 11 (F2) | TAT Exon 2, 78 | CAC Exon 2, 257 |
| 12 (E2) | AGG Exon 2, 28 | CAC Exon 2, 257 |
| 13 (D2) | CCG Intron 1 | CCT Exon 2, 173 |
| 14 (C2) | TCC Exon 2, 112 | CCA Exon 2, 257 |
| 15 (B2) | CGT Exon 2, 37 | GAT Exon 2, 199 |
| 16 (A2) | GTT Exon 2, 46 CGT Exon 2, 37 | CTC Exon 2, 211 GAA Exon 2, 199 |
| 17 (H3) | CGT Exon 2, 37 | CGG Exon 2, 173 CAG Exon 2, 173 |
| 18 (G3) | GTT Exon 2, 46 GTC Exon 2, 38 | TCC Exon 2, 209 |
| 19 (F3) | CAG Exon 2, 102 | TCT Exon 2, 209 CTT Exon 2, 211 |
| 20 (E3) | CGT Exon 2, 37 GTT Exon 2, 110 | CAG Exon 2, 220 GAT Exon 2, 199 |
| 21 (D3) | CGT Exon 2, 37 ACA Exon 2, 38 | CAG Exon 2, 170 |
| 22 (C3) | GCT Exon 2, 29 | AAT Exon 2, 230 AGT Exon 2, 230 |
| 23 (B3) | ACC Exon 2, 130 | CCA Exon 2, 257 |
| 24 (A3) | GTA Exon 2, 38 | TCT Exon 2, 260 |