

DRB SSP Kit

Primer positions/ Primerpositionen

REF

826215

LOT

B203151-db2020-10V1

187837/03-08/2012

PRIMERMIX	5' PRIMER	3' PRIMER
1 (H1)	Negative	control
2 (G1)	GAA Exon 2, 42	GAG Exon 2, 199
3 (F1)	GAA Exon 2, 42	GAT Exon 2, 199
4 (E1)	AGG Exon 2, 39	CGC Exon 2, 211 GAT Exon 2, 199
5 (D1)	ACT Exon 2, 91	CAA Intron 2
6 (C1)	CGT Exon 2, 37	GCT Exon 2, 212
7 (B1)	GCG Intron 1	CTC Exon 2, 84
8 (A1)	ATC Exon 2, 97	TCT Exon 2, 260
9 (H2)	GGG Exon 2, 32	CAC Exon 2, 257
10 (G2)	GTT Exon 2, 46	CAG Exon 2, 220 CAC Exon 2, 257 CCT Exon 2, 173
11 (F2)	TAT Exon 2, 78	CAC Exon 2, 257
12 (E2)	AGG Exon 2, 28	CAC Exon 2, 257
13 (D2)	CCG Intron 1	CCT Exon 2, 173
14 (C2)	TCC Exon 2, 112	CCA Exon 2, 257
15 (B2)	CGT Exon 2, 37	GAT Exon 2, 199
16 (A2)	GTT Exon 2, 46 CGT Exon 2, 37	CTC Exon 2, 211 GAA Exon 2, 199
17 (H3)	CGT Exon 2, 37	CGG Exon 2, 173 CAG Exon 2, 173
18 (G3)	GTT Exon 2, 46 GTC Exon 2, 38	TCC Exon 2, 209
19 (F3)	CAG Exon 2, 102	TCT Exon 2, 209 CTT Exon 2, 211
20 (E3)	CGT Exon 2, 37 GTT Exon 2, 110	CAG Exon 2, 220 GAT Exon 2, 199
21 (D3)	CGT Exon 2, 37 ACA Exon 2, 38	CAG Exon 2, 170
22 (C3)	GCT Exon 2, 29	AAT Exon 2, 230 AGT Exon 2, 230
23 (B3)	ACC Exon 2, 130	CCA Exon 2, 257
24 (A3)	GTA Exon 2, 38	TCT Exon 2, 260