

DQB SSP Kit

Primer positions/ Primerpositionen

REF

826220

LOT

B301201-db2020-10V1

187852/03-08/2012

PRIMERMIX	5' PRIMER	3' PRIMER
1 (H1)	GGG Exon 2, 77	GGT Exon 2, 260
2 (G1)	TTA Exon 2, 77 TCT Exon 2, 77	ACG Exon 2, 257 ACC Exon 2, 257
3 (F1)	AAG Exon 2, 89	GCT Exon 2, 257
4 (E1)	TCT Exon 2, 77 ACC Exon 2, 63	GGC Exon 2, 169
5 (D1)	TTA Exon 2, 77	GCG Exon 2, 164
6 (C1)	TCT Exon 2, 77 ACC Exon 2, 63	GCG Exon 2, 164
7 (B1)	TGA Exon 2, 170	GCT Exon 2, 257
8 (A1)	ACC Exon 2, 63	CTA Exon 2, 186