

HLA-A SSP Kit

Primer positions/ Primerpositionen

REF

826201

LOT

B928012-db2020-10V1

187890/03-08/2012

Primermix	5' PRIMER	3' PRIMER
1 (H1)	CCG Exon 3, 70	GGA Exon 3, 178
2 (G1)	GCA Exon 2, 71 TCA Exon 2, 144	GGA Exon 2, 167
3 (F1)	GAC Intron 1	CCA Exon 2, 226
4 (E1)	ACC Exon 2, 29	GCC Intron 2
5 (D1)	TTA Intron 2	CAC Exon 3, 181
6 (C1)	CAT Exon 3, 182	CCG Intron 3
7 (B1)	GGA Exon 2, 187	CCG Exon 3, 216
8 (A1)	GCC Exon 2, 71	GGC Exon 2, 229
9 (H2)	CCT Intron 1	TCG Exon 2, 226 GTG Exon 2, 219
10 (G2)	GCT Exon 3, 80	CGT Exon 3, 216
11 (F2)	GCT Exon 3, 80	CTC Exon 3, 216
12 (E2)	ACC Exon 2, 29	GAC Exon 2, 209
13 (D2)	GGT Exon 3, 178	CAC Exon 3, 228
14 (C2)	ACC Exon 2, 29	GCA Exon 2, 184
15 (B2)	CCA Exon 3, 12	TGT Exon 3, 163
16 (A2)	CCG Exon 2, 184	GAC Exon 3, 12
17 (H3)	ATG Exon 3, 20	CAA Exon 3, 105
18 (G3)	CAC Exon 2, 25	GCA Exon 2, 184
19 (F3)	ATA Exon 3, 20	TTC Exon 3, 69
20 (E3)	TCA Exon 2, 24	CCT Exon 2, 165 CTC Exon 2, 183
21 (D3)	GCT Exon 2, 244	CGT Intron 2
22 (C3)	ACA Exon 2, 29	GTT Exon 2, 186
23 (B3)	TTG Exon 2, 108	CCA Exon 2, 226
24 (A3)	ATA Exon 3, 20	CTC Exon 3, 216