

HLA-B SSP Kit

Primer positions/ Primerpositionen

REF

826206

LOT

B016062-db2020-10V1

187947/03-08/2012

| PRIMERMIX | 5' PRIMER | 3' PRIMER |
|-----------|---|------------------------------------|
| 1 (1H) | TGA Exon 3, 184 | AAC Intron 3 |
| 2 (1G) | TCG Exon 2, 24 | TGG Exon 2, 199 TGA Exon 2, 199 |
| 3 (1F) | ACC Exon 3, 36 | CTA Intron 3 |
| 4 (1E) | GTT Intron 1 | TGC Exon 2, 199 |
| 5 (1D) | CCG Exon 2, 30 | TAT Exon 2, 173 |
| 6 (1C) | GAG Exon 2, 188 | GTT Intron 2 |
| 7 (1B) | CTT Exon 2, 219 | GTT Intron 2 |
| 8 (1A) | GGA Exon 1, 45 | TCC Exon 2, 193 |
| 9 (2H) | TCT Exon 2, 69 GAA Exon 2, 133 | GTT Intron 2 |
| 10 (2G) | GGT Exon 3, 66 CTA Exon 2, 199 | ATC Intron 3 GTT Intron 2 |
| 11 (2F) | CTT Exon 2, 199 | GTT Intron 2 |
| 12 (2E) | ACG Exon 3, 228 GTC Exon 3, 229 CAC Exon 3, 217 | CAC Intron 3 |
| 13 (2D) | CTT Intron 1 | TCT Exon 2, 161 |
| 14 (2C) | TCC Exon 2, 24 | TCT Exon 2, 131 |
| 15 (2B) | TGA Intron 1 | GGC Exon 2, 229 |
| 16 (2A) | TGA Intron 1 | TGA Exon 2, 199 |
| 17 (3H) | GAC Exon 2, 133 | GGC Exon 2, 229 |
| 18 (3G) | TCT Exon 2, 69 | GTC Exon 2, 228 |
| 19 (3F) | GAA Exon 2, 173 | GCG Exon 2, 246 |
| 20 (3E) | GAA Exon 3, 201 | AGA Intron 3 |
| 21 (3D) | GTT Exon 3, 169 | GCT Exon 3, 262 |
| 22 (3C) | GCT Exon 1, 72 | CGT Exon 2, 120 |
| 23 (3B) | ATA Exon 3, 69 | GAT Intron 3 |
| 24 (3A) | TCA Exon 3, 12 | GCG Exon 3, 229 |
| 25 (4H) | TGG Exon 3, 14 | GCG Exon 3, 229 |
| 26 (4G) | CCA Exon 2, 33 | TAC Exon 2, 197 |
| 27 (4F) | GTT Exon 3, 25 | CTC Exon 3, 216 |
| 28 (4E) | GAA Exon 2, 133 | CGA Exon 2, 238 |
| 29 (4D) | TGG Exon 3, 14 | CCT Exon 3, 184 |
| 30 (4C) | GAA Exon 2, 133 | GGC Exon 2, 229 |
| 31 (4B) | TGG Exon 3, 14 | TCT Exon 3, 92 |
| 32 (4A) | ACC Exon 2, 92 | CGA Exon 2, 238 |
| 33 (5H) | GAC Exon 2, 206 | CTC Exon 2, 259 |
| 34 (5G) | GGC Exon 2, 126 | CTC Exon 2, 259 |
| 35 (5F) | TGA Intron 1 | CGA Exon 2, 238 GGG Exon 2, 238 |
| 36 (5E) | GCA Intron 1 | GCC Exon 2, 209 |
| 37 (5D) | GCA Intron 1 | CCG Exon 2, 102 |
| 38 (5C) | AGA Exon 2, 134 | CTG Exon 2, 207 |
| 39 (5B) | TGG Exon 3, 14 | CGT Exon 3, 216 |
| 40 (5A) | GCC Exon 3, 216 | AGA Intron 3 |
| 41 (6H) | GCA Intron 1 | CCC Exon 2, 183 |
| 42 (6G) | TGA Intron 1 | TCC Exon 2, 193 |
| 43 (6F) | TCT Intron 1 | GGT Exon 2, 204 |
| 44 (6E) | GAC Exon 3, 19 | GAT Intron 3 |
| 45 (6D) | CCT Exon 2, 241 | TCT Exon 3, 92 |
| 46 (6C) | GTT Exon 3, 25 | GCG Exon 3, 229 |
| 47 (6B) | AGG Intron 1 | GCG Exon 2, 246 GGG Exon 2, 246 |
| 48 (6A) | GCA Intron 1 | GCG Exon 2, 244 |