

# HLA-B SSP Kit

Primer positions/ Primerpositionen

REF

826206

LOT

B114061-db2020-10V1

187947/03-08/2012

PRIMERMIX	5' PRIMER	3' PRIMER
1 (1H)	TGA Exon 3, 184	AAC Intron 3
2 (1G)	TCG Exon 2, 24	TGG Exon 2, 199 TGA Exon 2, 199
3 (1F)	ACC Exon 3, 36	CTA Intron 3
4 (1E)	GTT Intron 1	TGC Exon 2, 199
5 (1D)	CCG Exon 2, 30	TAT Exon 2, 173
6 (1C)	GAG Exon 2, 188	GTT Intron 2
7 (1B)	CTT Exon 2, 219	GTT Intron 2
8 (1A)	GGA Exon 1, 45	TCC Exon 2, 193
9 (2H)	TCT Exon 2, 69 GAA Exon 2, 133	GTT Intron 2
10 (2G)	GGT Exon 3, 66 CTA Exon 2, 199	ATC Intron 3 GTT Intron 2
11 (2F)	CTT Exon 2, 199	GTT Intron 2
12 (2E)	ACG Exon 3, 228 GTC Exon 3, 229 CAC Exon 3, 217	CAC Intron 3
13 (2D)	CTT Intron 1	TCT Exon 2, 161
14 (2C)	TCC Exon 2, 24	TCT Exon 2, 131
15 (2B)	TGA Intron 1	GGC Exon 2, 229
16 (2A)	TGA Intron 1	TGA Exon 2, 199
17 (3H)	GAC Exon 2, 133	GGC Exon 2, 229
18 (3G)	TCT Exon 2, 69	GTC Exon 2, 228
19 (3F)	GAA Exon 2, 173	GCG Exon 2, 246
20 (3E)	GAA Exon 3, 201	AGA Intron 3
21 (3D)	GTT Exon 3, 169	GCT Exon 3, 262
22 (3C)	GCT Exon 1, 72	CGT Exon 2, 120
23 (3B)	ATA Exon 3, 69	GAT Intron 3
24 (3A)	TCA Exon 3, 12	GCG Exon 3, 229
25 (4H)	TGG Exon 3, 14	GCG Exon 3, 229
26 (4G)	CCA Exon 2, 33	TAC Exon 2, 197
27 (4F)	GTT Exon 3, 25	CTC Exon 3, 216
28 (4E)	GAA Exon 2, 133	CGA Exon 2, 238
29 (4D)	TGG Exon 3, 14	CCT Exon 3, 184
30 (4C)	GAA Exon 2, 133	GGC Exon 2, 229
31 (4B)	TGG Exon 3, 14	TCT Exon 3, 92
32 (4A)	ACC Exon 2, 92	CGA Exon 2, 238
33 (5H)	GAC Exon 2, 206	CTC Exon 2, 259
34 (5G)	GGC Exon 2, 126	CTC Exon 2, 259
35 (5F)	TGA Intron 1	CGA Exon 2, 238 GGG Exon 2, 238
36 (5E)	GCA Intron 1	GCC Exon 2, 209
37 (5D)	GCA Intron 1	CCG Exon 2, 102
38 (5C)	AGA Exon 2, 134	CTG Exon 2, 207
39 (5B)	TGG Exon 3, 14	CGT Exon 3, 216
40 (5A)	GCC Exon 3, 216	AGA Intron 3
41 (6H)	GCA Intron 1	CCC Exon 2, 183
42 (6G)	TGA Intron 1	TCC Exon 2, 193
43 (6F)	TCT Intron 1	GGT Exon 2, 204
44 (6E)	GAC Exon 3, 19	GAT Intron 3
45 (6D)	CCT Exon 2, 241	TCT Exon 3, 92
46 (6C)	GTT Exon 3, 25	GCG Exon 3, 229
47 (6B)	AGG Intron 1	GCG Exon 2, 246 GGG Exon 2, 246
48 (6A)	GCA Intron 1	GCG Exon 2, 244