

# HLA-C SSP Kit

Primer positions/ Primerpositionen

REF

826212

LOT

B114124-db2020-10V1

187949/03-08/2012

PRIMERMIX	5' PRIMER	3' PRIMER
1 (H1)	Neg. Kontr./contr.	
2 (G1)	GGT Intron 1	GGC Exon 2, 229
3 (F1)	CCA Exon 2, 45	AGT Exon 2, 216
4 (E1)	CAT Exon 3, 10 AGC Exon 3, 20	CTT Exon 3, 246
5 (D1)	GCT Exon 2, 32	ACT Exon 2, 270
6 (C1)	CCG Exon 4, 33	TGA Exon 4, 236
7 (B1)	CGA Intron 1	GGT Exon 2, 229
8 (A1)	AGC Exon 1, 47	CCA Intron 1
9 (H2)	CCC, Exon 2, 140	AGC Exon 2, 216
10 (G2)	CCG Intron 3	GGC Exon 4, 155
11 (F2)	GCA Exon 2, 103	GGC Exon 2, 229
12 (E2)	AGT Exon 3, 18 TAC Exon 3, 26	GCA Exon 3, 131 GCA Exon 3, 134
13 (D2)	AGT Exon 3, 18 TAC Exon 3, 26	TTT Intron 3
14 (C2)	TTA Intron 1	TGG, Exon 2, 25
15 (B2)	TCA Exon 3, 12	ATA Intron 3
16 (A2)	AGT Exon 3, 18	CCG Exon 3, 184
17 (H3)	CCA Exon 3, 12	TCC Exon 3, 237
18 (G3)	AAG Intron 1	GGT Exon 2, 229