

Erytype S Rh-Pheno double

REF 806150

187755 / 30

Ergebnisprotokoll / Worksheet / Feuille de travail /

Protocollo risultati / Hoja de control

Mikrotestplatte (Streifen) zur Bestimmung der Rh-Formel von 12 Probanden in Doppelbestimmung

Microtestplate (strips) for Rh-subtyping of 12 specimens in duplicate

Microplaque (barrettes) pour la détermination des sous-groupes Rh de 12 patients en double

Micrometodo (strips) su piastra per la determinazione dei gruppi sanguigni Rh di 12 campioni in doppio

Microplaca (tiras) para la tipificación del grupo Rh en 12 probandos para determinación doble

Ch.-B. / Lot-No. / Lot Nr. / Lotto / Lote:	
Datum / Date / Data / Fecha:	
Platte Nr. / Plate No. / Plaque No. / Piastra No. / Placa No.:	
Unterschrift / Signature / Firma:	

		1	2	3	4	5	6	7	8	9	10	11	12
A	Reagent I	Anti-C (Clones: MS24/P3x25513 G8)											
		Anti-c (Clone: MS33)											
		Anti-E (Clone: 906)											
		Anti-e (Clones: MS16/MS21/MS63)											
E	Reagent II	Anti-C (Clone: MS273)											
		Anti-c (Clone: MS35)											
		Anti-E (Clones: MS260/MS12)											
		Anti-e (Clones: MS62/MS69)											
Proband/Specimen/Echantillon/Campione/Probando													
Ergebnis/Result/Resultat/Risultato/Resultado													

