

Erytype S CDE+K

REF 806155

187756 / 25

Ergebnisprotokoll / Worksheet / Feuille de travail /

Protocollo risultati / Hoja de control

Mikrotestplatte (Streifen) zur CDE- und Kell-Bestimmung von 24 Probanden (mit Kell in Doppelbestimmung)

Microtestplate (strips) for CDE and Kell typing of 24 specimens (with Kell in duplicate)

Microplaque (bandes) pour la détermination CDE et KEL1 de 24 patients (avec KEL1 en double)

Micrometodo (strips) su piastra per la determinazione dei gruppi sanguigni CDE e Kell di 24 campioni (con Kell in doppio)

Microplaca (tiras) para la tipificación del grupo CDE y Kell en 24 probandos (con Kell determinación doble)

Ch.-B. / Lot-No. / Lot Nr. / Lotto / Lote:	
Datum / Date / Data / Fecha:	
Platte Nr. / Plate No. / Plaque No. / Piastra No. / Placa No.:	
Unterschrift / Signature / Firma:	

		1	2	3	4	5	6	7	8	9	10	11	12
A	Anti-CDE (Clones: P3x25513 G8/P3x61/P3x234)												
B	Anti-K (Clone: MS56)												
C	Anti-K (Clone: 601)												
D	Neg.-Kontr./Neg.-Contr./Contr. Neg.												
Proband/Specimen/Echantillon/Campione/Probando													
Ergebnis/Result/Resultat/Risultato/Resultado													
E	Anti-CDE (Clones: P3x25513 G8/P3x61/P3x234)												
F	Anti-K (Clone: MS56)												
G	Anti-K (Clone: 601)												
H	Neg.-Kontr./Neg.-Contr./Contr. Neg.												
Proband/Specimen/Echantillon/Campione/Probando													
Ergebnis/Result/Resultat/Risultato/Resultado													

