

# Erytype S AB0

REF 806115

187783 / 25

Ergebnisprotokoll / Worksheet / Feuille de travail /

Protocollo risultati / Hoja de control

Mikrotestplatte (Streifen) zur AB0-Bestimmung von 12 Probanden

Microtestplate (strips) for AB0 typing of 12 specimens

Microplaque (barrettes) pour la détermination AB0 de 12 échantillons

Micropiastre (strips) per la determinazione dei gruppi sanguigni AB0 di 12 campioni

Microplaca (tiras) para la tipificación del grupo AB0 en 12 probandos

Ch.-B. / Lot-No. / Lot Nr. / Lotto / Lote:	
Datum / Date / Data / Fecha:	
Platte Nr. / Plate No. / Plaque No. / Piastra No. / Placa No.:	
Unterschrift / Signature / Firma:	

		1	2	3	4	5	6	7	8	9	10	11	12
A	Anti-A (Clone: A003)												
B	Anti-B (Clone: B005)												
C	Anti-AB (Clone: BS63/BS85)												
D	Neg.-Kontr./Neg.-Contr./Contr. Neg.												
E	Serumgegenprobe/ Reverse/ Contre epreuve/ Test indiretto/ Test indirecto	A <sub>1</sub>											
F		A <sub>2</sub>											
G		B											
H		0											
Proband/Specimen/Echantillon/Campione/Probando													
Ergebnis/Result/Resultat/Risultato/Resultado													