

Erytype S ABD+Rev. A₁, B

REF 806127

187976 / 25

Ergebnisprotokoll / Worksheet / Feuille de travail /
 Protocollo risultati / Hoja de control

Mikrotestplatte (Streifen) zur AB0- und D-Bestimmung von 12 Probanden
 Microtestplate (strips) for AB0 and D typing of 12 specimens
 Microplaque (bandes) pour la détermination AB0 et D de 12 patients
 Micrometodo (strips) su piastra per la determinazione dei gruppi sanguigni AB0 e D di 12 campioni
 Microplaca (tiras) para la tipificación del grupo AB0 y D en 12 probandos

Ch.-B. / Lot-No. / Lot Nr. / Lotto / Lote:	
Datum / Date / Data / Fecha:	
Platte Nr. / Plate No. / Plaque No. / Piastra No. / Placa No.:	
Unterschrift / Signature / Firma:	

		1	2	3	4	5	6	7	8	9	10	11	12
A	Anti-A (Clone: A003)												
B	Anti-B (Clone: B005)												
C	Anti-AB (Clones: BS63/BS85)												
D	Anti-D (Clone: BS226)												
E	Anti-D (Clone: BS232)												
F	Neg.-Kontr./Neg.-Contr./Contr. Neg.												
G	Serumgegenprobe/ Reverse/ Contre epreuve/ Test indiretto/ Test indirecto	A ₁											
H		B											
Proband/Specimen/Echantillon/Campione/Probando													
Ergebnis/Result/Resultat/Risultato/Resultado													

