


Biotestcell - Di^a

REF 816031
LOT 3412011-00
 2024-05-20

187501 / 10

Last Name		First Name		Born				
Bloodgroup		Identified Antibodies		ID				
Hospital / Physician				Date		Signature		

			Rh-hr						Kell						Duffy		Luth-eran		Kidd				MNS				Lewis		P	Xg	Colton		Diego		Special Type**		Test Results			
	Donor- No.	Rh- Typ	D	C	E	c	e	C ^w	K	k	Kp ^a	Kp ^b	Js ^a	Js ^b	Fy ^a	Fy ^b	Lu ^a	Lu ^b	JK ^a	JK ^b	M	N	S	s	Le ^a	Le ^b	P ₁	Xg ^a	Co ^a	Co ^b	Di ^a	Di ^b								
Di ^a	550011	R ₁ r	+	+	o	+	+	o	o	+	o	+	o	+	+	o	o	+	o	+	+	+	+	o	+	+	+	+	o	o	+	nt								
Autocontrol																																								

+ = positive 0 = negative nt = not tested w = weak s = strong

** Donors are tested for HLA class I antigens according to Bio-Rad test method using a monoclonal anti-HLA class I reagent. Positive reactions are reported accordingly.

