


# Biotestcell - Di<sup>a</sup>

<b>REF</b>	816031
<b>LOT</b>	3436011-00
	2024-11-04

187501 / 10


Last Name		First Name		Born	
Bloodgroup		Identified Antibodies		ID	
Hospital / Physician				Date	Signature

		Rh-hr					Kell					Duffy		Luth-eran		Kidd		MNS				Lewis		P	Xg	Colton		Diego		Special Type**	Test Results				
Donor-No.	Rh-Typ	D	C	E	c	e	C <sup>w</sup>	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Js <sup>a</sup>	Js <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Lu <sup>a</sup>	Lu <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	M	N	S	s	Le <sup>a</sup>	Le <sup>b</sup>	P <sub>1</sub>	Xg <sup>a</sup>	Co <sup>a</sup>	Co <sup>b</sup>	Di <sup>a</sup>	Di <sup>b</sup>				
Di <sup>a</sup>	990017	R <sub>2</sub> r	+	o	+	+	+	o	o	+	o	+	o	+	+	o	o	+	+	o	+	+	o	+	o	o	nt	nt	+	nt					

Autocontrol																																		
-------------	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--

+= positive    0 = negative    nt = not tested    w = weak    s = strong

\*\* Donors are tested for HLA class I antigens according to Bio-Rad test method using a monoclonal anti-HLA class I reagent. Positive reactions are reported accordingly.

 **Bio-Rad Medical Diagnostics GmbH**  
 Industriestr. 1, D-63303 Dreieich, Germany  
 Tel.: +49-6103-3130-0, Fax: +49-6103-3130-724  
 www.bio-rad.com, techsupport.bmd@bio-rad.com

