

Biotestcell - Di^a

REF

816031

LOT

3520011-00

2025-07-14

187501 / 10

Last Name		First Name		Born			
Bloodgroup		Identified Antibodies		ID			
Hospital / Physician				Date		Signature	

			Rh-hr						Kell						Duffy		Luth-eran		Kidd		MNS				Lewis		P	Xg	Colton		Diego		Special Type**			Test Results								
	Donor-No.	Rh-Typ	D	C	E	c	e	C ^w	K	k	Kp ^a	Kp ^b	Js ^a	Js ^b	Fy ^a	Fy ^b	Lu ^a	Lu ^b	JK ^a	JK ^b	M	N	S	s	Le ^a	Le ^b	P ₁	Xg ^d	Co ^a	Co ^b	Di ^a	Di ^b												
Di ^a	550014	R ₁ r	+	+	o	+	+	o	o	+	o	+	o	+	+	o	o	+	+	+	+	+	o	+	+	+	+	+	+	o	+	nt												
Autocontrol																																												

+= positive

0= negative

nt= not tested

w= weak

s= strong

** Donors are tested for HLA class I antigens according to Bio-Rad test method using a monoclonal anti-HLA class I reagent. Positive reactions are reported accordingly.

